

# Origin Medical EXAM ASSISTANT™ (OMEA™) Cuts **Scan Time by 33%**, Improves **Completion Rate** and **Capacity**, while providing **Ergonomic Relief for Sonographers**



## About Trial

A two-phase prospective clinical trial evaluating the use of OMEA™, an AI based software system for the automated obstetric ultrasound imaging for first-trimester NT exams (13 diagnostic views, head-to-toe). The task involved capturing 13 diagnostic views, assessing 77 quality criteria, and obtaining 2 measurements (NT and CRL). Conducted across two clinical sites in the US, the trial analyzed 228 examinations and 3,420 acquisitions performed by 8 experienced maternal-fetal medicine sonographers (4-22 years of routine OB imaging experience, currently practicing at OB or MFM offices). Phase 1 compared manual versus OMEA™ assisted workflows to determine time savings and potential capacity gains on using AI; Phase 2 evaluated the benefits of using OMEA™ to boost capacity, revenue, and also assessed the feedback of the sonographers involved in the study.

## Key Findings

- 01 **Performance and Workflow Impact:** Sonographers assisted by OMEA™ achieved **33.6% shorter median scan times (8.17 vs. 12.30 minutes)** while simultaneously improving their **exam completeness from 78.6% to 92.05%**. Crucially, they also achieved **>90% completeness across all BMI ranges (17-60) and scan durations (determined by fetal position, movement, body habitus)**, ensuring imaging consistency and speed regardless of workflow challenges. Notably, there was also a median **reduction of 75% in keystrokes** for each exam (150 to 38).
- 02 **Operational and Capacity Impact:** Standardized scan times of the OMEA™ assisted workflow eliminate scheduling unpredictability, enabling tighter appointment slots, fewer cascading delays, and optimized room/machine utilization. This consistency **allows practices to schedule complex cases (e.g., high-BMI patients) with the same time-certainty as routine cases**, while adding **1-3 additional exam slots (20-30 mins) per sonographer daily within existing shift hours**.
- 03 **Economic Impact:** With additional daily capacity of a median of 2 scans per day per sonographer with OMEA™ assisted workflows, **practices can unlock around 500 additional exams per year per sonographer**, resulting in **additional revenue opportunity of \$125k per full-time equivalent (FTE) sonographer**, without additional spend on staffing, equipment, or facility costs.

# Trial Overview



## Study Design



**Phase 1 (Performance & Workflow Study)** - Manual and OMEA™ assisted workflows for each patient with randomized sequence assignment



**Phase 2 (Capacity and Economic Impact Study)** - Operator focused only on OMEA™ assisted workflow for each patient



## Scanning Protocol (AIUM)

13 diagnostic views, 15 acquisitions per exam, 2 CRLs, 2 NTs, and one of each diagnostic view from head to toe (First trimester NT Scan)



## Sample size

228 patients across 2 sites, 3,420 acquisitions



## Team

8 MFM sonographers (4-22 years of experience)




## BMI of Patients


Ranges between 17-60 (n= 228 patients)

# Key Findings: Performance and Workflow Study


OMEA™ assistance shortened exam times by 33.6% while maintaining an image quality acceptance rate of 92.05%, and with 75% lesser keystrokes for each exam (150 to 38).




Sonographers achieved consistently high completion rates while using OMEA™ within the timeframe, preventing complex exams from causing schedule overruns or requiring repeat visits due to incomplete exams.



**Sonographers consistently achieved >90% exam completion rates, regardless of BMI, while using OMEA™ assistance.**



**Sonographers using OMEA™ experienced shorter exam times than with manual workflows across the full BMI spectrum, supporting efficient, predictable workflows even in high-BMI exams.**







## Performance and Workflow Study (Phase 1) Outcomes

With a median time **savings of 35 to 45 minutes per sonographer per shift** with OMEA™ assisted workflows, **1-2 additional exam slots (20-30 minutes each)** can be added during the same shift, increasing each sonographer's daily capacity by up to 25%.

For Phase 2 of the study, **each sonographer now performed 10 scans/shift (vs. 8 scans/shift in Phase 1)**, with OMEA™ assistance, to assess the impact of the additional capacity.

# Key Findings: Capacity and Economic Impact Study

With a median of 2 additional scans/day/FTE sonographer, a site can perform 500 additional scans annually, unlocking up to \$125K/year in additional revenue per FTE sonographer.

Capacity and Economic Impact Per FTE Sonographer

Number of scans with manual workflow	Number of scans with OMEA™ assisted workflow	Scan time gained per day	Additional slots/day	Additional scans/year*	Additional revenue/year**
8 exams/day	Conservative - 9	25 min to 35 mins	1	250	\$62.5K
	Median - 10	35 mins to 45 mins	2	500	\$125.5K
	Optimistic - 11	45 mins to 60 mins	3	750	\$187.5K

\* 250 working days, \*\*Assuming revenue per scan is \$250

## Sonographer's Experience

### Ergonomic & Productivity Gains

100% of sonographers reported significant time savings, and 75% reported ergonomic relief. Reducing physical strain and cognitive load allowed sonographers to feel less stressed and focus more on the patient experience during and after exams.

- **Reduced Repetition:** A 75% reduction in keystrokes **lowered repetitive clicking and typing.**
- **Less Physical Strain:** Automated image capture **reduced the need to twist, turn, or hold awkward positions to get the perfect image,** even in high-BMI patients.
- **Lower Cognitive Load:** OMEA™ **automatically tracks the protocol, removing the mental burden of maintaining a checklist** while scanning.

## Sonographer's Feedback

“

...Saves clicks, measurements and annotations...



**Sonographer 6**  
(OB - 15 years)

”

“

...It captured anatomy as I was scanning other anatomy that I wasn't necessarily going for at the time...



**Sonographer 2**  
(OB - 4 years)

”

“

...Reduces need for typing or clicking and Time savings...



**Sonographer 8**  
(MFM - 12 years)

”

“

...Capturing images without me needing to annotate or freeze...



**Sonographer 4**  
(OB - 10 years)

”

“

...Capturing structures as I was looking at something else...



**Sonographer 5**  
(MFM - 10 years)

”

“

...Saved time by capturing images quickly...



**Sonographer 3**  
(OB - 22 years)

”

# Explore Collaborative Opportunities with Origin Medical

---

Scan to book a  
meeting with us



Access our research  
and clinical insights



Learn more about  
Origin Medical



## Contact Us



[www.originmedical.ai](http://www.originmedical.ai)



[marketing@originmedical.ai](mailto:marketing@originmedical.ai)



One Broadway, 14th Floor,  
Cambridge, MA 02142