

Total control.


Absolute positioning freedom and ultra-high definition intraoral imaging. RX DC - eXTend technology simplifies your work. Total, wireless control ensures fast installation and adaptation to all possible space requirements.

MyRay, just right for you.



TECHNICAL DATA	
Generator	Constant potential, microprocessor-controlled
Working frequency	145 - 230 KHz with self-adjustment (typically 175 KHz)
Focal spot	0.4 mm (IEC 336)
Total filtration	2 mm @ 60 kV / 2 mm @ 65 kV / 2 mm @ 70 kV (*)
Anode current	4 / 8 mA
Voltage at X-ray tube	60 / 65 / 70 kV (*)
Exposure times	0.020 – 1.000 seconds, R'10 and R'20 scale
Source-skin distance	20 and 30 cm
Irradiated field	Ø 60 mm and Ø 55 mm (with round cone)
Additional collimators	35 x 45 mm (with rectangular cone for size 2 sensors), 31 x 41 mm and 22 x 35 mm, for size 1 and size 0 sensors
Power supply	50/60 Hz, 115-120 V AC ±10% or 230-240 V AC ±10%
Duty Cycle	Continuous operation with self-adjustment up to 1s/80s total
Arms (for Standard version only)	Available in 3 lengths: 40 cm – 60 cm – 90 cm
Max. arm extension	230 cm, from wall
Versions	Standard (wall mounted) or Mobile (on portable cart)
Dose delivered	Viewing on a handheld device with option of digital archive on PC via iRYS software which can be automated via the "RX DC connect" (optional) accessory
PC connection cable	Serial with USB adapter available in various lengths
SOFTWARE	
Acquisition software (for PC)	iCapture for automatic saving of RX DC exposure parameters on PC
Image management software (for PC)	iRYS (compliant with ISDP®10003:2020 in accordance with EN ISO/IEC 17065:2012 - certificate number 2019003109-2) and iPad iRYS viewer App (free)
Protocols supported in iRYS	ICOM 3.0, TWAIN, VDDS
DICOM Node Connectivity	iRYS - IHE compliant (Print; Storage Commitment, SR document; WorkList; MPPS; Query/Retrieve)
X-ray log	iRYS feature to associate exposure parameters with the X-ray images of each examination (exportable in PDF or CSV format)
MINIMUM SYSTEM REQUISITES	
Supported operating systems	Microsoft® Windows® 10, 11 Professional 64 bit
Processor	Intel Core i3 or higher
Hard Disk	100 GB SSD (250 GB recommended)
RAM	4 GB (8 GB recommended)
Graphics card	Discrete 3D Video Card or integrated GPU
Display settings	1920x1080 pixel 24-bit RGB Full HD
Power supply	Use a power adapter of a power suitable for the video card in use
Port	USB 2.0 or later versions

(*) values depend on the country where the product is marketed.



myray
new comfort
in digital imaging

www.my-ray.com



BU Medical Equipment

Plant - Via Bicocca, 14/c - 40026 Imola - Bo (Italy) tel. +39 0542 653441 - fax +39 0542 653555

Headquarters - Cefla s.c. Via Selice Provinciale, 23/a - 40026 Imola - Bo (Italy) tel. +39 0542 653111 - fax +39 0542 653344

Cefla Medical North America

6125 Harris Technology Blvd. - Charlotte, NC 28269 - Ph: 704 598 0020 - www.ceflamedicalna.com - info@cefladental.com

01/2024
MRXEGB171S00
Data may be subject to change without notice.

According to the regulations in force, some products and/or features may have different availability and characteristics in areas outside of the European Union. Please contact your local distributor.



RX DC
X-ray unit
with eXTend technology




myray
new comfort
in digital imaging


Precision diagnosis.

Maximum image quality, minimum dose for the patient.
RX DC - eXTend technology provides always-sharp images,
a full configuration range and the exclusive flexibility of
wireless technology.

IMAGE BLURRING
FOCAL SPOT
0,8 mm




OPTIMAL DEFINITION
FOCAL SPOT
0,4 mm



The DC generator in the head tube is high-frequency and constant-potential. This technology gives sharp images with greater detail and lower exposure times than would be attainable with AC X-ray units, which are characterised by variable emissions. Moreover, constant-potential design ensures image generation is unaffected by power fluctuations. RX DC - eXTend technology is reliable for all diagnostic needs and always provides high-definition images by adapting to the sensor type.

Efficient and reliable real-time imaging.



0,8 mm

0,4 mm

ALWAYS-SHARP IMAGES


The focal spot of just 0.4 mm is placed in the tube head in such a way as to obtain a source-to-skin gap of 30 cm (total bulk remaining equal). In this way RX DC - eXTend technology implements extensive internal collimation of the X-rays and gives an extremely small focal spot, producing ever-sharper images and ever-more precise detail.

MAXIMUM PRECISION
Focal spot 0.4 mm
and power 70 kV, 8 mA.



FAST INSTALLATION AND WIRELESS CONTROL

The efficiency of wireless technology with maximum simplicity of use. The wireless controller frees users from the limits posed by on-machine control panels or wall-mounted controls. It is equipped with a button for ultra-fast shooting (fraction of a second) and two simple settings which make it easy to select the most suitable X-ray acquisition programme.




MULTI-MODE

Automatic parameter modulation ensures the best exposure power/time selection: parameters are automatically determined on the basis of the patient's build and the specific region of investigation. With 28 selectable sensitivity levels, sharp images are guaranteed with any sensor.




MINIMUM X-RAY DOSE

Attention to patient health is meticulous: a high frequency, constant potential generator minimises exposure times and reduces harmful radiation. Where deemed appropriate, the 4 mA mode halves the amount of X-rays. The interchangeable rectangular collimator cone (at 30 cm) further reduces the irradiated body surface area by adapting it to the effective surface area of the sensor.



SEQUENTIAL EXPOSURE

No downtimes as a result of tube overheating, not even when repeated use is required. The fast dynamic duty cycle allows, in fact, sequential exposures by keeping tube temperature under constant control on the large hand-held unit display.



RX DC CONNECT (optional)

The RX DC X-ray unit can easily be connected to your PC via RX DC CONNECT. Via the USB port, you can log the X-ray exposure dose data in digital format. With iRYS you can add the image to the patient's record and the relative X-ray log. Monitor the dose value over time, display and export to other applications via shareable file.




SIMPLE INSTALLATION, VERSATILITY, RELIABILITY

The solid arms are made of high quality materials that ensure strength and durability while reducing the risk of accidental vibration during acquisition. They are available in lengths of 40 cm, 60 cm and 90 cm and can be pointed in 6 directions to provide maximum adaptability and simplicity of installation.

TECHNICAL DATA	
Generator	Constant potential, microprocessor-controlled
Working frequency	145-230 KHz (typically 175 KHz)
Focal spot	0.4 mm (IEC 336)
Anode current	4 / 8 mA
Voltage at X-ray tube	60 / 65 / 70 kV (*)
Exposure time	0.020 – 1.000 seconds, R'10 and R'20 scale
Source-skin distance	20 and 30 cm
Irradiated field	35 x 45 mm (with rectangular cone for size 2 sensors), Ø 60 mm and Ø 55 mm (with round cone)
Additional collimators	31 x 41 mm and 22 x 35 mm, for size 1 and size 0 sensors
Total filtration	2 mm @ 60 kV / 2 mm @ 65 kV / 2 mm @ 70 kV (*)
Power supply	50/60 Hz, 115-120 V AC ±10% or 230-240 V AC ±10%
Duty Cycle	Continuous operation with self-adjustment up to 1s/80s total
Stability	Automatic lock/release, with touch-sensitive activation (HyperSphere technology)
Arms	Available in 3 lengths: 40 cm - 60 cm - 90 cm
Maximum arm extension	230 cm, from wall
Dose delivered	Viewing on a handheld device with option of digital archive on PC via iRYS software which can be automated via the "RX DC connect" (optional) accessory
PC connection cable	Serial with USB adapter available in various lengths
SOFTWARE	
Acquisition software (for PC)	iCapture for automatic saving of RX DC exposure parameters on PC
Image management software (for PC)	iRYS (compliant with ISDP®10003:2020 in accordance with EN ISO/IEC 17065:2012 - certificate number 2019003109-2) and iPad iRYS viewer App (free)
Protocols supported in iRYS	ICOM 3.0, TWAIN, VDD5
DICOM Node Connectivity	iRYS - IHE compliant (Print; Storage Commitment, SR document; WorkList; MPPS; Query/Retrieve)
X-ray log	iRYS feature to associate exposure parameters with the X-ray images of each examination (exportable in PDF or CSV format)
MINIMUM SYSTEM REQUISITES	
Supported operating systems	Microsoft® Windows® 10, 11 Professional 64 bit
Processor	Intel Core i3 or higher
Hard Disk	100 GB SSD (250 GB recommended)
RAM	4 GB (8 GB recommended)
Graphics card	Discrete 3D Video Card or integrated GPU
Display settings	1920x1080 pixel 24-bit RGB Full HD
Power supply	Use a power adapter of a power suitable for the video card in use
Port	USB 2.0 or later versions

(*) values depend on the country where the product is marketed.



www.myray.com



BU Medical Equipment

Plant - Via Bicocca, 14/c - 40026 Imola - Bo (Italy) tel. +39 0542 653441 - fax +39 0542 653555

Headquarters - Cefla s.c. Via Selice Provinciale, 23/a - 40026 Imola - Bo (Italy) tel. +39 0542 653111 - fax +39 0542 653344

Cefla Medical North America

6125 Harris Technology Blvd. - Charlotte, NC 28269 - Ph: 704 598 0020 - www.ceflamedicalna.com - info@cefladental.com

Data may be subject to change without notice. 01/2024 MRX-HGB171500
According to the regulations in force, some products and/or features may have different availability and characteristics in areas outside of the European Union. Please contact your local distributor.



RX DC
X-ray unit
with HyperSphere technology



**FREE TO MOVE**

Maximum freedom of movement with innovative ball joint (Patented).

RX DC

Hyper Technology.

Innovative design, revolutionary ergonomics, advanced technology. RX DC - HyperSphere technology brings the best of DC X-ray units into your surgery.

MyRay, just right for you.



MAXIMUM FLEXIBILITY

The wireless remote controller, the multi-mode option and the 28 adjustment levels (depending on sensor sensitivity) ensure full adaptability whatever your operating requirements.



PRECISION X-RAY IMAGING

A constant potential head tube (8 mA) with a tiny focal spot (0.4 mm at 30 cm) produces optimal images under all circumstances.



INNOVATIVE ERGONOMICS

The RX DC unit features HyperSphere technology which, thanks to the full-swivel ball joint, can reach any position with ease.



COMPLETE RELIABILITY

Built from high quality materials and featuring a comprehensive array of equipment. Versatile and easy to install, this X-ray unit is reliable whatever the situation.

**SIMPLE AND IMMEDIATE: WIRELESS**

The wireless remote controller lets the user control the device (by communicating with the X-ray tube) while enjoying full freedom of movement. Access to exposure programmes is provided via two simple settings. The large display shows the sequential exposure monitor and the patient exposure dose; moreover, the controller has a wireless X-ray **snapshot** button. Wireless device control allows fast, easy installation: no fixed control

panels are required, thus providing greater freedom when positioning the X-ray unit.

**MECHANICAL RELIABILITY**

The solid, light arms feature an effective, integrated self-balancing system that reduces any risk of tube head vibration during image acquisition.

Hyper Ergonomy.

RX DC - HyperSphere technology allows attainment of any position with ease thanks to the revolutionary ball joint. Outstanding ergonomics ensures all your diagnostic needs are met effortlessly.

HyperSphere technology gives the RX DC unit full rotation capability. The tube revolves freely around the joint, allow it to reach practically any position, including the vertical.

RX DC - HyperSphere technology also features an automatic touch-sensitive device for simple, efficient locking/release of the X-ray head tube so it can be repositioned effortlessly between one exposure and the next. Ergonomic zones on the sides of the head provide a firm grip for effective positioning.

Extensive positioning.

- electro-brake with touch-sensitive control
- infinite position range
- maximum versatility
- complete reliability

EASY TOUCH

Effortless locking, release and repositioning.

COMFORT POSITIONING

Two comfortable handles for a secure grip and the best possible positioning.



INFINITE POSITIONS, INFINITE DIAGNOSTIC CAPABILITY

Diagnosis with unlimited movement thanks to the revolutionary ball joint which allows simple yet precise head repositioning and effortless attainment of even the trickiest positions.




INSTALLATION VERSATILITY

The extruded aluminium arms - available in lengths of 40 cm, 60 cm and 90 cm to ensure outstanding installation versatility - are equipped with an integrated self-balancing system. Solid and light, they can be pointed in any direction and reduce any risk of tube head vibration during image acquisition.



Hyper Performance.

In RX DC - HyperSphere technology, advanced ergonomics, technological innovation and revolutionary design merge to provide users with ultra-sharp images.





Sharp images at all times RX DC - HyperSphere technology provides your surgery with optimum X-ray quality whatever the type of sensor connected. Now even more powerful, with 70 kV and 8 mA, even more flexible and suitable for all commercially available sensors. The constant potential head tube, associated with the smallest intraoral imaging focal spot available (0.4 mm), ensures the best images whatever your diagnostic needs.

High definition diagnostic.


MAXIMUM QUALITY

With a tiny focal spot of 0.4 mm (at 30 cm), RX DC - HyperSphere technology produces sharp images under any condition. The tube head is now even more powerful as it operates at 70 kV, 8 mA. RX DC - HyperSphere technology gives your surgery the precision and quality of cutting-edge know-how.




MULTI-MODE

Maximum flexibility to meet your diagnostic needs. Automatic parameter modulation ensures exposure power and time are always selected according to the patient's build and the specific region of investigation.



SEQUENCED EXPOSURE

The dynamic service cycle allows uninterrupted use of the RX DC, as in the case of systematic examinations, and real-time monitoring of tube head temperature on the large wireless controller display.



RX DC CONNECT (optional)

The RX DC X-ray unit can easily be connected to your PC via RX DC CONNECT.

Via the USB port, you can log the X-ray exposure dose data in digital format.

With iRYS you can add the image to the patient's record and the relative X-ray log.

Monitor the dose value over time, display and export to other applications via shareable file.

MINIMUM DOSE

The constant potential high frequency (DC) generator reduces the most harmful low energy radiation that is characteristic of analogue (AC) generators: current is adjustable (from 8 mA to 4 mA), as are exposure times. Moreover, the long cone (30 cm) with incorporated rectangular collimator reduces the exposed surface area. This maximises image quality and safeguards patient and worker health.



INCLUDES

RECTANGULAR COLLIMATOR

Reduces the exposed body area and allows current to be adjusted from 8 mA to 4 mA.

Versatile and adaptable.


Wall-mounted with variable positions or in a mobile cart-mounted version (to be shared among multiple workstations), RX DC is extremely versatile and easily adapts to all your working needs.

MyRay, just right for you.



TECHNICAL DATA	
Generator	Constant potential, microprocessor-controlled
Working frequency	145 - 230 KHz with self-adjustment (typically 175 KHz)
Focal spot	0.4 mm (IEC 336)
Total filtration	2 mm @ 60 kV / 2 mm @ 65 kV / 2 mm @ 70 kV (*)
Anode current	4 / 8 mA
Voltage at X-ray tube	60 / 65 / 70 kV (*)
Exposure times	0.020 – 1.000 seconds, R'10 and R'20 scale
Source-skin distance	20 and 30 cm
Irradiated field	Ø 60 mm and Ø 55 mm (with round cone)
Additional collimators	35 x 45 mm (with rectangular cone for size 2 sensors), 31 x 41 mm and 22 x 35 mm, for size 1 and size 0 sensors
Power supply	50/60 Hz, 115-120 V AC ±10% or 230-240 V AC ±10%
Duty Cycle	Continuous operation with self-adjustment up to 1s/80s total
Arms (for Standard version only)	Available in 3 lengths: 40 cm – 60 cm – 90 cm
Max. arm extension	230 cm, from wall
Versions	Standard (wall mounted) or Mobile (on portable cart)
Dose delivered	Viewing on a handheld device with option of digital archive on PC via iRYS software which can be automated via the "RX DC connect" (optional) accessory
PC connection cable	Serial with USB adapter available in various lengths
SOFTWARE	
Acquisition software (for PC)	iCapture for automatic saving of RX DC exposure parameters on PC
Image management software (for PC)	iRYS (compliant with ISDP®10003:2020 in accordance with EN ISO/IEC17065:2012 - certificate number 2019003109-2) and iPad iRYS viewer App (free)
Protocols supported in iRYS	ICOM 3.0, TWAIN, VDDS
DICOM Node Connectivity	iRYS - IHE compliant (Print; Storage Commitment, SR document; WorkList; MPPS; Query/Retrieve)
X-ray log	iRYS feature to associate exposure parameters with the X-ray images of each examination (exportable in PDF or CSV format)
MINIMUM SYSTEM REQUISITES	
Supported operating systems	Microsoft® Windows® 10, 11 Professional 64 bit
Processor	Intel Core i3 or higher
Hard Disk	100 GB SSD (250 GB recommended)
RAM	4 GB (8 GB recommended)
Graphics card	Discrete 3D Video Card or integrated GPU
Display settings	1920x1080 pixel 24-bit RGB Full HD
Power supply	Use a power adapter of a power suitable for the video card in use
Port	USB 2.0 or later versions

(*) values depend on the country where the product is marketed.



www.my-ray.com

cefla **BU Medical Equipment**
Plant - Via Bicocca, 14/c - 40026 Imola - Bo (Italy) tel. +39 0542 653441 - fax +39 0542 653555
Headquarters - Cefla s.c. Via Selice Provinciale, 23/a - 40026 Imola - Bo (Italy) tel. +39 0542 653111 - fax +39 0542 653344
Cefla Medical North America
6125 Harris Technology Blvd. - Charlotte, NC 28269 - Ph: 704 598 0020 - www.ceflamedicalna.com - info@cefladental.com

Data may be subject to change without notice. 01/2024 MRXDCB171500
According to the regulations in force, some products and/or features may have different availability and characteristics in areas outside of the European Union. Please contact your local distributor.



RX DC
High frequency
X-ray unit






MAXIMUM PRECISION
Focal spot 0.4 mm
and power 70 kV, 8 mA

ALWAYS THE BEST ACQUISITION PROGRAMME

Sharp images with simple and immediate configuration. The controller has two settings with which to select the programme most suitable for optimal X-ray image capture. The large display lets users monitor the temperature of the head tube which, thanks to the dynamic duty cycle, allows sequential exposure.


Highest quality with lowest exposure.

Always-sharp images, versatility and meticulous attention to patient health. With RX DC you get the best DC technology with the lowest X-ray dose.



The constant potential high frequency generator (DC) provides sharp images with the very highest level of detail. Compared to AC systems, they also reduce exposure times and the amount of harmful radiation by containing the dose administered to the patient. A focal spot of just 0.4 mm - one of the smallest available - ensures images are always sharp and of the highest quality.

High definition real-time imaging.



IMMEDIATE CONFIGURATION

Multi-Mode automatic exposure parameter modulation always ensures optimal time and power selection. Parameters are, in fact, adjusted automatically according to patient build and the region under investigation.



MINIMUM IRRADIATION

Attention to patient health is meticulous thanks to the constant potential DC generator with adjustable power (from 8 to 4 mA). Moreover, rectangular collimators can be used: these reduce the irradiated body area and thus lower the dose received by the patient.





OPTIMAL ERGONOMICS

The ergonomic handle is designed to maximise grip comfort and ensure easy, stable positioning of arms and tube head. A protractor with a graduated scale allows optimal repositioning of the tube head.


RX DC CONNECT (optional)

The RX DC X-ray unit can easily be connected to your PC via RX DC CONNECT. Via the USB port, you can log the X-ray exposure dose data in digital format. With iRYS you can add the image to the patient's record and the relative X-ray log. Monitor the dose value over time, display and export to other applications via shareable file.



PRECISION DIAGNOSTICS

Superb image definition: sharp edges and excellent detail. An embedded collimator cone gives a source-to-skin distance of 30 cm. This increases X-ray parallelism, providing more precise images, lower doses and ensuring greater attention to patient health. RX DC offers maximum flexibility and optimum X-ray quality whatever the type of sensor connected.



SIMPLE INSTALLATION, VERSATILITY, RELIABILITY

RX DC provides outstanding adaptability and simplicity of installation thanks to arms with an integrated self-balancing system that can be pointed in 6 directions - available in lengths of 40 cm, 60 cm and 90 cm. All parts are made from materials of only the finest quality to minimise maintenance costs and reduce the risk of accidental vibration during acquisition.