

Solubility measurement of crystalline and LLPS forms by extrapolation from values measured in water-methanol mixtures

R6182

Kyoowon Baek ^[1], John Comer ^[1], Karl Box ^[1], Jon Mole ^[2], Jaechun Woo ^[3]

[1] Sirius Analytical Ltd. Forest Row, RH18 5DW, UK. [2] Sirius Analytical Inc. Sirius Analytical Inc., Beverly, MA, USA. [3] Korea Research Institute of Chemical Technology, Daejeon, Korea

Purpose

To verify that solubility of ionizable drugs measured in aqueous solution is correctly predicted from solubility measured in water-methanol mixtures.

Methods

Four ionizable drugs were chosen for this study: bendroflumethiazide (acid) and cinnarizine, loratidine and clotrimazole (bases). Their aqueous pK_a values were determined by Yasuda-Shedlovsky extrapolation from p_sK_a values measured in water-methanol mixtures. Their aqueous solubilities were measured by pH-metric techniques (CheqSol or Curve Fitting) in which the aqueous pK_a was used in the calculation of results from the experimental data. Their solubilities were also measured in water-methanol mixtures, in which p_sK_a values interpolated from the Yasuda-Shedlovsky graphs were used in the calculation of results from the experimental data. The pH-metric solubility experiments lasted typically for 1 – 2 hours.

Results

Bendroflumethiazide crystallizes soon after precipitation. Cinnarizine initially precipitates in a non-crystalline form but converts within minutes to a crystalline form. It was the solubility of the crystalline forms of cinnarizine and bendroflumethiazide that the pH-metric experiments measured. Loratidine and clotrimazole are known to undergo Liquid-Liquid Phase Separation (LLPS) and to precipitate initially in a glassy, non-crystalline form that can persist for several hours, and it was the solubility of this non-crystalline form that the pH-metric experiments measured. Measured aqueous solubilities agreed well with results obtained by extrapolation. For all four compounds the graphs of log (pH-metric solubility) vs. methanol weight % were linear with R^2 of 0.986 or better


<i>Results at 25 °C</i>	pK_a (1)	pK_a (2)	Extrapolated solubility			Crystalline solubility		
			M	log S (M)	µg/mL	M	log S (M)	µg/mL
Bendroflumethiazide	8.33	9.86	0.000029	-4.54	12.0			
Cinnarizine	2.41	7.53	0.000001	-6.09	0.3			
Clotrimazole	5.80		0.000039	-4.41	13.4	0.000001	-5.93	0.4 [1]
Loratidine	5.26		0.000036	-4.44	13.7	0.000004	-5.38	1.6 [1]

Conclusion


Solubilities for the four compounds measured in aqueous solution were predicted by extrapolation from solubilities measured in water-methanol mixtures. The pH-metric experiments revealed information about the transition of precipitated compounds between solid state forms.

Reference

[1] Hsieh, Y. L.; Ilievbare, G. A.; Van Eerdenbrugh, B.; Box, K. J.; Sanchez-Felix, M. V.; Taylor, L. S., pH-Induced Precipitation Behavior of Weakly Basic Compounds: Determination of Extent and Duration of Supersaturation Using Potentiometric Titration and Correlation to Solid State Properties. Pharm Res 2012, 29 (10), 2738-2753.



Clotrimazole is an antifungal, and **loratidine** is an H1 histamine antagonist. When they precipitate in response to pH change, they undergo Liquid-Liquid Phase Separation (LLPS). The LLPS form comprises non-crystalline droplets of supercooled liquid. It is the solubility of the LLPS that has been measured here. Shake-flask measurements of the original solid indicate a much lower solubility corresponding to the crystalline form [1]. Further work will focus on the determination of the time taken to convert from LLPS to crystals using controlled supersaturation experiments



The pK_a values and solubilities of these samples were measured on SiriusT3 instruments. Bendroflumethiazide and cinnarizine were measured by Karl Box and co-workers at Sirius UK, and clotrimazole and loratidine were measured by Kyoowon Baek using an instrument in Dr. Jaechun Woo's laboratory at KRICT.



Photo: SiriusT3.

For more information please contact:

Sirius Analytical, Forest Row Business Park, Station Road, Forest Row, East Sussex, RH18 5DW, UK

T +44 (0)1342 820720 F +44 (0)1342 820725 E sales@sirius-analytical.com

Sirius Analytical - World-class provider of innovative instrumentation for chemical analysis

www.sirius-analytical.com