

Application of MacroFLUX™ apparatus for screening formulations before bioequivalence studies

Enikő Borbás¹, Konstantin Tsinman², György Marosi¹, Zsombor K. Nagy¹, Bálint Sinkó²

¹Budapest University of Technology and Economics, Budapest, Hungary, ²Pion Inc, 10 Cook Street, Billerica, MA 01821, USA



PURPOSE

For **generic drug development** traditional (USP) dissolution tests have been used in the pharmaceutical industry to compare performance of different drug product formulations before or instead of conducting bioequivalence studies. Although **dissolution tests** provide a simple way of testing formulations, the **in vivo predictive power** of these tests are questionable. Namely, when a poorly water-soluble API is formulated to enhance its dissolution, additives, such as surfactants and polymers have an effect not only on dissolution profile, but also on flux through the membrane. The aim of this study was to represent the importance of simultaneous dissolution-absorption studies using **MacroFLUX™** apparatus before conducting **bioequivalence studies**.

MATERIALS and METHODS




Figure 1. Schematic of Dissolution – Flux device used in this study.

Brand and generic formulations of Telmisartan, an antihypertensive drug, were tested using **MacroFLUX™**. Receiver chamber integrated with permeation membrane, overhead stirrer and fiber optic (FO) UV probe was inserted in the standard **900 mL** vessel of **USP II** apparatus. A filter-supported **artificial membrane** with **3.8 cm²** area was separating the dissolution (donor) compartment from the receiver compartment contained **15 mL** of pH 7.4 (Prisma™ HT, Pion Inc). The experiment began in **850 mL** at pH 1.6 simulating gastric conditions and then after 30 min media in the dissolution vessel was converted to **FaSSIF** by adding 212 mL of specially formulated concentrate containing SIF powder. The integrated fiber-optic UV probes were positioned in the donor and receiver compartments allowing **real time concentration monitoring** in both chambers. Concentration monitoring was enabled through **fiber optic UV probes** connected to the Rainbow Dynamic Dissolution Monitor instrument (Pion Inc). **Flux (J)** of a drug through a membrane is defined as the amount of drug (m) crossing a unit area (A) perpendicular to its flow per unit time (t).

$$J = \frac{dm}{A \cdot dt} \quad (1)$$

RESULTS

I. Brand name compared to generic formulation containing the same excipients




Figure 2. Dissolution profile (on the left) and appearance profile (on the right) of telmisartan from different formulation and pure API.

The dissolution and flux results of the **brand name** (Micardis) and **generic** (Pritor) Telmisartan 40 mg tablets were compared (Figure 2). Both formulations showed slow release kinetics in SGF and instant dissolution after media conversion to FaSSIF full with the final concentration around **35 µg/mL** (more than 90% of the API dissolved). In the first 30 minutes of the experiments the concentration of the API in the acceptor chamber was under the detection limit ($\sim 0.1 \mu\text{g/mL}$). After media change in the time interval of 50-120 minutes the flux through membrane was found to be **$0.337 \pm 0.028 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$** in case of the brand and **$0.308 \pm 0.014 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$** in case of the generic product. In comparison the initial flux from unformulated API at the same load in the donor compartment was **$0.120 \pm 0.016 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$** and it was decreasing due to **precipitation** of API in FaSSIF.

II. Brand name compared to generic formulation containing different excipients




Figure 3. Dissolution profile (on the left), in vitro appearance profile (in the middle) and in vivo appearance (on the right) of telmisartan from different formulations.

The dissolution and flux results of the **brand name** (Micardis) and **generic** (Actavis) Telmisartan 40 mg tablets were compared (Figure 3). Actavis showed a slower release kinetics than Micardis, though reached the same maximum concentration after 110 min. After media change the flux from the generic product was found to be **$0.240 \pm 0.011 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$** , which is only **71%** of the flux of the brand name (**$0.337 \pm 0.028 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$**). This *in vitro* result showed excellent correlation with the *in vivo* data from bioequivalence studies¹, where the appearance rate or the drug in blood from Actavis was **72 %** of the rate from Micardis.

RESULTS

III. Brand name with additional solubility modifier formulation additives




Figure 4. Dissolution profile (on the left) and flux (in the middle) and solubility (on the right) of telmisartan from Micardis tablet with additional solubility modifiers

The effect of solubility modifier formulation additives, such as **Soluplus**, **PVP VA 64** were studied by adding excipients (80 mg) to the donor chamber. Although the addition of excipients to the donor media had **no significant effect on the dissolution** kinetics (Figure 4) of the API resulting in more than 90% dissolved drug in FaSSIF media, the **flux** through membrane was significantly **decreasing**. In case of Soluplus containing media the flux was found to be $0.252 \pm 0.017 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$ while in case of PVP VA 64 containing media $0.274 \pm 0.003 \mu\text{g}/(\text{cm}^2 \cdot \text{min})$. Flux results can be better understood considering that flux not only depends on the concentration gradient between the two chambers, but also **inversely proportional to the solubility²** of the API as shown in Figure 4.

CONCLUSION

The results showed that by adding **solubility modifier formulation excipients** the **dissolution** of Telmisartan was not altered, but the **flux** through the membrane was found to be **decreasing significantly**. These results point out the limitations of traditional (USP) dissolution tests and emphasize the importance of **simultaneous dissolution-absorption** studies, where the effect of formulation excipients on dissolution and also on membrane transport can be measured. This knowledge is essential for **generic formulation development** before **bioequivalence studies** are conducted.

The **in vivo predictive power** of the **simultaneous dissolution-absorption** test was demonstrated by comparing the *in vitro* fluxes to *in vivo* rate of appearance in blood of brand name and generic formulation of telmisartan. The performance of the generic product was 71% of the brand *in vitro*, while 72 % *in vivo*, showing excellent **in vitro-in vivo correlation**.

REFERENCE

- ¹Public Assessment Report of Actavis
- ²Borbás, et al. "Investigation and mathematical description of the real driving force of passive transport of drug molecules from supersaturated solutions." Mol. Pharm. (2016)