

# A novel liquid biopsy assay with 88% detection rate of genomic alterations across CNS tumors

Matthew G. Varga\*, Christian Klosowski\*, Xavier Bower\*, Naomi Searle\*, David Tsao\*, Gary Palmer\* and Santosh Kesari\*\*

\* BillionToOne, Inc., Menlo Park, CA  
\*\* Pacific Neuroscience Institute, Santa Monica, CA

#2077



NORTHSTAR

BY BILLIONTOONE

## BACKGROUND

Non-invasive detection of challenging tumors like those in the CNS is needed

Genomic assessment of central nervous system (CNS) tumors, including glioblastomas (GBM), is necessary to inform clinical decision making, especially for genes which have prognostic and treatment implications such as *IDH1/2*.

Tissue biopsy remains the current standard for histological and genomic analyses, however there are limitations:

- **Tumor location** renders sample extraction infeasible, or when available, limited in size
- **High heterogeneity** in tumors such as GBMs makes it difficult to fully capturing the diversity of its genomic profile

Plasma-based liquid biopsies have emerged as a complement or alternative to tissue. However, ctDNA from CNS tumors can be difficult to detect due to:

- **Low ctDNA shed rates** inherent to CNS tumor biology
- **Blood brain barrier** impedes ctDNA permeability into peripheral blood

These factors make CNS tumor-derived ctDNA difficult to detect, with current literature suggesting 27 to 55% detection for any alteration, including variants of uncertain significance (VUS):

| Study                                | Patients (n) | ctDNA Detection Rate* |
|--------------------------------------|--------------|-----------------------|
| Schwaederle et al.<br>PMID: 26848768 | 33           | 27%                   |
| Piccioni et al.<br>PMID: 30855176    | 222          | 55%                   |
| Zill et al.<br>PMID: 29776953        | 107          | 51%                   |
| Bagley et al.<br>PMID: 31666247      | 20           | 55%                   |

**Table 1.** Studies reporting ctDNA sequenced by next-generation sequencing (NGS) from plasma obtained from GBM patients.  
\*Includes all alterations, including VUS, pathogenic, and actionable.

Currently, there are no assays on the market approved for use in CNS tumors because of this poor detection rate. Therefore, there is unmet clinical need for a highly sensitive liquid biopsy assay to genomically profile CNS tumors to overcome the above limitations.


## REFERENCES

1. Tsao, D.S., Silas, S., Landry, B.P. et al. A novel high-throughput molecular counting method with single base-pair resolution enables accurate single-gene NIPT. Sci Rep 9, 14382 (2019). <https://doi.org/10.1038/s41598-019-50378-8>
2. Deveson IW, Gong B, Lai K, et al. Evaluating the analytical validity of circulating tumor DNA sequencing assays for precision oncology. Nat Biotechnol. 2021;39(9):1115-1128. doi:10.1038/s41587-021-00857-z
3. LOD defined as lowest concentration of analyte that can be detected >95% probability.

## METHODS

Highly sensitive, plasma-based genomic profiling assay

Northstar Select is a plasma-based, comprehensive genomic profiling test that leverages Quantitative Counting Template™ (QCT) technology<sup>1</sup>, optimized chemistry and panel design, and bioinformatics innovations to increase sensitivity. From January to December 2023, 62 patients with diagnosed CNS tumors were submitted commercially to BillionToOne for Northstar Select analysis.



**Figure 1:** Cited single-nucleotide variant (SNV) limits of detection (LODs) of common assays in the market range from 0.25% to 0.50% VAF<sup>2</sup>. Northstar Select LOD ranges from 0.13-0.16% VAF<sup>3</sup>.

| SNVs / Indels 62 genes |               |       |        |        |
|------------------------|---------------|-------|--------|--------|
| AKT1                   | CCNE1         | EZH2  | JAK2   | MYC    |
| AKT2                   | CD274 (PD-L1) | FANCA | JAK3   | NF1    |
| ALK                    | CDH1          | FBXW7 | KIT    | NOTCH1 |
| APC                    | CDK12         | FGFR1 | KRAS   | NPW1   |
| AR                     | CDK4          | FGFR2 | MAP2K1 | NRAS   |
| ARAF                   | CDK6          | FGFR3 | MAPK1  | NTRK1  |
| ARID1A                 | CDKN2A        | FGFR4 | MAP2K2 | PALB2  |
| ATM                    | CDKN2B        | GATA3 | MEK2   | PDGFRA |
| BRAF                   | CHEK2         | GNAI1 | MET    | PIK3CA |
| BRCA1                  | CTNNB1        | GNAQ  | MLH1   | PMS2   |
| BRCA2                  | DDR2          | GNAS  | MPL    | PTEN   |
| BRIP1                  | EGFR          | HRAS  | MSH2   | PTPRN1 |
| CCND1                  | ERBB2 (HER2)  | IDH1  | MSH6   | RAD51C |
| CCND2                  | ESR1          | IDH2  | MTOR   | RAD51D |


| CNAs: Amplifications 19 genes |       |        | Fusions 9 genes |       |
|-------------------------------|-------|--------|-----------------|-------|
| AR                            | ERBB2 | MET    | ALK             | NTRK1 |
| BRAF                          | ESR1  | MYC    | BRAF            | NTRK2 |
| CCNE1                         | FGFR1 | PDGFRA | FGFR2           | NTRK3 |
| CD274 (PD-L1)                 | FGFR2 | PIK3CA | RET             | ROS1  |
| CDK4                          | KIT   | RAF1   |                 |       |
| CDK6                          | KRAS  | RET    |                 |       |
| EGFR                          |       |        |                 |       |

| CNAs: Losses 5 genes |        | Biomarker |
|----------------------|--------|-----------|
| ATM                  | CDKN2A |           |
| BRCA1                | PTEN   | MSI       |
| BRCA2                |        |           |


**Table 2.** Northstar Select gene panel

## DATA

Increased sensitivity results in high detection rate of clinically actionable alterations in CNS tumors




**Figure 2:** Summary of alterations detected by Northstar Select in CNS tumors



**Figure 3.** (A) Distribution of all variants detected in 62 CNS patient cases (B) Distribution of clinically actionable alterations only among all 62 CNS patient cases.

## CLINICAL CASE STUDIES

Clinically informative alterations of select genes



**Figure 4.** (A) Pathogenic / likely pathogenic alterations detected by Northstar Select (B) Patient-level details on clinically relevant alterations of 3 select genes: *TERT* and *IDH1/2*.

## RESULTS

88% detection rate in CNS tumors, 91% in GBMs only

In this study, Northstar Select's detection rate of any alteration was 88% in CNS and 92% for GBMs only.

When limiting to clinically actionable alterations (i.e. there is a drug, treatment implication, or clinical trial available), 36% of our cases had at least one alteration detected, and in GBMs, this number increased to 49% of cases.

Majority of the alterations found were below 0.5% VAF, which is the reliability threshold of current assays. This is due to the assay's uniquely low LOD of 0.13% to 0.16% VAF, which is below the 0.17% median VAF demonstrated across CNS tumors in this study.

## CONCLUSIONS

Higher detection rates uncover more low-abundance alterations

Current literature indicates an unmet clinical need for a highly sensitive liquid biopsy assay to determine the genomic landscape of CNS tumors and to aid clinical decision making.

Northstar Select has a detection rate of 88.7% in all cases and 91.4% in GBMs, nearly doubling current rates.

Nearly 50% of detected alterations in GBMs were actionable, supporting the clinical utility of Northstar Select to inform clinical decisions either as a complement to imaging and tissue analyses, or independently when tissue biopsies are infeasible.

These data further suggest that the clinical utility of Northstar Select may be expanded into metastatic brain disease and other low ctDNA shedding tumors.