

# Raspberry Pi (RP) and Huawei Modem Quick Guide

## 1. Connect the RP to the ECG machine

- Plug one end of the USB cable into the RP, and the other end into the USB port on the ECG machine (Figure 1).


## 2. Connect the Ethernet cable

- Insert one end of the Ethernet cable into the RP and the other end into the Ethernet port on the ECG machine (Figure 1).

## 3. Connect the modem to the RP

- Plug one end of the short USB cable into the modem and the other end into the RP (Figure 1).
- Alternatively, you can connect the RP to a power source using the power adapter (Figure 1).

**1** Note: RP does not have any physical button for turning on or off.




## 4. Powering on the RP and modem

- First, ensure the ECG machine is turned on.
- Then, connect the RP to a power source. Once powered, the RP will take approximately 2 minutes to fully initialize.
- After the RP is ready, **press and hold the power button on the modem**. You should see a green light (battery status) and a green light (Wi-Fi signal), indicating the modem is functioning correctly.

## 5. ECG machine is ready

- Once the modem is turned ON (Figure 2), please look for the flashing green light in the front of the RP (Figure 3). The machine is now ready to perform a patient's ECG.



# Huawei Modem Quick Guide



## 1. Battery indicator

- **Green Light:**  
Sufficient battery
- **Blinking green:**  
Charging the battery
- **Red:** Low battery

## 2. Power button

## 3. Micro USB port

## 4. Signal indicator

- **Green:** Strong signal
- **Orange:** Weak signal
- **Red:** No signal

## Turning on the Mobile Wi-Fi

Press and hold the power button until indicator lights up, indicating that the Mobile Wi-Fi is turned on.

**Note:** Press and hold the power button to also turn off the Mobile Wi-Fi.

## Battery saving mode

If the Mobile Wi-Fi (Huawei) has been inactive for some time, it enters standby mode and all the indicators are off. Pressing the power button causes the indicators to light up again.

## Charging the battery

If the Mobile Wi-Fi has not been used for a long time or if the battery is fully depleted, the Mobile Wi-Fi may not power on immediately after the charger is plugged in. Please charge the battery for a while before attempting to power on the Mobile Wi-Fi.