

Is the Epitomee® Capsule Safe? Absolutely.



Clinically tested: Proven results for both safety and effectiveness. Most side effects were mild and similar to placebo. Visit epitomeecapsule.com for more safety information.



Non-systemic: It doesn't enter the bloodstream



Drug-free: No hormones, stimulants, or appetite suppressants



Prescribed by your healthcare provider: Ensures it's right for you

Are You a Candidate for the Epitomee® Capsule?

- 1 Indicated for overweight individuals who are 18+ years old, with a minimum BMI of 25 kg/m². The Epitomee® capsule may be considered for adults with or without comorbidities, subject to the judgment of a healthcare provider and in accordance with the instructions for use.
- 2 Individuals who want a non-drug, non-surgical approach.
- 3 Individuals who experienced GLP-1 intolerance.
- 4 Individuals seeking other weight management options after completing medication-based interventions.
- 5 Individuals who struggle with portion control or early satiety and have had little success with diet and exercise.

Less Struggle, More Success

With the Epitomee® Capsule, weight management becomes an intuitive part of daily life—allowing people to feel in control of their health and confident in their bodies for the long run.

Talk to Your Healthcare Provider Today

Learn how the Epitomee® Capsule can help you stay on track safely.



Discover More About the Epitomee® Capsule

Prescriber/Clinic Info

Name
000-000-000

epitomeecapsule.com

Epitomee® Capsule, FDA cleared, Rx Only, is indicated to aid weight management in adults with excess weight or obesity, Body Mass Index (BMI) of 25-40 kg/m², when used in conjunction with diet and exercise. You cannot take Epitomee® if you are pregnant or have had allergic reaction to one or more of the Epitomee® components. Please read the patient information leaflet before use. For more information, please visit our website at epitomeecapsule.com. Manufacturer: Epitomee Medical, 17 Hatochen St. Caesarea 3079892, Israel.

1 Kushner RF, Ard JD, Wadden TA, et al. Quality of Life Improvements Associated with Weight Loss Using a Novel Shape-Shifting Hydrogel Capsule: RESET Study Results. International Journal of Obesity. 2025; In Review. 2 Kamar M, Ryan DH, Leonard S, et al. The safety and efficacy of extended use of an oral shape-shifting superabsorbent hydrogel capsule for weight loss: the ELECT extension study. Obesity Science and Practice. 2025; In Review. 3 Bays HE, Ard JD, O'Neill PM, et al. Weight and cardiometabolic effects of a novel oral shape-shifting superabsorbent hydrogel capsule: Prespecified and exploratory analysis of the Epitomee capsule RESET Study. Obes Pillars. 2025. 4 Jakicic JM, Ryan DH, Ard JD, et al. Association of the Early Response to an Oral Shape-Shifting Superabsorbent Hydrogel Capsule with Weight Loss. Clin Obes. 2025; In Review. 5 Shirin H, Neeland IJ, Ryan DH et al. Effects of an oral biodegradable device used for 12 weeks on weight reduction, cardiovascular risk factors, satiety, snacking, and meal size. Obesity Pillars. 2023. 6 Ard JD, Ryan DH, O'Neill PM, et al. Efficacy and safety of a novel oral hydrogel capsule in adults with overweight or obesity: the pivotal randomized RESET study. Obesity. 2025.

*Refers to the device's mechanism of action.

epitomee®
weight management capsule

The Brilliant Little Weight Management Capsule



The Epitomee® Capsule: Drug-Free, Prescription-Only, Support for Weight Management

It helps you feel full sooner, eat less, and stay on track—without stimulants, hormones, or invasive treatments.

Designed by Science, Made for Real Life

The Epitomee® Capsule is a safe and drug-free weight loss and management tool. Its innovative capsule expands in the stomach, mimicking food to create a feeling of fullness—helping you eat less while still enjoying your meals.

Combined with a healthy lifestyle program consists of diet and physical activity, it supports sustainable weight management and a better quality of life.¹

- Clinically proven and **FDA-Cleared**
- **Works naturally*** in the stomach
- Designed for **use over time**: The safety and effectiveness of Epitomee were studied in the 24 weeks RESET study⁶
- Companion with **diet and physical activity**
- Helps build **healthy, sustainable habits**
- **Easy to use**, take before lunch & dinner




Proven by Clinical Results

Reduce Waistline by 2 Jeans Sizes Over 24 Weeks³

Participants treated with the Epitomee® Capsule reduced their waistline by 2–3 inches (6 cm), which is equivalent to 2 jeans sizes after 24 weeks.

Predictable Weight Loss⁴

Early weight loss with Epitomee® treatment (losing at least 2% at week 8) resulted in an average weight loss of 9.3% within 24 weeks.



Boosts Satiety⁵

Participants steadily lost weight week by week, and showed they felt:

- Fuller sooner
- Less need to snack between meals
- More satisfied with smaller portions

Low Risk⁶

Epitomee® Capsule had zero related serious adverse events. The Epitomee® Capsule is designed to provide effective results without safety risk. Adverse events reported by participants were similar between the Epitomee® and placebo groups, indicating that the Epitomee® treatment has no apparent side effects.

How the Epitomee® Capsule Works

Each Epitomee® Capsule contains a proprietary polymer blend that gently expands in the stomach after being taken with water, forming a soft, gel-like structure. This structure creates a natural* feeling of fullness, signaling the brain to help you feel satiety—so you can eat smaller portions without feeling deprived.



INGESTING

Swallow one capsule with two 8-ounce (250 ml) glasses of water 30 minutes before lunch and dinner.



EXPANDING

Once in your stomach, the device inside the capsule absorbs the water and gently expands into a triangle shape.



SATIATING

The expanded triangle gently presses on your stomach walls, signaling fullness to your brain.



DISINTEGRATING

After several hours, the capsule breaks down completely in the small intestine and is naturally* cleared out of the body.



Tips for Success

- Drink plenty of water and continue healthy eating and movement
- Use as part of your healthy routine
- Be mindful of hunger cues and portion sizes
- Track your progress and stay in touch with your healthcare provider
- Consistency is key—use daily as prescribed