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8.Explanation Of Symbols

ISO 15233-1

Symbols to be used with medical device labels, labelling, and information to be supplied

 REF	Catalog number		Caution
 LOT	Batch code		Temperature limit
	Use-by date		Consult instructions for use
 R_x only	Prescription use only		Do not use if package is damaged

License Holder and Manufacturer:

Epitomee® Medical Ltd.

17 Hatochen St. Caesarea 3079892, Israel.

Telephone number: 074-7400799

Distributor: Rafa Laboratories, Ltd.,
POB 405, Jerusalem 9100301, Israel.

Registration number: 35410002

Patent: www.epitomeemedical.com/patents

7. Glossary

- **Abdomen** – The part of the body between the chest and the pelvis that contains the organs of the digestive system.
- **Abdominal adhesions** – Abnormal bands of fiber that form in the abdomen, and connect organs and tissues, causing them to stick together.
- **Adhesive peritonitis** – Inflammation of the lining of the abdomen caused by inflammation of abdominal organs that are sticking together.
- **Body mass index (BMI)** – A common measure that uses height and weight to determine whether a person's weight is in the healthy range.
 - 18.5-24.9 kg/m²: healthy
 - 25-29.9 kg/m²: overweight
 - 30 kg/m² or higher: obese
- **Crohn's disease** – A disease that causes painful inflammation in the digestive system.
- **Diverticula** – Pouches of tissue from the large intestine that can become inflamed and painful and may cause bleeding.
- **Esophagus** – The tube that allows the passage of food, liquids, and saliva from the mouth to the stomach.
- **Esophageal rings or webs** – Folds that can form in the esophagus and block food from passing through easily to the stomach.
- **Ileus** – Lack of normal movement somewhere in the intestines that leads to a buildup and potential blockage of food material.
- **Intestinal obstruction** – A condition in which food and waste products are blocked from passing through the intestines.

- Epitee® should be taken under the direction of a health care provider, as part of a structured weight loss program. Failure to adhere to prescribed dietary and exercise instructions may result in failure to lose weight.

6.2 Possible Risks Related to Taking Epitee®

In the clinical study performed, side effects with Epitee® were similar in number to placebo (dummy capsule). The most common side effects were infections (the study was performed during the COVID-19 pandemic) abdominal discomfort (pain, cramping), constipation, headache and nausea.

**Potential GI Adverse Events
(rates observed compared to placebo)***

Greater than Placebo	Lower than Placebo	Equal to Placebo	Not observed
• Abdominal discomfort (pain, cramping)	• Constipation • Nausea • Abdominal distension (bloating) • Flatulence • GERD • Vomiting	• All GI adverse events combined (93.5% mild) • Diarrhea	• Negative health effects due to weight loss • Allergic Reaction • Bowel obstruction • Choking • Death • Slowed or stopped movement of food through the GI system • Effects on absorption of medication • Need for emergency surgery

*Rates observed in the RESET pivotal study

On an empty stomach Epitomee® can affect how your body absorbs medications you take at the same time. Please take Epitomee® according to “**Directions for Use**” section.

6.1 Precautions

- Contact a doctor right away if problems occur, if you have a severe or continued adverse event. If you have a severe allergic reaction or severe stomach pain, stop using Epitomee® until you can speak with your attending doctor.
- If you have a medical condition that makes it hard to swallow, you may have trouble swallowing Epitomee®.
- Consult with your attending doctor if you are taking medications causing disturbances in the bowel movements (ileus) (such as opioids and tricyclic antidepressants).
- Use Epitomee® with caution if you have active gastrointestinal conditions such as gastro-esophageal reflux disease (GERD), ulcers or heartburn.
- Avoid taking Epitomee® if you:
 - have problems with your esophagus, including webs and/or rings.
 - have a diverticula.
 - are likely to have strictures (a narrowing of the intestine such as patients with Crohn’s disease may have).
 - have had gastrointestinal surgery with complications (for example obstruction and/or adhesive peritonitis or known abdominal adhesions).
- Do not consume Epitomee® if blister pack or package are damaged.
- Do not take any capsule which is broken, crushed, or damaged. Discard this capsule.
- Epitomee® is NOT a food substitute. It is not absorbed by the body and therefore has no nutritional value.

For each dose, follow the steps below:

1. Ensure that the capsule is not broken, crushed, or damaged (throw away any damaged capsule)
2. Swallow one capsule with two glasses of water (8 fl oz/250 ml each). Taking the capsule with beverages other than water could lead to ineffective treatment.
3. After taking Epitee[®], wait about 30 minutes before starting the meal.

It is recommended to drink water regularly throughout the day, in addition to the water you drink with the capsule.

If a pre-meal dose is missed, take the capsule during that meal according to the directions for use, and swallow one capsule with two glasses of water.

To avoid impact on the absorption of medications:

- For all medications that should be taken with food, take them after starting a meal.
- For all medications that are taken once daily, take them either before taking the Epitee[®] capsule in the morning (fasting or with breakfast) or after taking the Epitee[®] capsule, at bedtime, and as prescribed by your doctor.

6. Warnings and Safety Information

Read this package insert in its entirety before using Epitee[®].

Keep out of reach of children.

Epitee[®] may alter the absorption of medications. Please review "Directions for Use" section carefully.

Do not use Epitee[®] capsules after the expiration date printed on the blister pack or package.

When combined with diet and exercise, Epitomee® can help manage your body weight. In the RESET clinical trial, 56% of the individuals taking Epitomee® for 6 months lost 5% or more of body weight and lost on average, 6.6% of their body weight.

2. Who Can Take Epitomee®?

Epitomee® is indicated to aid weight management in adults (18+) with excess weight or obesity, Body Mass Index (BMI) of 25-40 kg/m², when used in conjunction with diet and exercise.

3. Who Cannot Take Epitomee®?

Do not take **Epitomee®**:

- If you are pregnant.
- have had allergic reaction to one or more of the **Epitomee® components**: cross-linked polyacrylic acid sodium salt, hydroxy propyl methyl cellulose (HPMC), cellulose acetate phthalate, hydroxypropyl cellulose, chitosan (vegan), diethyl phthalate, polyvinyl alcohol polyethylene copolymer (Kollicoat), triethyl citrate, cellulose acetate, acetyl tributyl citrate.

4. Storage

- Epitomee® capsules should be kept closed in the blister pack and stored at room temperature (between 59°F - 77°F / 15°C - 25°C).
- Epitomee® capsules should be kept in their original package until use.

5. Directions for Use

Epitomee® capsule should be taken with water twice a day, 30 minutes before lunch and 30 minutes before dinner. Each dose includes 1 capsule of Epitomee® provided in a blister pack.


1.Introduction

Epitomee[®] is an oral capsule that helps you manage your weight. Take the capsule with water 2 times a day, 30 minutes before lunch and dinner. Epitomee[®] works in the stomach to help you feel full and eat less food throughout the day.

When taken with water, Epitomee[®] releases compartments containing gel particles that form a triangular shape in the stomach (see Figure). The particles absorb water, and the compartments swell - thereby creating a feeling of early satiety leading to reduced eating.

Epitomee[®] passes through the stomach and breaks down in the small intestine into small particles. Then the particles are secreted out of the body in normal bowel movements (they are not absorbed in the body).

Epitomee[®] is supplied in a carton package containing either 1 blister pack or 4 blister packs. The blister pack contains 14 doses, each dose includes one capsule, to be taken with water before lunch and dinner.

Epitomee[®] ingestion and passage through GI tract

epitomee Patient Information

Rx only

Please talk to your healthcare provider if you have any questions or do not understand any of this information.

Content

1. Introduction	2
2. Who Can Take Epitomee [®] ?	3
3. Who Cannot Take Epitomee [®] ?	3
4. Storage	3
5. Directions for Use	3
6. Warnings and Safety Information	4
6.1 Precautions	5
6.2 Possible Risks Related to Taking Epitomee [®]	6
7. Glossary	7
8. Explanation Of Symbols	8

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