



WAVEWRITER ALPHA™ IMAGEREADY MRI FULL BODY PATIENT ELIGIBILITY

This form provides information about the patient's implanted WaveWriter Alpha or WaveWriter Alpha Prime Spinal Cord Stimulator System and MRI scan eligibility. It may be provided to the radiologist to support the confirmation of the patient's scan.

- Prior to performing an MRI Scan, confirm that the patient's stimulator is placed in MRI Mode.
- Refer to www.bostonscientific.com/imageready for labeling and safety conditions

Patient's Name: _____ Date: _____

Physician's Name: _____ Phone: _____

Physician's Signature: _____

A. MR CONDITIONAL WAVEWRITER ALPHA SYSTEM INFORMATION	MODEL #	MRI FULL BODY ELIGIBLE	NOT MRI ELIGIBLE
1. Implantable Pulse Generator (IPG)			
WaveWriter Alpha IPG (32 contact or 16 contact)	SC-1232 or SC-1216	<input type="checkbox"/>	
WaveWriter Alpha Prime IPG (32 contact or 16 contact)	SC-1432 or SC-1416	<input type="checkbox"/>	
Note: if you have another model number IPG, please refer to the labeling specific to your IPG model number.			
2. Percutaneous and/or surgical paddle leads (check all that apply)			
Avista™ MRI Percutaneous Lead, 8-contact lead (56 cm or 74 cm)	SC-2408-56 or SC-2408-74	<input type="checkbox"/>	
Linear™ Percutaneous Leads (50 cm or 70 cm)	SC-2158-50, SC-2138-50 or SC-2158-70, SC-2138-70	<input type="checkbox"/>	
Linear ST Percutaneous Leads (50 cm or 70 cm)	SC-2218-50, SC-2208-50 or SC-2218-70, SC-2208-70	<input type="checkbox"/>	
Linear 3-4 Percutaneous Leads (50 cm or 70 cm)	SC-2352-50 or SC-2352-70	<input type="checkbox"/>	
Linear 3-6 Percutaneous Leads (50 cm or 70 cm)	SC-2366-50 or SC-2366-70	<input type="checkbox"/>	
Infinion™ CX Percutaneous Leads (50 cm or 70 cm)	SC-2317-50 or SC-2317-70	<input type="checkbox"/>	
Artisan™ or Artisan MRI Surgical Leads (50 cm or 70 cm)	SC-8416-50, SC-8416-70, or SC-8216-50, SC-8216-70	<input type="checkbox"/>	
CoverEdge™ 32 or CoverEdge 32 MRI Surgical Leads (50 cm or 70 cm)	SC-8336-50, SC-8336-70, SC-8436-50, SC-8436-70	<input type="checkbox"/>	
CoverEdge X 32 or CoverEdge X 32 MRI Surgical Leads (50 cm or 70 cm)	SC-8352-50, SC-8352-70, SC-8452-50, SC-8452-70	<input type="checkbox"/>	
Other Lead(s):			<input type="checkbox"/>
Extensions or Splitters:			<input type="checkbox"/>
3. Adapters			
Precision™ Adapter M8, 15cm	SC-9218-15	<input type="checkbox"/>	
Precision Adapter M8, 55cm	SC-9218-55	<input type="checkbox"/>	
Precision Adapter S8, 15cm	SC-9208-15	<input type="checkbox"/>	
Precision Adapter S8, 55cm	SC-9208-55	<input type="checkbox"/>	
4. Medtronic Lead connected to Precision Adapter M8			
Vectris™ SureScan™ MRI 1x8 subcompact Lead, 60cm Vectris	977A160	<input type="checkbox"/>	
SureScan MRI 1x8 subcompact Lead, 75cm	977A175	<input type="checkbox"/>	
Vectris SureScan MRI 1x8 subcompact Lead, 90cm	977A190	<input type="checkbox"/>	
Vectris SureScan MRI 1x8 compact Lead (60cm / 75cm / 90cm)	977A260 / 977A275 / 977A290	<input type="checkbox"/>	
Specify™ SureScan MRI 2x8 paddle	977C265, 977C290	<input type="checkbox"/>	
Specify SureScan MRI 5-6-5 paddle	977C165, 977C190	<input type="checkbox"/>	
5. St. Jude (Abbott) Lead connected to Precision Adapter			
Octrode™, 60cm	3186	<input type="checkbox"/>	
Penta™, 60cm	3228	<input type="checkbox"/>	
<p>Note: Leads should be connected directly into the IPG, Patient should not be implanted with lead extensions or splitters. The Precision Adapter M8, 55cm (SC-9218-55) with 75cm or 90cm Medtronic Leads (977A175, 977A275, 977A190, and 977A290) was not tested for MRI conditionality for the WaveWriter Alpha System.</p> <p>When connecting to the Adapter M8 or Adapter S8, use one 16 Contact Lead, one 8 Contact Lead, or two 8 Contact Leads, All Leads must be from the same manufacturer. Do not add Boston Scientific Leads. Having a combination of Leads from different manufacturers was not tested for MRI conditionality for the WaveWriter Alpha System.</p> <p>Do not use Adapters of different lengths (for example, one 15cm and one 55cm) when connecting two of the Adapter M8 or Adapter S8 to the WaveWriter Alpha System. Connecting Adapters of different lengths was not tested for MRI conditionality for this system.</p>			

B. PATIENT IMPLANT CONFIGURATION INFORMATION (ALL QUESTIONS MUST BE ANSWERED)	MRI FULL BODY ELIGIBLE	NOT MRI ELIGIBLE
1. The lead implant location is epidural	Yes	No
2. The IPG is implanted in the upper buttock or lower flank	Yes	No
3. Patient has no abandoned leads or IPGs (i.e. leads or IPGs that are not connected to the functioning WaveWriter Alpha or WaveWriter Alpha Prime System)	Yes	No
4. No evidence can be found of fractured leads or compromised IPG-lead system integrity	Yes	No

C. INSTRUCTIONS FOR THE PATIENT PRIOR TO THE MRI EXAM

1. For patients with a rechargeable IPG (SC-1232 or SC-1216), instruct the patient to fully charge their IPG and bring the Charger to the MRI Center (in case charging is necessary)
2. Instruct the patient to bring their Remote Control to the MRI exam and to enable MRI mode before the MRI Scan
Note: The Charger and Remote Control are MR Unsafe and must not be brought into the MRI Scanner Room.

Boston
Scientific

Advancing science for life™

Europe, Middle East and
Africa (EMEA) Headquarters
2 Rue René Caudron, 78960
Voisins-le-Bretonneux, France
[bostonscientific.com/en-EU](https://www.bostonscientific.com/en-EU)

Boston Scientific – S.A. 8 Anslow
Crescent Anslow Office Park,
Bryanston, Sandton, 2021,
South Africa

© 2024 by Boston Scientific
Corporation or its affiliates. All rights
reserved.

NM-855510-AB

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.