

AKADEMIKA



Product Catalogue



## FEATURES

- This laboratory is designed to study the physics of fiber optics.
- It deals with the study of various important characteristics of optical fiber cable like numerical aperture, V number, attenuation, coupling loss and mode field diameter

## SPECIFICATIONS

## Laser

- Source Type : He-Ne Laser
- Wavelength : 632.8nm
- Light : Red Visible
- O/P Power : 2 mW
- Complete system with power supply

### Laser to fiber coupler

- Laser to fiber coupler with Lens adjustment facility to adjust beam into core of fiber
- Coupling Efficiency : >70% for SM fibers  
> 90% for MM fibers
- Operation Wavelength : 180 to 2000nm
- Power Handling Capacity : > 1watt

### Laser power meter

- Separate Sensor Unit with Stand
- Separate Display Unit
- Power Measurement :  $2\mu\text{W} \sim 20\text{W}$
- Resolution :  $1\text{nW}$
- Power Density :  $15\text{KW}/\text{cm}$
- Sensor Diameter :  $20\text{mm}$
- Operation Wavelength :  $400 \sim 1100\text{nm}$
- Calibrated to  $633\text{nm}$

### Optical Patchcord ST-ST


- Single Mode : 9/125 um 1mtr.
- Multi Mode : 62.5/ 125 um 1 mtr.
- Multi Mode Fiber Pigtail : 62.5/ 125um 3 mtr.
- Multi Mode : 100/140um 1mtr.

## Spool

- 3 in 1 Multi Mode Spool(100m,500m,1000m)

## MEASUREMENT OF COUPLING LOSS


- In order to achieve maximum coupling efficiency, the fiber core diameter has to be bigger than the spot diameter
- Another factor which determines the coupling efficiency is the matching of fiber numerical aperture to numerical aperture of the focused rays.
- The size of the single mode fiber is very small, around 4-9 microns.



- For minimum coupling loss couple alignment should be proper and LASER beam characteristics should be matched properly

## MEASUREMENT OF NUMERICAL APERTURE

- Numerical Aperture is range of acceptance angles rotated around axis, where light can be coupled into a waveguide such as the core of an optical fiber.




### NUMERICAL APERTURE

- The blue cone is known as the cone of acceptance
- Light waves within the acceptance cone can be collected in a small area which can then be sent into the optical fiber
- Numerical aperture is measure of light gathering capacity of optical fiber Light outside the acceptance cone cannot be coupled to onto optical fiber
- It can be calculated quantitatively using formula

$$NA = n \sin \theta$$

### MEASUREMENT OF ATTENUATION OF A OPTICAL FIBER BY CUTBACK METHOD:

- The cutback method is used for measuring the total attenuation of an optical fibre.
- The cutback method begins by measuring the output power  $P_y$  of the test fibre of known length 'L'. Without disturbing the input conditions, the test fibre is cut back to a shorter length






- The output power  $P_x$  of the short test fiber is then measured and the fibre attenuation 'A' and the attenuation coefficient 'a' are calculated


### MEASUREMENT OF MODE FIELD DIAMETER OF AN OPTICAL FIBER

- The mode field diameter (MFD) represents a "measure of the transverse extent of the electromagnetic field intensity of a mode of light in a fiber cross section".
- In optical fiber, this typically is larger than the fiber core, since a portion of the light propagates through the cladding.
- MFD is determined using a Gaussian approximation of the intensity distribution with the MFD defined as the width of the curve at the  $1/e^2$  power level.
- MFD is a function of wavelength, the core diameter, and the refractive-index.




### MODE OBSERVATION

- Light can only be launched at certain discrete angles within the N.A. cone which can propagate through fiber core successfully
- Acceptance cone can not be assumed as a solid cone of rays, must be viewed as composed of discrete annular conical rings of rays.
- Discrete values of launching angles lead to formation modes



- In optical fiber, modes are different patterns of light intensity distribution around the core axis



NUMERICAL APERTURE

### ACCESSORIES

- Optical Bread Board : 01 No
- X-Y-Z Translation Stage : 01 No
- Gimbal Mount for He-Ne Laser : 01 No
- Rotary Stage : 01 No
- Laser Fiber Coupler : 01 No
- He-Ne Laser Source With Power Supply : 01 No
- Laser Power Meter : 01 No
- Laser to Fiber Coupler : 01 No
- Coupler Holder : 01 No
- Coupler Holder Base : 01 No
- Rotary Stage Block : 01 No
- Rotary Fiber Holder With Receptacle : 01 No
- Photo Detector Holder : 01 No
- Photo Detector Cap With Receptacle : 01 No
- Photo Detector Cap With Knife-edge Slit : 01 No
- Fiber Holder Stand : 01 No
- Fiber Holder Stand Base : 01 No
- Fiber Holder Stand Horizontal Plate : 01 No
- Fiber Holder Stand Vertical Plate : 01 No
- Photo Detector Holder Adjustable Stand : 01 No
- Gimbal Mount Block : 01 No
- Fiber Coupler Block : 01 No
- Front Plate : 01 No

### Fiber Cable

- 9/125 micron 1 Meter SMPatch Cord : 01 No
- 62.5/125 micron MMPatch Cord : 02 No
- 100/140 micron 1 meter MM Pigtail : 01 No
- 3 in 1 Multi Mode Spool (100,500,1000) : 01 No
- Allen Key Set
- Allen Bolt
- Acrylic Plate 4.5" X 4.5" X 6 mm : 01 No
- Power Cord : 01 No

# Our Products



Basic Electronics



Analog Communication



Digital Communication



RF/Antenna & Microwave



Fiber Optics



Radar



Wireless Communication



Network Laboratory



Processor



Controls & Instrumentation



Test & Measuring Instruments



Drone Technology Laboratory

**AKADEMIKA**

Unit No 128/129, Hema Industrial Estate  
Sarvodaya Nagar, Jogeshwari ( E )  
Mumbai – 400060

+91 9004904462

[www.akademika.in](http://www.akademika.in)

[info@akademika.in](mailto:info@akademika.in)

• DISTRIBUTOR •

