
P2 Nasal Pillows Interface User Manual


CE 0123

COMPONENTS OF THE MASK


- | | |
|------------|--|
| 1 Headgear | 2 Nasal Pillows (available in S, M, L) |
| 3 Panel | 4 Soft sleeves |
| 5 Clip | 6 Elbow |
| 7 Frame | 8 Short tube |
| 9 Swivel | |


Fitting

Disassembly

Reassembly

P2 Nasal Pillows Interface

Thank you for choosing BMC's P2 Nasal Pillows Interface. The nasal pillows interface is light, flexible and quiet. It is designed to minimize contact with your face, thus ensuring that you feel comfortable during therapy. This user manual provides you with the information you need for the correct use of your mask. The P2 Nasal Pillows Interface and accessories are not made with natural rubber latex or DEHP.

Intended use

The P2 Mask is intended to provide an interface for patients who have been prescribed CPAP or bi-level therapy.

- To be used by adult patients (>66 lb / 30 kg) for whom positive airway pressure therapy has been prescribed.
- Intended for single-patient re-use in the home environment and multi-patient re-use in the hospital/institutional environment.

CAUTION: In the US, Federal law restricts this device to sale by or on the order of a physician.

Contraindications

- This mask should be used with caution for patients with following conditions:
- Facial deformity.
 - Injured facial skin which will contact with the mask during use.

Before using the mask

⚠️ WARNINGS

- The vent holes must be kept undamaged and free of blockages.
- The mask is not suitable for patients who need ventilation to support life.
- This mask should only be used with CPAP or bi-level devices recommended by a physician or respiratory therapist.
- The mask should not be used unless the device is turned on. Once the mask is fitted, ensure the device is blowing air.
- Explanation of the warning:** CPAP and bi-level devices are intended to be used with special masks (or connectors) which have vent holes to allow continuous flow of air out of the mask. When the device is turned on and operating properly, the fresh air from the device flushes the exhaled air out through the attached mask exhalation port. However, when the device is not operating, the fresh air supplied through the mask is insufficient, and exhaled air may be re-inhaled. Rebreathing of exhaled air for more than several minutes may, in some circumstances, lead to suffocation. This applies to most models of CPAP or bi-level devices.

- At low CPAP pressures the flow through the exhaled port may be inadequate to clear all exhaled gas from the tube. Some rebreathing may occur.
- To minimize the risk of vomiting during sleep, the patient should avoid eating or drinking three hours before using the mask. This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Not applicable for silicone-allergy user.

- This mask should not be used on patients who are uncooperative, unconscious, unresponsive, or unable to remove the mask.
- Using the mask may cause pain in the teeth, gums or jaw or aggravate existing dental conditions. If these conditions occur, consult your health care professional.
- When the patient is unable to remove the mask himself, the mask must be used under qualified supervision.
- Follow all precautions when using supplemental oxygen.
- Oxygen flow must be turned off when the CPAP or bi-level device is not operating, so that unused oxygen does not accumulate within the device enclosure and create a risk of fire.

- Oxygen supports combustion. Oxygen must not be used when smoking or in the presence of an open flame. Only use oxygen in a well-ventilated room.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration varies, depending on the pressure settings, patient breathing pattern, mask, point of application and leak rate.

- The technical specifications of the mask are provided for your clinician to check if it is compatible with the CPAP or bi-level device. If used outside specification or if used with incompatible devices, the seal and comfort of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect the CPAP or bi-level device function.
- Stop using the mask and consult your physician or sleep therapist, if you have ANY adverse reaction to the use of the mask.
- Refer to your CPAP or bi-level device manual for details on settings and operational information.
- Remove all packaging before using the mask.
- Images shown here are indicative only. If there is inconsistency between the image and actual product, the actual product shall govern.
- Any serious incidents relating to the product should be reported to BMC and to the competent authorities in your country.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents, other than those specified in the accompanying documents, or exceeding the number of processing cycles can have an adverse effect on the mask and consequently the safety or the quality of the therapy.

Fitting

- Thread the two side straps of the headgear through the gaps in panels separately, then fold back over to secure.
- Put the headgear over your head by holding the back strap and stretching it around the back of your head. The top strap should sit comfortably on top of your head.
- Adjust each end of the side strap evenly until the mask sits comfortably just under the nose. Reattach the fastening tabs.
- Move the mask away from your face then place it back to get a comfortable fit. Ensure the nasal pillows fit securely into your nostrils with both pillows positioned upright.
- Finally connect the end of the mask tubing into your device air tubing. Your mask is now fitted and ready to use.

Notes:

- Refer to your CPAP or bi-level device manual for details on setting and operational information.
- Three nasal pillows sizes (S, M, L) are available, hold the pillows on your nostrils ensuring that the pillows size (e.g. M) correctly fits for you.
- If you feel air leak, readjust the angle of the nasal pillows to get a better seal. If you are unable to resolve any leaks, you may need an alternative nasal pillows size.

Removal

To remove the P2, pull the nasal pillows assembly away from your nose and lift the headgear off your head. You do not need to loosen the headgear straps.

Disassembly

- Holding the short tube on the mask, squeeze and pull the nasal pillows away from the clip in the mask frame.
- Remove the headgear from the mask.

Notes:

- The short tube assembly cannot be taken apart.
- Do not pull on the short tube or the device air tubing as they may get damaged.

Reassembly

- Attach the nasal pillows to the clip in the mask frame. Make sure the Nasal Pillows triangular opening and the frame platform opening is to point down.
- The Nasal Pillows is in the correct position. Make sure the frame and Nasal Pillows will fit tightly together without any gaps.
- Thread the two side straps of the headgear through the gaps in panels separately, then fold back over to secure.

Cleaning the mask at home

Notes:

- Wash your hands before cleaning the mask.
- Wash the headgear before first time use. The headgear may fade slightly that is normal. It is important to always handwash the headgear.

⚠️ CAUTIONS

- Do not use solutions containing bleach, chlorine, alcohol, aromatics, moisturizers, antibacterial agents, or scented oils to clean any part of the mask. These solutions may cause damage and reduce the life of the product.
- Exposing any part of the mask to direct sunlight or heat may cause deterioration.
- Visual criteria for product inspection: If there is any visible deterioration of a mask component (cracking, crazing, tears, etc.), the component should be discarded and replaced.
- Replacement of parts must be carried out by professional service personnel.

Preparation before cleaning

- Prepare tools: soft bristle brush, drinking quality water (10 litres per mask), mild liquid detergent (e.g. Alconox diluted to 1%).
- Disassemble the mask components according to the disassembly instructions.

Daily / After each use

- Soak the components in a solution of mild liquid detergent (e.g. Alconox diluted to 1%) with drinking quality water (25°C to 35°C) for 10 minutes.
- While soaking, clean each component with a soft bristle brush for 1 minute. Pay special attention to all crevices and cavities.
- Agitate the components vigorously twice in drinking quality water. 5 litres of water per mask at a time.
- Allow the components to air dry and do not allow direct sunlight.
- When all components are dry, reassemble the mask according to the reassembly instructions.

Note:

- Cleaned once a day or after each use, the mask can be washed up to 365 times.

Weekly

Handwash the headgear and its components in a solution of mild liquid detergent (e.g. Alconox diluted to 1%) with drinking quality water (25°C to 35°C) for 10 minutes. Rinse the headgear and its components well with drinking quality water and allow them to air dry out of direct sunlight.

Note:

- Cleaned once a week, the headgear can be washed up to 52 times.

Reprocessing the mask between patients

This mask should be reprocessed when used between patients. Cleaning and disinfection instructions are available from the BMC website, www.bmc-medical.com. If you do not have internet access, please contact your BMC representative.

Note:

- The maximum number of disinfection cycles is 20.

Intended clinical benefits

- Improve patients' quality of sleep.
- Improve patients' compliance with therapy.

Storage

Ensure that the mask is thoroughly clean before storing it for any length of time. Store the mask in a dry place out of direct sunlight.


Disposal

The mask does not contain any hazardous substances and may be disposed of with your normal household waste. For patients with transmissible diseases, the product shall be discarded as medical waste, otherwise it may cause the spread of disease.

Service life

The service life of the P2 Nasal Pillows Interface depends on the intensity of usage, maintenance, and environmental conditions to which the mask is used or stored. As this mask system and its components are modular in nature, it is recommended that the user maintain and inspect it on a regular basis, and replace the mask system or any components if deemed necessary or according to the "Visual criteria for product inspection" in the "Cleaning the mask at home" section of this manual. The mask should be replaced within 12 months from the date of first use.

Technical specifications

Pressure (cmH ₂ O)	4	8	12	16	20
Flow rate (L/min)	20	29	38	45	50

Dead space information	Dead space is the empty volume of the nasal pillows and short tube assembly up to the swivel. The dead space of the P2 for Large pillows is 92 mL.	
Therapy pressure	4 to 20 cmH ₂ O	
Resistance	Drop in pressure measured (average for three nasal pillows sizes) at 50 L/min: 0.5 cmH ₂ O; at 100 L/min: 1.6 cmH ₂ O	
Sound	DECLARED DUAL-NUMBER NOISE EMISSION VALUES in accordance with ISO 4871. The A-weighted sound power level of the mask is 26 dBA, with an uncertainty of 3 dBA. The A-weighted sound pressure level of the mask at a distance of 1 m is 18 dBA, with an uncertainty of 3 dBA.	
Environmental conditions	Operating temperature: +5°C to +40°C (41°F to 104°F) Operating humidity: 0% to 93% relative humidity non-condensing Storage and transport temperature: -20°C to +55°C (-4°F to 131°F) Storage and transport humidity: 0% to 93% relative humidity non-condensing	
Date of manufacture	See product label	

Symbols

	Caution		Use-by date
	Authorized Representative in the European Community		Not made with natural rubber latex
	Batch code		Importer
	Catalogue number		Model number
	Temperature limit		Made in China
	Humidity limitation		Medical device
	Manufacturer		Unique device identifier

Limited warranty

It is warranted that P2 Nasal Pillows Interface, including its components, frame, nasal pillows, headgear, short tube and other accessories shall be free from defects in material and workmanship for a period of ninety (90) days from the date of purchase by the initial consumer. To exercise the rights under this warranty, contact the local authorized dealers.

MANUFACTURER:
BMC MEDICAL CO., LTD.
 Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA
 URL: en.bmc-medical.com
 E-mail: intl@bmc-medical.com
 Tel: +86-10-51663880

EU AUTHORISED REPRESENTATIVE:
Shanghai International Holding Corp. GmbH (Europe)
 Effestraße 80, 20537 Hamburg, Germany
 Tel: 0049-40-2513175
 Fax: 0049-40-255726

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