




# Ion-Exchange Chromatography



Charge-based resolution from capture to polishing, covering cation and anion exchange, strong and weak functional groups, and bind/elute versus flow-through modes.

# 01 What IEX chromatography does

Ion exchange chromatography (IEX) separates molecules according to differences in surface charge under defined buffer conditions. The stationary phase carries fixed ionic functional groups, either negatively charged (cation exchanger) or positively charged (anion exchanger), each associated with mobile counterions. Biomolecules carrying the appropriate opposite charge displace these counterions and bind to the resin; those with insufficient or opposite charge pass through.

For proteins, the operative concepts are pH and isoelectric point. A protein carries a net positive charge at pH values below its isoelectric point (pI) and a net negative charge at pH values above its pI. The relationship between operating pH and pI determines which exchanger type is appropriate. The general rule is: the exchanger is named for what it exchanges, not for what it retains. A cation exchanger retains cations.


| Exchanger        | Properties                                                                     |
|------------------|--------------------------------------------------------------------------------|
| Cation exchanger | carries a negative charge and retains positively charged (cationic) molecules. |
| Anion exchange   | carries a positive charge and retains negatively charged (anionic) molecules.  |

Net charge at a given pH gives the first approximation of IEX behaviour, but it is only an approximation. Actual chromatographic retention depends on the spatial distribution of charges across the protein surface. Proteins with identical pI values can show different retention on the same column because their surface charge patterns differ. Buffer species, ionic strength, resin chemistry, and pore geometry all modulate the interaction further. Method development must therefore be empirical as well as principle-guided.

IEX is applicable across a wide range of targets: proteins, peptides, nucleic acids, and other charged biomolecules. It can function as a capture step, an intermediate purification step, a polishing step, or an analytical tool for charge-variant characterisation.

# 02 The separation process

The generalized IEX process is:



Equilibration associates the resin functional groups with their mobile counterions under the binding buffer conditions. When the sample is loaded, target molecules with sufficiently strong electrostatic complementarity to the resin displace the counterions and are retained. Molecules with insufficient charge, the wrong charge sign, or strong competition from high-ionic-strength conditions pass through in the unbound fraction.

Binding requires low ionic strength. High concentrations of salt ions compete with the target molecule for electrostatic interaction with the resin surface, suppressing binding. Conductivity is therefore the operational lever for controlling loading conditions: samples are adjusted to a sufficiently low conductivity before or during loading, and purity is improved by washing at the same or slightly elevated conductivity before elution.

Elution is most commonly achieved with a salt gradient. Increasing counterion concentration progressively displaces bound species in order of their electrostatic affinity for the resin: weakly retained components elute first, and more strongly retained species elute later. A well-resolved salt gradient can separate proteins with small differences in charge; a step elution defines discrete pools for process applications..

pH-gradient elution is an alternative. For a cationic protein bound to a cation exchanger, increasing pH toward and above its pI reduces net positive charge and weakens binding. For an anionic protein bound to an anion exchanger, decreasing pH toward the pI reduces net negative charge. pH elution can provide different selectivity than salt elution and may be advantageous when salt sensitivity is a concern, but it requires careful buffer design to maintain a defined, smooth pH gradient.

The distinction between strong and weak ion exchangers affects the pH dependence of the functional group. Strong exchangers, such as sulfonate (SP) for cation exchange and quaternary ammonium (Q) for anion exchange, retain their charge across a broad pH range. Weak exchangers, such as carboxymethyl (CM) and diethylaminoethyl (DEAE), have ionization states that are more sensitive to pH. Neither type is universally superior; the choice is a selectivity and process-design decision made in the context of the specific application.

## 03 Applications and method development

Ion exchange is the most broadly applied chromatographic mode in protein purification. Its versatility spans every stage of the purification workflow. At capture, IEX can concentrate and partially purify protein from clarified cell culture or lysate. At intermediate stages, it separates product from process-related impurities that survived an earlier capture operation. At polishing, it resolves product-related variants, including charge isoforms, deamidated species, and glycoforms, at the resolution level required for biopharmaceutical characterization and release.

Anion exchange in flow-through mode warrants particular attention as a bioprocess polishing strategy. Under appropriately chosen pH and conductivity conditions, the product of interest passes through the column while negatively charged impurities, including host-cell DNA, endotoxins, host-cell proteins, and some viral particles, bind to the positively charged resin. This configuration can deliver high impurity clearance at high throughput, since the product does not need to be retained and eluted.

Method development should systematically evaluate:

- Target pI and expected net charge at candidate operating pH values.
- Target stability across the required pH range.
- Choice of cation versus anion exchanger, based on the pH window available and expected impurity charge profiles.
- Strong versus weak exchanger, based on the required pH operating range and selectivity requirements.
- Loading conductivity and ionic strength.
- Buffer species and concentration: the ionisation state of the buffer affects background conductivity and can interact with resin selectivity.
- Salt-gradient versus pH-gradient elution.
- Dynamic binding capacity at the intended residence time.
- Target recovery and impurity clearance at each candidate condition.

A practical starting point is to select a pH that places the target sufficiently far from its pI to generate productive binding at accessible ionic strength, while remaining within the protein's stability range. Small pH and conductivity screens then identify the binding window. Once binding is confirmed, gradient scouting reveals whether the target and principal impurities have sufficiently different retention to support a step elution in a process setting.

Because IEX responds to both pH and ionic strength, buffer preparation and pH control are critical determinants of reproducibility. Variation in conductivity at loading or in gradient slope at elution can shift retention times and affect purity and yield. Process robustness assessment should include the sensitivity of the separation to these variables.

## 04 Strengths, limitations, applications

### Strengths

**High resolving power** — for proteins and other charged biomolecules, including closely related charge variants.

**Applicable at every stage** — from capture to analytical characterisation.

**Variable selectivity** — strong and weak exchanger formats provide different selectivity and pH operating behaviour.

**Scalable** — from analytical to large-scale preparative chromatography with straightforward linear scaling.

**Compatible with sensitive proteins** — salt-based elution and flow-through anion exchange both operate under mild aqueous conditions

### Limitations

**Ionic only** — non-ionic molecules are not retained by electrostatic interaction and cannot be separated by this mechanism.

**pH selection** — charge, and therefore retention, changes significantly with pH

**High conductivity** — suppresses binding and must be controlled at loading.

**Tolerance** — of the selected pH and ionic strength by the target must be ensured throughout the process


**pI alone does not predict retention** — surface charge distribution and resin chemistry both contribute.

**Buffer conductivity** — must be measured and controlled for reproducibility.

## Best-fit applications

IEX is the appropriate first-choice screening mode when the target or impurities differ in charge and a suitable pH operating window exists. It is especially valuable for protein capture and polishing, charge-variant analysis, and the removal of charged process impurities including DNA, endotoxin, and host-cell proteins. Its broad applicability makes it the foundational chromatographic mode against which other separation strategies are often benchmarked.

Product-specific operating ranges, capacities, cleaning conditions, and regulatory claims should be taken from current LT Biotech resin documentation and product-specific validation data.



# LT Biotech

UAB LT Biotech

Phone: +370 5 2160227

E-mail: [info@ltbiotech.lt](mailto:info@ltbiotech.lt)

Address: Mokslininku 2, Vilnius, Lithuania, LT-08412

Resource Series

No. 01