

Hydrophobic Interaction Chromatography



How to use salt concentration to turn surface hydrophobicity into a separation handle for proteins and variants.

01 What HIC separates

Hydrophobic interaction chromatography (HIC) separates proteins according to differences in their exposed hydrophobic surface area. All proteins present both hydrophilic and hydrophobic regions; the balance between them, and the degree to which hydrophobic patches are accessible to the surrounding solvent, determines how a protein behaves on an HIC column under a given set of conditions.


Because this selectivity mechanism operates independently of net molecular charge, HIC can resolve proteins and protein variants that appear almost identical by ion exchange chromatography. Two molecules with the same pI and similar size but different surface hydrophobicity can be cleanly separated by HIC. This complementarity makes HIC a powerful addition to any multi-step purification process.

HIC is distinct from reversed-phase chromatography and must not be conflated with it. Both modes exploit hydrophobic interactions, but they do so under fundamentally different conditions and with different consequences for the biomolecule. Reversed-phase chromatography uses organic solvents and typically acidic mobile phases that can unfold proteins. HIC operates entirely in aqueous buffers, uses salt to modulate binding rather than organic solvent to drive elution, and is routinely applied to proteins for which preservation of native structure and biological activity is essential.

The stationary phase consists of porous hydrophilic beads carrying hydrophobic ligands of defined chemistry. Common ligands include butyl, phenyl, and octyl groups, each presenting a different character and density of hydrophobic surface. Ligand choice affects selectivity significantly: a phenyl group offers mild hydrophobicity with some aromatic character, while a longer alkyl chain such as octyl provides stronger retention. Higher ligand density increases retention but can also increase nonspecific binding. Matrix, ligand, and operating conditions must therefore be selected together.

02 The separation process

The process follows a salt-in, salt-out logic that is the opposite of what many practitioners first expect:



Binding is promoted by high ionic strength. Kosmotropic salts, which enhance the ordering of water molecules around solutes, are particularly effective. Ammonium sulfate is the classical choice; sodium sulfate, sodium citrate, and potassium phosphate are alternatives. At high salt concentration, water molecules are preferentially recruited around the dissolved ions, reducing the hydration shell available to hydrophobic protein surfaces. This thermodynamically destabilizes exposed hydrophobic patches in the aqueous environment and drives their association with the hydrophobic ligands on the stationary phase.

Elution is achieved by reducing salt concentration. As ionic strength decreases, the thermodynamic penalty for hydrophobic exposure in solution is reduced, and proteins progressively desorb. Proteins with weaker or less exposed hydrophobic regions leave the column first; those with more extensive or deeply hydrophobic patches require a greater reduction in salt or additional elution conditions. A gradient provides higher resolution between closely related species; a step elution provides a more robust, defined pool for process applications.

The exact response is not simply a linear function of salt concentration. Salt identity, protein conformation, ligand chemistry, and temperature all influence the interaction. Temperature is worth noting: hydrophobic interactions generally strengthen with increasing temperature, which can be exploited or must be controlled. This complexity is the reason results should be interpreted as chromatographic selectivity under the tested conditions, not as an absolute ranking of molecular hydrophobicity.

HIC integrates naturally after upstream operations that leave the sample at elevated ionic strength. An ammonium sulfate precipitation, an IEX elution at high conductivity, or a similar step may produce a sample that requires little additional conditioning before HIC loading, a practical workflow advantage worth planning for.

03 Applications and method development

HIC is used across the purification workflow: as a capture step when appropriate binding capacity exists, as an intermediate purification step when it complements a preceding affinity or IEX step, and as a polishing step when hydrophobic differences between product and impurity are sufficient for resolution. It is particularly valuable for resolving protein aggregates from monomers, for separating product-related variants with different hydrophobic exposure, and for reducing host-cell protein populations that survive affinity or charge-based capture.

Method development typically begins with a small-scale resin screen across multiple ligand chemistries and initial salt conditions. Key variables to evaluate include:

- Salt identity and starting concentration.
- pH, which affects protein conformation and charge and can alter hydrophobic exposure.
- Ligand chemistry and density.
- Gradient slope versus step-elution conditions.
- Residence time and resin loading.
- Temperature, particularly for sensitive molecules or when large-scale heat management is relevant.

A gradient scouting experiment at several salt concentrations is the most efficient way to locate the binding and elution window before committing to step development. Once the approximate conditions are known, the step height, wash stringency, and loading limit can be systematically optimised. The starting salt concentration must be high enough to drive productive binding without precipitating the protein or causing aggregation. For proteins prone to aggregation at high ionic strength, a lower-hydrophobicity ligand or a reduced starting salt concentration may be required. Recovery and product quality, including aggregate content, activity, and downstream compatibility, should be assessed alongside yield at each stage of development.

04 Strengths, limitations, applications

Strengths

Selectivity — based on surface hydrophobicity is orthogonal to charge-based IEX, making the modes complementary in some processes.

Mild conditions — operates entirely under aqueous, non-denaturing conditions.

Versatile — useful across capture, intermediate purification, and polishing stages.

Integrates efficiently — after upstream steps that elevate sample ionic strength.

Multiple ligand chemistries — provide a tuneable range of hydrophobic retention. Can resolve product-related variants, aggregates, and impurities that charge-based modes do not distinguish.

Limitations

Precipitation — high salt loading conditions can cause precipitation for proteins with limited solubility at elevated ionic strength.

Method transfer — salt identity and concentration interact strongly with selectivity. Method transfer between different salt systems requires re-evaluation.


Irreversible binding — very hydrophobic proteins or aggregates may bind irreversibly under standard conditions.

Resin performance — depends jointly on ligand chemistry, density, pore geometry, and operating conditions; no single variable predicts outcome independently.

Best-fit applications

Choose HIC when hydrophobic surface differences are expected to provide useful selectivity and when a fully aqueous, non-denaturing process is required. It is a strong complementary mode after affinity or IEX capture and is particularly valuable for resolving aggregates, product variants, and impurities that remain similar in charge. For applications where organic solvents are acceptable and the highest possible resolution is needed, reversed-phase chromatography may be the more appropriate choice.

Product-specific operating ranges, capacities, cleaning conditions, and regulatory claims should be taken from current LT Biotech resin documentation and product-specific validation data.



LT Biotech

UAB LT Biotech

Phone: +370 5 2160227

E-mail: info@ltbiotech.lt

Address: Mokslininku 2, Vilnius, Lithuania, LT-08412

Resource Series

No. 01