Instructions for Use

Scope
Non sterile Contrast Hand Controller (PN 12101030) for sterilization and subsequent sterile single use.

Intended Purpose
The Contrast Hand Controller is used for real-time and variable-flow of the contrast injection rate for precise and consistent contrast administration on a Cvi@ Contrast Delivery System during angiography or diagnostic and interventional PTCA (percutaneous transluminal coronary angioplasty) procedures.

Technical description
The Hand Controller activates dispense valves on a Cvi@ Contrast System. The Hand Controller has two buttons, one for contrast marked 'C' and one for saline marked 'S'. Each button controls a pressure operated valve on the pump unit. The Hand Controller works with air only, no fluid is present. Dispense volumes of saline and contrast are programmed in the pump operating screen. Each time the button is pressed the programmed volume of saline or contrast is delivered.

Sterilization Instructions
To maintain the sterility the products must be packaged in a primary package according to EN 868 / EN ISO 11607. It is the responsibility of the processor that packaging is carried out under suitable conditions concerning cleanroom. Therefore process validation and routine process monitoring are necessary.

EO sterilization parameters to be used:
EO concentration: 550 - 700 mg/l
Relative Humidity: 30 - 85 % RH
Temperature: 49 - 55 °C
Air removal: -70 kPa
Conditioning: 30±10 min
EO Injection: 40 ml, maximum pressure 35 - 45 kPa
EO Exposure Time: 360 min
EO Removal: -4 kPa
Flushing: 4 air change cycles
Aeration: 7 - 9 h at 30 - 58 °C

Instructions for Use

These sterilization instructions above are capable for preparing the sterile product. It remains the responsibility of the processor to ensure that the sterilization under the processor-specific conditions, as actually performed, using equipment, materials, sterilizer load, and product density in the processing facility, achieves the desired result. This requires validation and routine monitoring of the process.

Disclaimer:
1. For use by physicians and healthcare professionals only.
2. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Setup Instructions:
1. Connect the Hand Controller to the Cvi@ panel, ensuring correct connection of the Male and Female Luer.
2. Complete the fluid set-up as per the manufacturer's instructions for the machine.
3. After completing the fluid set-up process press the OK button on the Cvi@ panel to begin the calibration of the handset.
4. Press and hold the (C) button on the Hand Controller, the machine will display a number on the main screen and finish the step. This process calibrates the user's push to the handset.

Note: A calibration value of 10 or less will result in a warning message on the Cvi@ panel that the calibration command has failed.

Instructions during Use:
1. The harder you push the faster the contrast medium will flow.
2. The actual flow rates seen during injection will be displayed on the Cvi@ panel.
Note: The (S) button for saline is on ON/OFF action.
1. Press once to begin saline flush.
2. Press once to stop saline flush.
3. If a stop action push is not seen the saline flush will time out after 10 seconds.

Warnings
1. Read any instructions before use.
2. Do not use if the sterile packaging is damaged.
3. The product is intended for SINGLE USE ONLY.
4. Re-use of a single-use device creates a potential risk to the patient or user. This may lead to a contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
5. The product should not be used for anything other than the intended use.

Storage
Store in the original packaging under clean and dry conditions.

Precision
1. All connections should be finger tight.
2. Overtightening may lead to cracks and leaks. There are no safety-related aspects to be considered during disposal.

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