1. Symbols / Symbole / Symboles

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Warnings

1. Read any instructions before use.
2. Do not use if the sterile packaging is damaged.
3. The product is intended for SINGLE USE ONLY.
4. Reuse of a single-use device creates a potential risk to the patient or user. This may lead to a contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
5. The product should not be used for anything other than the intended use.
6. The product is not to be used in connection with an active device (e.g. pump).

Precaution

1. Use proper aseptic techniques when handling this device.
2. The presence of air in the system may lead to embolism and in rare instances to death.
3. There are no safety-related aspects to be considered during disposal.

Storage

Store in the original packaging, under clean and dry conditions.

Instructions for Use / Gebrauchsanweisung
Instructions d'Utilisation
Y-PTCA Connector / Y-PTCA Konnektor
Connecteur Y-PTCA

MEDNET

Instructions for Use / Gebrauchsanweisung
Instructions d'Utilisation
Y-PTCA Connector / Y-PTCA Konnektor
Connecteur Y-PTCA

CE

EN

Instructions for Use / Gebrauchsanweisung
Instructions d'Utilisation
Y-PTCA Connector / Y-PTCA Konnektor
Connecteur Y-PTCA

CE

REV 04 - 2022 - 06

EN / DE / FR
Part 1 of 4

Scope
Non sterile Y-PTCA connector (PN 1 020X00X) for sterilization and subsequent sterile single use.

Intended Use
The Y-PTCA connector serves as a valve to avoid blood loss during the introduction of a guide wire or catheter into the vascular system during a PTA-/PTCA-procedure.

Technical Description
PTCA valve: Tuohy Bors, Push Pull or Push Click Bors, 9 French (OD 3.00 mm) for Push Click Y-PTCA connector
Maximum size of the catheter: 8 French (OD 2.67 mm) for Push Pull, 8.5 French (OD 2.83 mm) for Tuohy Bors, 9 French (OD 3.00 mm) for Push Click Y-PTCA connector

Sterilization Instructions
To maintain the sterility the products must be packaged in a primary package according to EN 868 / EN ISO 11607. It is the responsibility of the processor that packaging is carried out under suitable conditions concerning cleanroom.
Therefore process validation and routine process monitoring are necessary.
The following EO sterilization process parameters have been validated:
EO concentration: 800 mg/l
Relative Humidity: 30 - 90 % RH
Temperature: 48 - 58 °C
Air removal: max. vacuum -60 kPa
Conditioning: 10 min
EO Injection: ≥ 30 min, max. pressure 40 kPa
EO Exposure Time: 480 min
EO Removal: -30 kPa
Fushing: 5 air change cycles
Aeration: 168 h at room temperature

The sterilization instructions above are capable of preparing the sterile product. It remains the responsibility of the processor to ensure that the sterilization under the processor-specific conditions as actually performed, using equipment, materials, sterilizer food, and product density in the processing facility, achieves the desired result. This requires validation and routine monitoring of the process.

Instructions for Use:

1. The product may only be used in a sterile condition.
2. Connect the side port of the Y-Connector to fluid source. Then flush and fill the assembly with the fluid. To fill the valve section, open the valve, place one finger over the opening and continue to fill the assembly.
3. Connect the Y-Connector to the guiding catheter. Aspirate the valve to remove any trapped air and flush thoroughly with fluid. Purge any blood by opening the valve whilst continuing to flush the assembly. Close the valve once the blood has been purged. Inspect carefully for any air bubbles and reflux if necessary.
4. Open the valve and insert the guide wire or catheter. (If desired, insert the guidewire introducer and advance it through the valve). Advance the guide wire or catheter to the required distance into the vasculature. Close the valve around the shaft of the guide wire or catheter.

Note: It is important for the Y-Connector Valve to be closed tightly to prevent blood leakage but not too tight to restrict function.

5. Refer to the dilatation catheter labelling for intended use, contra-indications and potential complications associated with the use of dilatation catheters.
6. Withdraw the catheter to the point that the distal tip remains 10 to 20 cm inside the vasculature. Open the Y-Connector Valve and withdraw the catheter completely. Disconnect the Y-Connector Valve from the guiding catheter.

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