EN

MECNET

Instruction for Use / Gebrauchsanweisung / Instructions d'Utilisation
Fluid Extension Lines / Flüssigkeitsmanagement-Verlängerungen /
Lignes de prolongement pour fluides

CE

0482

1. Symbols / Symbole / Symboles

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0482

REV 01 - 2025 - 11

Scope
Non sterile Fluid Extension Lines (PN 1051 XXXXXX, 1052 XXXXXX, 1054 XXXXXX) for sterilization and subsequent sterile single use.

Intended Purpose
The Fluid Extension Lines are intended as a kind of multifunction connection tube and can be widely used for patients who require the administration or infusion of liquids/contrast media as part of infusion and/or angiography, or for the passage of liquids for pressure measurements. They can be connected to an active medical device (e.g. pressure monitor) with pressure resistance according to label.

Technical Description
Fluid Extension Lines maintain an air-tight and liquid-tight connection to the environment between the tube ends to allow the channelling of liquids.
Lengths: 15 cm – 200 cm
Materials: DEHP-free Polyvinylchloride (PVC) or Polyurethane reinforced (PUR)
Rated burst pressures 500 psi, 1200 psi according to label
Luer connection: Luer Lock or rotating adapter

Packaging and Sterilization Instructions
Do not process products with damaged inner packaging. To maintain the sterility the products must be packaged in a primary package according to EN 868 / EN ISO 11607. It is the responsibility of the processor that packaging is carried out under suitable conditions concerning cleanroom, contamination and endotoxins in accordance with the intended use. Therefore, process validation and routine process monitoring are necessary.

EO Sterilization parameters to be used:
EO concentration: 800 mg/l
Relative Humidity: 30 - 90 % RH
Temperature: 48 - 58 °C
Air removal: 10 min
Conditioning: 30 min
EO Injection: max. pressure 40 kPa
Exposure Time: 480 min
EO Removal: -30 kPa
Flushing: 5 air change cycles
Aeration: 168 h at 0 – 40 °C

Warnings:
1. Instructions must be read carefully prior to use.
2. The product should not be used for anything other than the intended use.
3. Use proper aseptic techniques when handling.
4. The product is intended for SINGLE USE ONLY.

Contra-Indications: None known.
Societies must be taken into account accordingly.
Indications, contraindications and use of the product must always be based on current medical standards. The guidelines and recommendations of the respective medical societies must be taken into account accordingly.

Warnings:
1. Instructions must be read carefully prior to use.
2. The product should not be used for anything other than the intended use.
3. Use proper aseptic techniques when handling.
4. The product is intended for SINGLE USE ONLY.

Storage and Transport:
Store and transport in the original packaging under clean and dry conditions away from sunlight.

Instructions for use:
1. Before use, the label must be checked for the maximum allowable pressure conditions and ensure that these are sufficient for the intended application.
2. Aspirate and set up the line according to hospital protocols for the procedure.
3. The presence of air in the system may lead to embolism and in rare instances to death.
4. All connections should be finger tight. Over tightening may lead to cracks and leaks.
5. Once connected, verify there is no leakage. Leakage may lead to contact with biohazardous material.
6. Re-use of a single-use device creates a potential risk to the patient or user. This may lead to a contamination and/or impairment of functional capability.
7. Conformation and/or limited functionality of the device may lead to injury, illness or death of the patient. If this device comes into contact with medicines or disinfectants containing such similar ingredients as medium-chain triglycerides (MCTs), dimethylacetamide (DMAc), dimethylsulfoxide (DMSO), 2-Propanol (Isopropanol, IPA), allylated (DHA) or phenolaldehyde, cracks or leaks may occur.

Internal hospital guidelines must be observed for disposal, especially if there is a risk of infection.
Serious incidents:
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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