



# Scalable Oncology Insights with dWMS

5X Methylation Profiling and  
Structural Variant Detection

This case study highlights the utility of 5x direct Whole Methylome Sequencing (dWMS) for multi-omic biomarker discovery in oncology. By downsampling 30x ONT datasets from ovarian, lung, and lymphoma tumor samples, we demonstrate that low-pass dWMS enables cost-effective, high-confidence detection of genome-wide CpG methylation and large structural variants (SVs). These findings position dWMS as a powerful platform for translational research and clinical assay development.

## Study Objectives

To evaluate whether 5x dWMS provides sufficient resolution to detect genome-wide methylation and large structural variants in cancer tissue, using repurposed oncology datasets—supporting applications in early biomarker discovery, translational assay development, and low-input tumor profiling.



# Methods

## Samples

Flash-frozen human ovarian, lung, and lymphoma tumor tissues (n = 1 each) were selected according to the client's original study protocol.

## Sequencing Workflow

Genomic DNA (gDNA) was extracted and prepared using Renew Biotechnologies' (Renew) proprietary direct Whole Methylome Sequencing (dWMS) protocol. Native-read sequencing was performed on the Oxford Nanopore Technologies (ONT) PromethION platform, with reads aligned to the GRCh38 human reference genome using Minimap2.

## Downsampling & Methylation Analysis

Each 30x dataset was downsampled to 15x, 10x, and 5x genomic coverage using a custom pysam (<https://pysam.readthedocs.io/en/latest/faq.html>)-based pipeline. Aligned BAMs were converted to WBL's proprietary .CH3 format, which compresses file size by ~95% while retaining base-resolution CpG detail.

Methylation analysis was conducted using ModSeqR, Renew's native-read R package for ONT methylation data. The tool supports CpG summarization by site, genomic region, or annotation window, and includes built-in modules for quality control, sliding window analysis, and differential methylation detection.

High-confidence CpG calls were defined as those with model-predicted probabilities  $\geq 90\%$ , and high-quality reads as those with base quality scores  $\geq 10$ .

### Renew CH3 File Format

Proprietary file format that compresses whole genome ONT methylation data by ~95%, enabling rapid downstream analysis with full-resolution CpG.

### ModSeqR

- Renew's custom R package for native-read methylation analysis
- Summarizes CpGs by site, sliding window, or BED annotations
- Performs QC and regional differential methylation analysis
- Fully compatible with CH3 files generated from ONT dWMS



## Structural Variant Detection

Large insertions, deletions, duplications, and inversions were identified and visualized at each coverage level using Samplot (<https://github.com/ryanlayer/samplot>) and custom scripts optimized for ONT native-read alignments.

## Key Findings

### Genome-Wide Methylation Coverage and Quality

High-coverage dWMS was performed on ovarian, lung, and lymphoma tumor samples (n = 1 each) using ONT PromethION with one sample per flow cell. All samples achieved >30x whole-genome coverage with alignment rates averaging 99.5%, confirming high-quality data suitable for methylation analysis and SV detection (Table 1). These datasets were subsequently downsampled to 15x, 10x, and 5x to evaluate dWMS performance at lower depths.

Given the consistency in sequencing quality and coverage across tumor types, data from the three samples were averaged to illustrate general dWMS performance trends across multiple cancer tissues at varying depths.

Cancer Type	Output	Coverage (x)	% Alignment	Avg Read Length (bp)
Ovarian	107.3	33.5	99.4	8353
Lung	138.4	43.2	99.6	8137
Lymphoma	111.8	37.1	99.5	5288


**Table 1. Summary of sequencing performance metrics across tumor samples.** dWMS was conducted on three human tumor types using ONT PromethION. Data were downsampled from 30x to assess methylation and SV detection at lower genomic coverage.

This study demonstrates the power of nanopore-based dWS to deliver scalable, high-resolution methylation data without bisulfite conversion<sup>1,2</sup> or PCR amplification<sup>3</sup>. Even at 5x genomic coverage, dWMS achieved 100% 1x CpG coverage across the genome (Figure 1) -far exceeding the reach of traditional Illumina EPIC arrays<sup>4</sup>, which interrogate only ~ 3% of unique CpGs (~850,000 of ~28 million). In contrast, 5x dWMS provides native-read access to nearly all CpG sites, enabling base-resolution methylation analysis with full-molecule context.



CpG coverage depth scaled predictably with sequencing input: at 10x, ~90% of CpGs were covered at  $\geq 5x$ , and at 15x, nearly 80% reached  $\geq 10x$  coverage. Notably, even at 5x, dWMS captured more than 3x the number of CpGs covered at  $\geq 10x$  depth than the EPIC array.


Across all coverage levels, the proportion of CpG sites classified as high-confidence ( $\geq 90\%$  model probability) and high-quality (base quality  $\geq 10$ ) remained consistent, averaging 78.4% and 93.6%, respectively—highlighting the reproducibility and analytical integrity observed in this dataset using dWMS.



**Figure 1. Genome-wide CpG site coverage across downsampled dWMS sequencing depths.** Bar plots show the percentage of CpG sites covered at  $\geq 1x$ ,  $\geq 5x$ ,  $\geq 10x$ , and  $\geq 15x$  across 5x, 10x, 15x, and 30x genomic coverage using direct Whole Methylome Sequencing (dWMS). Data represent pooled results from three tumor types (ovarian, lung, lymphoma;  $n = 3$ ). Even at 5x genomic coverage, 100% of unique CpGs were covered at least once, with increased depth achieved at progressively higher sequencing coverages. These findings demonstrate the scalability of dWMS for flexible methylation profiling—from broad discovery to high-confidence regional analysis.

## Large Structural Variant Detection

To assess the performance of dWMS for structural variant (SV) detection we selected four large rearrangements spanning insertions, deletions, inversions, and duplications from the ovarian, lung, and lymphoma samples (Figure 2). Although not cancer-specific, these variants were chosen for their size and clarity to illustrate the resolution and consistency of SV detection across coverage levels. Each variant was confidently identified at 5x, 10x, 15x, and 30x, underscoring the power of native-read sequencing to detect complex rearrangements missed by short-read methods.



**Figure 2. Representative structural variants detected across 5x–30x dWMS coverage.** Variants include a 471 bp insertion (A, ovarian), 997 bp deletion (B, lymphoma), 9.9 kb inversion (C, lung), and 2.25 kb duplication (D, ovarian). All were confidently detected across depths, including 5X sequencing coverage, and visualized with Samplot, demonstrating the robustness of native-read SV detection at low-pass coverage.

## Conclusions

Renew's 5x dWMS platform delivers high-resolution methylation and SV profiling from a single native-read assay. Even at reduced coverage, dWMS achieves broad CpG detection and resolves large genomic rearrangements, offering a cost-effective, scalable solution for biomarker discovery and assay development in oncology research.



## Reference List

1. Tanaka K, Okamoto A. Degradation of DNA by bisulfite treatment. *Bioorganic & medicinal chemistry letters*. 2007;17(7):1912-1915.
2. Olova N, Krueger F, Andrews S, et al. Comparison of whole-genome bisulfite sequencing library preparation strategies identifies sources of biases affecting DNA methylation data. *Genome biology*. 2018;19(1):1-19.
3. Kanagawa T. Bias and artifacts in multitemplate polymerase chain reactions (PCR). *Journal of Bioscience and Bioengineering*. 2003/01/01/ 2003;96(4):317-323. doi:[https://doi.org/10.1016/S1389-1723\(03\)90130-7](https://doi.org/10.1016/S1389-1723(03)90130-7)
4. Pidsley R, Zotenko E, Peters TJ, et al. Critical evaluation of the Illumina MethylationEPIC BeadChip microarray for whole-genome DNA methylation profiling. *Genome biology*. 2016;17:1-17.

### Partnering with Renew Biotechnologies

Unlock The Full Potential Of Your Genomic  
Research With Renew Biotechnologies.

Learn more at [www.renewbt.com](http://www.renewbt.com)

Contact us at [info@renewbt.com](mailto:info@renewbt.com)