


## Introduction

- **Robust preclinical models** are required for the evaluation of novel molecular radiotherapies (MRT):
  - accurately mimicking human tumor microenvironments
  - Expressing optimal levels of target proteins.
- Two biomarkers of interest in the field of MRT and Oncology are:
  1. **Fibroblast activation protein (FAP)**: Minimally expressed in regular tissues, upregulated in cancer-associated fibroblasts within the tumor microenvironment.
  2. **Prostate-specific membrane antigen (PSMA)**: overexpressed in most prostate cancers.
- Oncodesign Services already had **PSMA**-overexpressing models available, but only on Balb/c Nude mice. This strain is radiosensitive due to a mutation and has become unavailable for purchase in the recent months.
- Furthermore, inherently expressing **FAP** models are available, but have minimal FAP-expression
- To close the gap in our catalogue, **Oncodesign Services has recently developed novel FAP and PSMA-overexpressing models in non-radiosensitive mouse strains.**

## Methods

### Development of two FAP-overexpressing mouse models

- Transduction of two cell lines using a lentiviral vector encoding for FAP and ZsGreen
  - HT-1080: FAP-negative, 0% FAP+ cells
  - U-87-MG: low FAP-positive, 15% FAP+ cells
- Confirmation of *in vitro* FAP expression using immunofluorescence and flow cytometry
- Subcutaneous engraftment of  $10 \times 10^6$  HT-1080-FAP or U-87-MG-FAP cells in Swiss Nude mice
  - Twice weekly monitoring of body weight and tumor volume
  - Harvesting of tumors at ethical endpoint and *ex vivo* confirmation of FAP expression




### Development of two PSMA-overexpressing mouse models

- Subcutaneous engraftment of  $10 \times 10^6$  22RV1-Luc-mCherry or Ln-Cap C4.2 cells (with matrigel) **in male BRGSF mice**
- Before the engraftment of 22RV1-Luc-mCherry, one group of mice was castrated
- Use of enriched diet (JL Mouse) to prevent known cachexia
- Twice weekly monitoring of body weight and tumor volume


### Development of two FAP-overexpressing cell lines

- HT-1080 and U-87-MG cells were successfully transduced with a lentiviral vector encoding for ZsGreen and FAP.
- Expression of both ZsGreen and FAP was confirmed by immunofluorescence and flow cytometry.
  - **HT-1080-FAP**:  $63.00 \pm 0.57$  % FAP-positive cells
  - **U-87-MG-FAP**:  $95.80 \pm 0.28$  % FAP-positive cells




### Development of FAP-overexpressing mouse models

- Swiss Nude mice subcutaneously engrafted with **HT-1080-FAP** cells presented with rapid tumor growth, which was associated with important body weight loss
- Swiss Nude mice subcutaneously engrafted with **U-87-MG-FAP** cells presented with a very slow, but homogenous tumor growth.



### Ex vivo confirmation of FAP expression using flow cytometry

- *Ex vivo* flow cytometry confirms FAP expression in both tumor models:
  - **HT-1080-FAP tumors**:  $46.85 \pm 12.11$ % FAP-positive cells
  - **U87-MG-FAP tumors**:  $69.88 \pm 17.20$ % FAP-positive cells
- *Ex vivo* FAP expression results correspond to the *in vitro* FACS results.




## Results

### Development of PSMA-overexpressing mouse models


#### 1. 22RV1-Luc-mCherry

- **Non-castrated** BRGSF mice subcutaneously engrafted **22RV1** cells present with rapid and homogenous tumor growth.
- Tumor growth in **castrated** BRGSF mice subcutaneously engrafted with **22RV1** cells was slightly slower and more heterogenous.



#### 2. Ln-Cap C4.2

- BRGSF mice subcutaneously engrafted with Ln Cap C4.2 cells presented a very heterogenous tumor growth, which was associated with a slow, but significant body weight loss despite the use of enriched diet



## Conclusion

- **Four novel preclinical mouse models were successfully developed at Oncodesign Services:**
  - Two FAP-overexpressing models with variable expression of FAP, mimicking clinical heterogeneity (negative HT-1080, slightly positive U-87-MG, positive HT-1080-FAP and overexpressing U-87-MG-FAP cells)
  - Two PSMA-overexpressing models, developed in non-radiosensitive, BRGSF mice
  - Further histological validation of the models is on-going.
- These models offer **new opportunities to validate novel imaging agents and molecular radiotherapies targeting FAP and PSMA**, facilitating their translation to clinical applications.