

LIBERTY STAIRLIFT



OWNER'S MANUAL



Peace of mind comes in many forms. It's knowing that you can go up and down the stairs on your own. It's regaining your independence. It's feeling safe every time you use your stairlift, knowing you don't have to worry about power outages or parking it in a specific spot on the track. It's knowing you've made the best decision and purchased the right stairlift. It's the feeling you get when you own an Atlas Liberty stairlift.

DANGER: Indicates a hazardous situation that will result in death or serious injury if not avoided.

WARNING: Indicates a hazardous situation that could result in death or serious injury.

CAUTION: Indicates a hazardous situation that could result in minor or moderate injury.

NOTE: Provides important information about the equipment or installation.

This energy-efficient stairlift is built with the "Atlas difference"—it's simple to use, durable, and lightweight.

This owner's manual must be read and understood in its entirety before attempting to operate the equipment. All operating procedures must be followed to avoid injury and/or property damage.

TABLE OF CONTENTS

SAFETY	4
Definitions of safety	4
Special precautions.....	4
INTRODUCTION.....	5
Device name: Liberty.....	5
Read and understand.....	5
Purchasing information	5
Technical specifications	5
Code declaration	5
Main features of the stairlift.....	6
Safety features.....	6
Controls.....	6
OPERATION.....	7
Stair seat	7
CARE AND MAINTENANCE.....	9
Battery	9
Care and cleaning.....	9
Maintenance/service sheet.....	10
TROUBLESHOOTING	12
Manual control operation.....	12
WARRANTY	14

SECTION 1

SAFETY

SAFETY DEFINITIONS



This safety alert symbol appears with safety statements. It means caution, be alert, your safety and that of others is at stake! Please read and follow the message that follows the safety alert symbol.

CAUTION

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE

Indicates a situation that could damage the device and/or the environment or cause the device to malfunction.

NOTE: Indicates a condition that must be met for the device to function as intended.

SPECIAL PRECAUTIONS



Pay attention to the following special precautions when using your stairlift.

- Do not exceed the maximum load capacity of 300 lb.
- Ensure there are no obstacles under the stairlift before using it.
- Do not disable the safety equipment or switches provided with the device.
- Do not attempt to repair the stairlift yourself. Contact your Atlas dealer for assistance.
- Do not allow children to use the stairlift or play near it.
- This product is designed solely to lift people within the specified load capacity. It must not be used to transport goods or for any other purpose.
- Read all instructions in this manual before using the stairlift.

SECTION 2

INTRODUCTION

NAME: LIBERTY

Indications for use: The Atlas Liberty is designed to facilitate the transfer of persons with reduced mobility up and down the levels of a residence or commercial application.

READ AND UNDERSTAND

This manual provides instructions for the proper use and maintenance of your Liberty stairlift. Please refer this user manual for operating instructions. Any modification of the equipment without the manufacturer's written authorization is prohibited and will void the warranty.

INFORMATION

If you have any questions about the operation and maintenance of your Liberty stairlift, please contact your dealer.

Be sure to complete the separate warranty form and return it within ten (10) days of purchase to register your unit.

SPECIFICATIONS

Visit [Atlas.com](https://atlas.com) for specifications on your stairlift model.

DECLARATION

Code requirements for the Liberty may vary depending on location. It is the installer's responsibility to contact the state, city, or local code enforcement office and determine all regulations to which the stairlift and installation are subject. This must be done before installing the Liberty. Certified by Intertek (ETL) to ASME A18.1-2017 Section 4 and Section 7.

INTRODUCTION



Safety features:

- Complies with ASME 18.1
- Integrated centrifuge Over-speed brake
- Swivel seat lock
- Chassis obstruction sensors
- Footrest obstruction sensors

Warranty:

- 1-year warranty on Drive train
- 1-year warranty on parts
- 1-year warranty on battery



Specifications

Maximum capacity	315 lbs
Minimum folded width	14.25
Lifting/lowering speed	15-20 fpm
Seat back height	18
Seat bottom	13" x 15.5"
Overall Seat Depth	13

SECTION 3

OPERATION

1. To start the stairlift, press the red "On/Off" switch located on top of the carriage to the "On" position and ensure that The key switch (optional) is in the "On" position. When the stairlift is turned on, a single beep indicates that the stairlift is ready for use.
2. To facilitate access to the chair, the armrest on the entry side can be raised. When you sit down on the seat at the top landing, make sure that the seat is fully swiveled, facing away from the stairs, and locked in the safety position. The seat is swiveled by pressing the swivel release handle located on the side of the seat. The footrest can be lowered to facilitate entry into the chair and to keep the feet in the correct position during transport.



CAUTION

Never attempt to use the device while standing. Always sit in the center of the seat with the safety belt fastened around your waist.

3. Make sure you are seated properly, your seatbelt is fastened, and the armrests are in the down position before using the stairlift
. To operate the stairlift, press and hold the hand control located on the armrest in the direction you wish to travel.



CAUTION

Keep your feet centered on the footrest when going up and down the stairs. Never let your feet hang off the footrest when using the stairlift.

4. If you release the switch, the stairlift will stop. Keep pressing the switch or button until you reach the landing. The stairlift will stop gently at the appropriate preset point.
5. At the upper landing, rotate the seat away from the stairs until it is fully rotated and locked, and remove your seat belt before standing up. To exit the stair seat safely, use the armrests for support, stand up, and step off the footrest. If you do not intend to use the device again immediately, the seat, footrest, and armrest can be folded away so as not to obstruct people walking up or down the stairs.



CAUTION

When transferring to or from the seat, use the armrest for support. Make sure the seat is securely locked before getting in or out.

6. To operate the stairlift remotely, press and hold the appropriate up or down button on the remote control call/send button on the remote control. Remote controls are usually located on the upper and lower landings of the staircase.

If you are operating the stairlift remotely for someone else, make sure they are seated and their safety belt is fastened.

NOTE: A single beep indicates that the stairlift has encountered an obstacle. Clear the obstacle before continuing.

SECTION 4

CARE AND MAINTENANCE

BATTERY

The stairlift must remain always charged. When the stairlift stops within three feet of an upper or lower landing, it automatically recharges. The power supply must remain always plugged in, and the stairlift can be left charging indefinitely. A fully charged battery will display a steady green light, while a battery that is charging will display a steady orange light on the charger.

If the stairlift emits intermittent beeps, it means that it has stopped far from the upper or lower charging stations or that the charger has been unplugged. Make sure the charger is plugged in and move the stairlift at each end of the rail to ensure that the batteries remain fully charged.

Note: The stairlift must be turned on to allow the batteries to charge. The batteries will not recharge when the stairlift is turned off.

Batteries should be replaced every 1 to 2 years on average, depending on the degree of use. Contact your local dealer for information on how to replace the batteries.



CAUTION

Non-compliant batteries may cause product malfunction.

Annual inspection is essential to identify the condition of the battery.

MAINTENANCE AND CLEANING

The stairlift should not require any technical maintenance to continue operating at full capacity. The drive system does not use grease or other lubricants, and no adjustments should be necessary.

However, the stairlift should be inspected annually. The stairlift and rail should be kept clean. Dust can accumulate on surfaces and dirt can appear. The seat, footrest, frame, and upholstery can be cleaned with a damp, soapy cloth, followed by another clean, damp cloth. Do not use abrasive cleaners, as they may damage the finish. For heavy soiling, clean by rubbing lightly with a soft brush.

SECTION 5

TROUBLESHOOTING

MANUAL CONTROL

If your stairlift is not working and the user cannot get off the device on the stairs, another person can use the optional manual control tool to lower or raise the stairlift to a landing. However, please follow the instructions on the lower safety flap of the stairlift and set the "On/Off" switch to "Off." Insert the manual control tool into the hole in the lower safety flap until it engages with the motor shaft, then turn in the desired direction.



CAUTION

Do not use the stairlift when the manual control tool is engaged.

WARRANTY

SECTION 6

ONE-YEAR LIMITED WARRANTY

PRODUCTS COVERED: LIBERTY

Your stairlift comes with a separate warranty page on page 15. You must return this page within ten days. (10) days after installing the device to register it. This warranty policy page must remain in your owner's manual for your records. Do not tear out this page.

Dealer/Installer Name: _____

Address _____

Phone: _____

Email: _____

If you need warranty service or repair under warranty, your dealer will need the information below to receive information from the factory or order parts for your stairlift.

Date of purchase:

Stairlift serial number

THIS EXPRESS WARRANTY REPLACES ALL OTHER WARRANTIES OR CONDITIONS, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE FOR A PARTICULAR PURPOSE, AND THERE IS NO WARRANTY THAT EXTEND BEYOND THE DESCRIPTION OF THE LIMITED WARRANTY SET FORTH IN THIS DOCUMENT. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER LEGAL RIGHTS THAT VARY FROM STATE TO STATE.

Atlas warrants to the original purchaser of a Liberty stairlift manufactured for us that it is free from material and mechanical defects in the motor, gearbox, worm gear, and rack for a period of one year. all other components and parts for a period of 1 year. Batteries supplied by Atlas are limited to one per unit. 1 year from the date of retail purchase. The

Liberty is not approved for outdoor use.

No warranty, express or implied, including merchantability or fitness for a particular purpose, shall survive the expiration of one year from the date of original purchase of the unit. The limit of liability of Atlas and its dealer hereunder shall be the purchase price of the device. Some states do not allow limitations on the duration of an implied warranty or the exclusion or limitation of incidental or consequential damages or legal remedies, so the above limitations may not apply to you. All warranty claims must be reported to the dealer from whom the stairlift was purchased, as they are responsible for handling your warranty claim. The dealer must contact Atlas Customer Service at 866-378-6648 and provide the stairlift's serial number and a description of the events leading to the warranty claim. Dealers may charge a labor fee.

labor, service, travel, or other costs associated with performing repairs, and these costs are not covered by this limited warranty. Repairs or replacement work may be performed by a repair shop.

The limited warranty does not cover repairs or replacements made by anyone other than the dealer because of a defect in material or workmanship. However, the limited warranty does not cover the costs or expenses assessed by any other person or company performing such repairs or replacement work. All parts used to replace defective materials must be genuine Atlas parts to be covered by this limited warranty.

PLEASE COMPLETE ALL FIELDS AND RETURN WITHIN TEN (10) DAYS OF PRODUCT INSTALLATION.

Fax the completed form to 1-866-234-5680 or mail it to: Atlas, ATTN Warranty Department, 2075^{47th} Street, Sarasota Florida 34234.

PRODUCT INFORMATION

Model: _____

Serial Number: _____

Date of purchase: _____

Date of installation: _____

INFORMATION ABOUT THE INSTALLER

Company name: _____

Contact name _____

Address _____

Phone _____

Fax _____

Email _____

INFORMATION ABOUT THE REQUEST

Scooter Electric wheelchair

Wheelchair N/A

Year _____

Manufacturer _____

Model: _____

PURCHASER INFORMATION

Name: _____

Address: _____

Phone: _____

Email: _____

ADDITIONAL INFORMATION

How did you hear about Atlas?

Atlas distributor

Internet

Magazine

Friend or acquaintance

I saw the Atlas product somewhere

Other _____

I purchased my Atlas lift for the following reasons

Style/appearance

Atlas representative

Previous experience

Ease of use

Recommendation

Price/Value

WARRANTY

Page intentionally left blank

Installation and use must comply with applicable codes and standards in force in the relevant jurisdiction.

