

PUSH-BUTTON SAFETY BLOOD COLLECTION SET Instructions for Use

Please take time to read the instructions before use.

NOTE: Gloves should be worn when using this device. For blood collection using a syringe, or for short-term IV administration, refer to the manufacturer's instructions included in the box.



STEP 1

Twist sheath to remove from Luer adapter and attach tube holder, transfer device, or syringe.



STEP 2

Twist and pull off needle cover and perform phlebotomy per your facility's procedure.



STEP 3

Hold the wings and insert needle into the vein. Observe for the presence of blood in the chamber.

STEP 4

Collect blood specimen(s) per your facility's procedure.



STEP 5

Cover the puncture site with a sterile gauze.



STEP 6

Depress the button to retract the needle in-vein and perform visual check, taking care not to impede device retraction. Apply pressure and/or bandage.

STEP 7

Dispose of all used materials in the nearest sharps container per your facility's procedure.



MYCO
MEDICAL
TheSharpChoice

NMSDC Certified Diversity Supplier

1.919.460.2535
2015 Production Drive
Apex, NC 27539 USA
www.mycomedical.com
info@mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2026 MYCO Medical, Inc. All rights reserved.

MYCO® is a registered trademark of MYCO Medical, Inc.

LIT-00981 Rev 1 Eff 08.19.2026