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Clauss, A. Acta Haemat. 17:237-246, 1957.

NCCLS: Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays. 2nd edition. Approved guideline. NCCLS Document H21-A3. Wayne, PA., 1998.

NCCLS: Procedure for determining fibrinogen in plasma. Approved guideline. NCCLS Document H30-A2. Wayne, PA., 2001.

Musgrave, K.A., Bick, R.L.: Quality Assurance in the Hemostasis Laboratory. In Bick, R.L., et al editors: Hematology: Clinical and Laboratory Practice. Vol. 2, pp 1309-1315. Mosby. St. Louis, MO., 1993

HEMOFIBRIN

Quantitative determination of Fibrinogen concentration

	Manufacturer: RAL TÉCNICA PARA EL LABORATORIO, S.A. Av. Mare de Déu de Montserrat 51 08970 Sant Joan Despí (Barcelona) - Spain
IVD	FOR IN VITRO DIAGNOSTIC USE
	STORE AT 2 - 8° C
	READ INSTRUCTIONS

REF	Cont.
Reference: GN 1045	Package: 5 x 2 mL bovine thrombin + 125 mL imidazole buffer + 1 mL calibrator

PRINCIPLE:

Gernon Hemofibrin kit is intended for use in the quantitative determination of fibrinogen in plasma samples.

The thrombin clotting time fibrinogen assay is based on the method originally described by Clauss. In the presence of high concentrations of thrombin, the time required for clot formation in dilute plasma is inversely proportional to the fibrinogen concentration.

REAGENTS:

Hemofibrin: Lyophilized buffered bovine thrombin.

Imidazole buffer: Imidazole buffer in saline, pH 7.4 ± 0.2, with 0.1% sodium azide as a preservative.

Fibrinogen calibrator: Human plasma collected with sodium citrate anticoagulant.

Reconstitute Hemofibrin reagent with 2 mL of distilled water. Agitate gently until solution is complete.

Reconstitute the calibrator plasma with 1.0 mL distilled water. Agitate gently until solution is complete.

STORAGE:

The unopened reagents are stable till the end of the expiry date if stored at 2-8°C.

The reconstituted Hemofibrin reagent is stable for 7 days at 2-8°C or 8 hours at 15-30°C or may be frozen within 4 hours for use within 30 days. Thaw rapidly at 37°C. Do not refreeze.

The reconstituted calibrator plasma is stable for 8 hours at 2-8°C.

WARNINGS AND PRECAUTIONS:

Imidazole buffer contains 0.1% sodium azide. Sodium azide under acid conditions yields hydrazoic acid, a toxic compound. Dilute with running water before discarding, flush with large volumes of water.

Each unit of source material used in the preparation of Fibrinogen Reference Plasma has been tested by an FDA licensed method and found non-reactive for HBsAg and negative for antibodies to HIV and HCV.

However, no known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS, or other infectious diseases. This product, like all materials of human origin, should be handled as potentially infectious biological material.

Take the necessary precautions for the use of laboratory reagents.

SPECIMEN:

Citrated plasma.

3.2% (0.109M) trisodium citrate anticoagulant is recommended for coagulation assays. Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90% of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22-24°C.

PROCEDURE:

Gernon Hemofibrin is suitable for use with manual, mechanical or photo-optical methods of end-point clot detection. For automated and semi-automated instruments, follow the instrument manufacturer's instructions.

For manual and mechanical test methods:

Prepare a minimum of three different dilutions of the reconstituted Fibrinogen calibrator plasma in imidazole buffer (for instance, 1:5, 1:10 and 1:20). The corresponding concentrations will be 2c, c and c/2 (c = fibrinogen concentration in calibrator). Two replicates of each dilution should be prepared.

Dilute quality control and patient samples 1:10 in imidazole and prewarm 0.2 mL of each dilution to 37°C for 4-6 minutes.

Add 0.1 mL of Bovine Thrombin reagent to prewarmed dilution and time clot formation. Do not prewarm thrombin.

Prior to use, mix well the reagents without shaking them.

Pipette directly into clean and dry cuvettes:

Add in cuvettes under stirring:	
Plasma (calibrator, control or sample diluted)	0.2 mL
Incubate at 37°C for 4-6 minutes	
Hemofibrin reagent	0.1 mL

The frequency of curve preparation is partially determined by the method of clot detection used. Always prepare a new curve with each change in reagent lots, instrumentation, or when controls no longer fall within established ranges.

CALCULATION:

Write down the clot time. Calculate the average value of calibrator replicates. Use all three calibrator points (dilutions 2c, c and c/2 of calibrator plasma) to construct a log-log curve that plots fibrinogen concentration vs. clotting time. Find the clotting time of quality control and patient samples on the curve and read the corresponding fibrinogen value. If clotting times for the 1:10 dilution fall outside the linear curve, prepare 1:5 or 1:20 dilutions as needed. If the sample is diluted 1:5, divide the result from the standard curve by 2; if the sample was diluted 1:20, multiply the curve result by 2 to get the final result.

REFERENCE RANGE:

The normal range for fibrinogen levels in human plasma is considered to be 150-350 mg/dL. Each laboratory should establish its own mean normal and normal range because of variances among different laboratories.

PERFORMANCE CHARACTERISTICS:

Accuracy:

A low, a normal, and a high fibrinogen plasma were tested in multiple laboratories using Hemofibrin reagent. The results were compared to results obtained using other manufacturer's reagents in multiple labs.

Sample	Hemofibrin	n	All reagents	n
Low	144 mg/dL	10	163 mg/dL	195
Normal	294 mg/dL	10	297 mg/dL	195
High	488 mg/dL	10	474 mg/dL	390

Precision:

A low, a normal, and a high fibrinogen plasma were tested on multiple days using Hemofibrin reagent on a photo-optical instrument. Ten standard curves were determined on each test day, for a total of 30 standard curves. The percent CV was determined to be 5.9% (low), 3.4% (normal), and 2.9% (high).

LIMITATIONS:

Blood must be immediately added to trisodium citrate anticoagulant (ratio 9:1) and gently mixed. EDTA and heparin are unsuitable anticoagulants. The citrate concentration must be adjusted for patients with hematocrits above 55%.

Hemolysis can cause clotting factor activation and end point detection interference.

Freezing and thawing of plasma that contains residual cells will generate damaged cell membranes that can affect results.

Acute inflammatory reactions can elevate circulating Factor I (fibrinogen).

High Fibrinogen Degradation Products (FDP) may prolong clotting times, especially when the fibrinogen level is below 150 mg/dL.

In patients with qualitative fibrinogen abnormalities, the thrombin clotting time assay may indicate decreased fibrinogen. The quantitative fibrinogen results may be normal on these same samples if tested by other methods.

Heparin does not interfere at therapeutic levels. However, very high heparin levels may cause low fibrinogen results.

High paraprotein levels, thrombin antibodies, and drugs that activate the fibrinolytic system can interfere with fibrinogen assays.

LITERATURE:

Clauss, A. Acta Haemat. 17:237-246, 1957.

NCCLS: Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays. 2nd edition. Approved guideline. NCCLS Document H21-A3. Wayne, PA., 1998.

NCCLS: Procedure for determining fibrinogen in plasma. Approved guideline. NCCLS Document H30-A2. Wayne, PA., 2001.

Musgrave, K.A., Bick, R.L.: Quality Assurance in the Hemostasis Laboratory. In Bick, R.L., et al editors: Hematology: Clinical and Laboratory Practice. Vol. 2, pp 1309-1315. Mosby. St. Louis, MO., 1993

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Clauss, A. Acta Haemat. 17:237-246, 1957.
NCCLS: Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays. 2nd edition. Approved guideline. NCCLS Document H21-A3. Wayne, PA., 1998.
NCCLS: Procedure for determining fibrinogen in plasma. Approved guideline. NCCLS Document H30-A2. Wayne, PA., 2001.
Musgrave, K.A., Bick, R.L.: Quality Assurance in the Hemostasis Laboratory. In Bick, R.L., et al editors: Hematology: Clinical and Laboratory Practice. Vol. 2, pp 1309-1315. Mosby. St. Louis, MO., 1993

IU/GN/1045/RO/0310/EF