

MRI Sacrum/SI Joints WO

Technical Parameters

Anatomic Coverage:	SI joints, sacrum, coccyx					
Coil:	Torso					
Scanner Preference:						
Scanning Plane:	Axial T1	Axial PDFS	Sagittal T2FS	Coronal Oblique T1	Coronal Oblique PDFS	Coronal STIR
Slice Thickness:	3 mm	3 mm	3 mm	3 mm	3 mm	4 mm
Gap:	0.5 mm	0.5 mm	0.5 mm	0.5 mm	0.5mm	1 mm
Weighting:	T1	PD Fat-Sat	T2 Fat-Sat	T1	PD Fat-Sat	STIR
Mode:	TSE	TSE	TSE	TSE	TSE	TSE
FAT SAT:	NO	YES	YES	NO	YES	YES
FOV:	18-22 cm	18-22 cm	18-22 cm	18-22 cm	18-22 cm	30-40 cm
MPR:	NO	NO	NO	NO	NO	NO
Notes:						

Prescribing the coronal oblique plane: parallel the sacrum on a sagittal localizer

