

Renew® Inserts

EN INSTRUCTIONS FOR USE

PLEASE READ CAREFULLY

GENERAL INFORMATION

Renew® Inserts are made of soft silicone and are safe and comfortable to use in the rectum. The inserts and the removable applicator are both single-use and non-sterile.

SIZING

Renew® Inserts come in two sizes: regular and large. Start with the regular size **Renew® Insert.** If it comes out unexpectedly or if you experience leakage, try using the large size.

INTENDED USE

The Insert is intended for the management of Accidental Bowel Leakage (ABL) in adult patients with faecal incontinence by providing a temporary rectal barrier.

INDICATIONS FOR USE

The **Renew® Insert** is indicated for use in adult patients experiencing faecal incontinence and is designed to help reduce or prevent involuntary leakage of stool when inserted into the rectum.

INTENDED USER

The device is intended for use by Healthcare Professionals, adult patients capable of self-insertion, and caregivers in home, clinical, and home care settings.

TARGET POPULATION

Adult patients experiencing faecal incontinence (FI) or Accidental Bowel Leakage (ABL).

ESSENTIAL PERFORMANCE

The essential performance of the device is its ability to remain securely positioned within the rectum and provide an effective barrier to help reduce or prevent involuntary stool leakage, while allowing safe insertion and removal without causing injury or device failure.

CLINICAL BENEFITS

Clinical study* data demonstrate that the Renew® Insert is safe and effective when used as intended.

Use of the device may provide the following benefits:
- Helps reduce the frequency of Accidental Bowel Leakage (ABL) episodes.
- A majority of users experienced a clinically meaningful reduction in leakage.
- Improvement in symptoms of faecal incontinence, as measured by Wexner score.
- High levels of user satisfaction and ease of use reported.

U.S. Clinical Study Renew 210 CLDI published in Dis Colon Rectum 2015

WARNINGS
- Do not use with suppositories.
- Discontinue use and consult your Healthcare Professional if you experience any of the following:
 - Bleeding haemorrhoids
 - Rectal bleeding or bloody stool
 - Anal/rectal pain, irritation, spasms, fissure or tearing
 - Undergone rectal surgery within the past 4 weeks
 - Displacement of the insert upward into the anus and it is not expelled with the next bowel movement.
- Do not use if the packaging has been deliberately or accidentally unsealed before intended use, or if the product is damaged.
- Do not re-use.
- Keep out of reach of children to prevent choking on small components.

CAUTIONS
- Consult instructions for use.
- Avoid multiple insertion attempts to limit the exposure of stool to other parts of the body; use a new insert if this happens.
- Do not use the device after the expiration date indicated on the packaging.

PRECAUTIONS
- Maintain good hygiene during insertion and removal to minimise the risk of infection.
- Consult a Healthcare Professional before use if you have a history of rectal or anal conditions (e.g. mild hemorrhoids or prior anorectal procedures).
- Use caution during insertion and removal to avoid irritation or minor injury.

CONTRAINDICATIONS
- Do not use the device if you have a known allergy or hypersensitivity to silicone or any of its components.
- Do not use in individuals with severe rectal or anal conditions (e.g. anal fissures, rectal prolapse, severe hemorrhoids, or infections) unless advised by a Healthcare Professional.
- Not recommended for use during active rectal bleeding or in the presence of acute or chronic anal infections.

SIDE EFFECTS
- False urge to empty your bowels.
- Mild irritation of anal area.
- Potential sensitivity to the materials used in the device.

COMPLICATIONS
- Allergic reaction to materials
- Device migration or incorrect positioning
- Incorrect size selection causing poor fit or irritation
- Failure to remove applicator

STORAGE CONDITIONS
- Keep dry
- Keep away from sunlight

REPORTING OF SERIOUS INCIDENTS

If a device malfunction or failure results in a serious incident, such as significant deterioration in health or death, you should contact Clinisupplies Ltd and your local Competent Authority as soon as possible. Details of your local Competent Authority can be found on your local government website.

Clinisupplies Ltd accepts no liability for any loss or harm that may arise if this product is not used in accordance with the Instructions for Use and the company's recommendations.

DISPOSAL

When used in a care facility, dispose of the used device in accordance with local clinical waste or healthcare waste policies.

When used at home, dispose of the device in normal household waste.

Please follow local regulations.

USER GUIDE STEP BY STEP

Preparation

1. Wash your hands, peel back the top of the package to fully expose the insert.
2. Place your index finger into the applicator and lift the insert out of its package.

Insertion

3. Sit comfortably on a toilet and insert from the front or back. Other common insertion positions are to raise yourself up onto one hip (while sitting on the toilet) or lie down on a bed (if this is more comfortable). Position the top disc of the insert so it's lightly resting on the anus. Note: If the insert falls off the applicator prior to insertion, discard it and use a new one.
4. Relax your muscles and gently slide the insert into the anus until you feel the applicator touches the anus.
5. Gently withdraw the applicator. The top disc of the insert should be resting inside the rectum, while the bottom disc remains just outside the anus.
6. Dispose of the applicator in a bin and wash your hands after use. The insert is now correctly positioned and works from the inside to help prevent ABL.

Removal

7. Simply pull on the bottom disc outside the anus to remove the insert and dispose of it. Wash your hands after direct contact with the used insert.

Replacement

8. For optimal protection, use a new Renew® Insert after each bowel movement.

Wearing instructions

The insert may be worn continuously, until your next bowel movement, or changed as advised by a Healthcare Professional.

You can bathe and be physically active while using a Renew® Insert.

Never use more than one insert at a time.

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