EN

Aquaflush® Sterile Stoma Stopper

Instructions for use

Aquaflush® Sterile Stoma Stopper is made from medical grade silicone that can be used to seal stoma sites at various parts of the body. The device is soft, low profile and comfortable and are unique as they are adjustable in length, they also come in a variety of shaft diameters depending on the size of the stoma.

A stoma connects a portion of a body cavity to the external environment. Parts of the gastroin-testinal tract which could be involved in the creation of a stoma are the oesophagus, stomach, duodenum, small bowel, and colon. Stomata in neonates are usually temporary and the most common conditions for which stomata may be fashioned in neonates are necrotizing entero-colitis, anorectal malformations and Hirschsprung's disease.

This guideline is applicable to neonates with stomata of the duodenum, jejunum, ileum, and colon.


Aquaflush® Stoma Stoppers were first indicated for patients with antegrade continence enemas (ACE) and then following this for patients with a Mitrofanoff stoma where they are required to seal the tract or maintain an open or unobstructed Stoma in between insertions of the enema catheter or drainage catheter (usually every 24 to 48 hours). They prevent Stomal Stenosis (closing of the stoma in-between catheterisations) and leakage from the stoma by forming a seal in between insertions of the catheter.

Indication:

It is indicated for neonates, children, and adult patients to manage antegrade continence enemas (ACE) and a Mitrofanoff stoma.

Preparing and inserting the Aquaflush® Stoma Stopper

- Wash your hands thoroughly before using this product
- ACE stoma site should be cleaned and dried before attempting to insert the stopper
- Using a clean sharp pair of scissors, cut the stopper to the length suggested by your Healthcare Professional. Cut through the indent along the shaft to ensure a smooth edge for insertion

- Gently insert the stopper into the tract, a small amount of water-based lubricant can be used to help
- The cap of the stopper should sit flush against your skin
- **Handy Hint:** make a note of the date you first used your **Aquaflush® Stoma Stopper** to ensure it is not used for an excess of 28 days

Applying Aquaflush® Dressing:

- The supplied extra thin sterile hydrocolloid dressing is designed to help keep the stopper in place
- Ensure the skin around the stopper cap is clean and dry before applying the dressing
- Remove the release film from the dressing and place centrally over the stopper cap
- Press gently to ensure it is secure
- **Aquaflush® hydrocolloid dressing** is single use only

Caring for your Aquaflush® Stoma Stopper:

- Wash with mild soap and water after each insertion and allow to dry before storage in a clean dry container
- If your stopper becomes damaged discard and replace with a new one
- **Handy hint:** Keep spare **Aquaflush® Stoma Stopper** and additional dressings in case of unforeseen event.

Storage Recommendations:

Store in a cool, dry place, away from direct sunlight and sources of heat.

Aquaflush® Stoma Stopper is available in four sizes in Fr gauge designated as nominal dimension to accommodate a range of anatomies and incontinence severities of intended user group.

The device is supplied sterile and preassembled with sterile **Stoma Stopper dressing** on individual, single-use blister packs.

Caution:

- Removal of the device for a prolonged period could result in stenosis (closing of the stoma) which may require urgent medical attention
- If irritation occurs under the hydrocolloid dressing, consult your Healthcare Professional for advice

Safety:

Must be discarded after a maximum of 28 days intermittent use and should be replaced.

If there is any damage to the device along the length of the stopper, do not insert to avoid dislodging a section. If you are concerned that a piece could have become dislodged, contact a Healthcare Professional immediately.

The **Aquaflush® Stoma Stopper** has a three-year shelf life from the date of manufacture.

- Do not use after the expiration date. For the expiration date check the box label
- Any serious incident that has occurred with this device should be reported to the manufacturer or competent authority
- Use a water-based lubrication for easier insertion

Contraindications:

Should not be used if there is a known silicone allergy. In case of irritation, consult your Healthcare Professional for advice. Do not remove unless advised to do so.

Disposal:

Stoma Stopper dressing and device can be removed and disposed of in line with local authority guidelines.

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Clinisupplies Ltd, 1 Blackmoor Lane, Croxley Park, Watford, Hertfordshire WD18 8GA, UK
Clinisupplies (Nanjing) Co., Ltd.,
EU REP **Advena Ltd.** Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta