SensaCath® Flex
Intermittente Katheter

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PLEASE READ CAREFULLY

GENERAL INFORMATION

SensaCath® Flex Intermittent Catheter has a sterile and activated hydrophilic surface. It is ready for insertion into the urethra without any delay. The catheter is stored in a foil pouch containing gas and is ready for use.

SensaCath® Flex intermittent catheter is a single use and latex free catheter which is made of phthalate free PVC (Polyvinyl Chloride) and coated with cross-linked PVP (Polyvinylpyrrolidone).

The surface layer of the catheter has a hydrophilic coating which binds with sterile water, creating a ready to use smooth slippery coating on the catheter.

SensaCath® Flex intermittent catheter is used for bladder drainage with individuals with bladder dysfunction. It is for single use only and must be discarded after use.

STORAGE

It should be stored in a cool and dry environment between 5º-35º C avoiding direct sunlight.

STERILISATION

Sterilised by irradiation.

WARNINGS

Intermittent self-catheterisation is a common procedure, however, you should always seek guidance from your healthcare professional and carefully read the instructions before use. Intermittent self-catheterisation using the **SensaCath® Flex intermittent catheter** should only be carried out under medical advice e.g. physician/nurse – and in accordance with the instructions given.

The sterile water in the catheter is harmless. Nevertheless, do not use the catheter if the package is open or torn, or has expired.

It is single use only. Reuse, reprocessing, cleaning, disinfection or resterilisation may compromise product characteristics which in turn may create a potential risk to the user. Catheterisation should be completed as soon as possible after the packaging is opened and the catheter should be removed once the bladder is empty.

User should be trained in the proper use of the **SensaCath® Flex Intermittent Catheter** which is ready to use. The user manual should be read carefully before use.

Clinisupplies accepts no liability or loss that may arise if this product is not used entirely according to the company's recommendations.

COMPLICATIONS

- Urinary tract infection
- Bleeding related to urethral trauma or infection
- Epididymitis, Urethritis, urethral stricture and bladder stones

SIDE EFFECTS

- Urinary tract infection/Sepsis
- Burning and/or pain in urethra
- Urethral bleeding
- Urinary system stone disease

CONTRAINDICATIONS

- Acute urethritis
- Acute prostatitis
- Acute epididymitis
- Patients with PVC-allergy
- Patients have calcareous urolithiasis

INSTRUCTION FOR USE

1. Begin by washing your hands and the area around the urethral orifice. It is important to avoid transferring bacteria onto the catheter from your hands or surrounding areas.
2. Attach the catheter package to a smooth surface using the sticky peel.
3. Peel the package open and the catheter can be removed and used immediately. Please take care not to drop or contaminate the catheter.

4. Preparation for insertion of the catheter:

 Loosen your non-touch insertion slider a bit and release it from the connector by pulling it downwards. Ten pinch the non-touch insertion slider again about 10 to 15 centimeters from the catheter tip.

5. Preparation for insertion of the catheter:

 a. Male:

 Straighten out urinary tract by lifting your penis with one hand. Avoid squeezing the penis, as this may make insertion of the catheter more difficult.

 b. Female:

 With one hand, part your labia to expose the urethral orifice.

6. a. Male:

 Without touching the clear catheter, use the non-touch slider to gently insert the catheter into the urethra approximately 2cm at a time.

 Just before the catheter reaches the bladder you may feel a slight resistance due to the sphincter.

 b. Female:

 Without touching the clear catheter, use the non-touch slider to end consecutive to gently insert the catheter into the urethra approximately 2cm at a time.

7. Continue to insert the catheter into the urethra until urine starts to flow.

 Note: The catheter can be connected to a urine drainage bag if there is no toilet nearby.

8. When the urine stops flowing, ease out the catheter 2 to 3cm. If the urine restarts, wait a few seconds before easing out the catheter another 1 to 2cm. Once the bladder is completely empty, remove the catheter slowly and dispose of it in a waste bin.

PACKAGING

Quantity per box: 30 pcs

Reference	Fr/Ch (2mm)	Length
PSCX10	10Fr/Ch (3.33mm)	40cm
PSCX12	12Fr/Ch (4mm)	40cm
PSCX14	14Fr/Ch (4.67mm)	40cm
PSCX16	16Fr/Ch (5.34mm)	40cm
PSCX18	18Fr/Ch (6mm)	40cm

Any serious incident that has occurred with this device should be reported to the manufacturer or competent authority.

CSP0737_0_v2 · 23-09-2026 · 11:44:11