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PLEASE READ CAREFULLY

GENERAL INFORMATION

SensaCath® Flex Compact Intermittent Catheter is a sterile, single use, latex-free medical device made of phthalate-free TPU (Thermoplastic Polyurethane). **SensaCath® Flex Compact** incorporates two smooth eyes for efficient drainage. The catheter is coated with a hydrophilic coating to reduce friction during its insertion and removal from the urethra or urinary bladder.

The **SensaCath® Flex Compact** is stored in a pouch with sterile water which keeps the catheter lubricated and smooth and activated, ready for use at any time. Once the **SensaCath® Flex Compact** is removed from the pouch, it can be immediately used.

INTENDED USE

SensaCath® Flex Compact are self-lubricating intermittent catheters intended to be inserted through the urinary meatus to the urinary bladder, by the patient/user or a healthcare professional, to provide urinary drainage, typically for an individual who is physiologically incapable of voiding. It is intended that each user waits until all the urine is drained completely from the bladder. Catheter should be removed once urine has stopped flowing.

INDICATIONS FOR USE

SensaCath® Flex Compact is used for bladder drainage with individuals who have bladder dysfunction. **SensaCath® Flex Compact** is effective in urinary retention and incontinence associated with neurological disorders, post-operative care, and other medical conditions. It is for single use only and must be discarded after use.

STORAGE

Do not store catheters where they are directly exposed to organic solvents, ionising radiation, UV light and high temperatures.

It should be stored in a cool and dry environment between 5° 35° C, avoiding direct sunlight.

STERILISATION

Sterilised by irradiation. Do not re-sterilise.

WARNINGS

Intermittent self-catheterisation is a common procedure, however, you should always seek advice from your healthcare professional and carefully read the instructions before use. Intermittent self-catheterisation using the **SensaCath® Flex Compact** should only be carried out under medical supervision by healthcare professionals – and in accordance with the instructions given.

The sterile water in the catheter is harmless. Nevertheless, do not use the catheter if the package is open or torn, or has passed its use-by date.

It is single use only. Reuse, reprocessing, cleaning, disinfection or sterilisation may compromise product characteristics, which in turn may create a potential risk to the user. Catheterisation should be completed as soon as possible after the package is opened and the catheter should be removed once the bladder is empty.

User should be trained in the proper use of the **SensaCath® Flex Compact** which is kindly to use. The user manual should be read carefully before applying.

Clinisupplies accepts no liability or loss that may arise if this product is not used entirely according to the company's recommendations.

COMPLICATIONS

- Urinary tract infection
- Bleeding linked to urethral trauma or ipsilateral
- Epididymitis, Urethritis, urethral obstruction and bladder stones – rare

POSSIBLE SIDE EFFECTS

- Urinary tract infection/Sepsis
- Burning and/or pain in urethra
- Urethral bleeding
- Urinary system stone formation

CONTRAINDICATIONS

- Patients with TPU allergy
- Patients have calcareous urolithiasis
- Acute urethritis
- Acute prostatitis
- Acute epididymitis
- Patients with PVC allergy

PRECAUTIONS

- If you feel any resistance during insertion, do not force the catheter. Pull the catheter back slightly, try and relax and then continue to introduce the catheter
- If you cannot advance the catheter into the bladder seek medical advice
- If urine is not flowing out, please check if the catheter is blocked or bent

INSTRUCTION FOR USE

SensaCath® Flex Compact is ready to use as soon as you open the packaging, so ensure you are ready to proceed with the catheterisation straight after opening it.

1. Begin by washing your hands and the area around the urethral orifice. It is important to avoid transferring bacteria to the catheter from your hands or surrounding areas.
2. Attach the catheter package to a smooth surface using the sticky pouch.
3. Peel the package open and the catheter can be removed and used immediately. Please take care not to drop or contaminate the catheter.

4. To loosen the non-touch insertion slider, release it from the connector by pulling it down. Then grip the non-touch insertion slider again about 10 to 15 centimeters from the catheter tip.

a. Male:

Straighten out your urethra by lifting your penis with one hand. Avoid squeezing the penis, as this may make insertion of the catheter more difficult.

b. Female:

With one hand, part your labia to expose the urethral orifice.

6. Without touching the clear catheter, use the non-touch slider to gently insert the catheter into the urethra approximately 2cm at a time. Just before the catheter reaches the bladder you may feel a slight resistance due to the sphincter.
7. Continue to insert the catheter into the urethra until urine starts to flow.

NOTE: The catheter can be connected to a urine drainage bag if there is no toilet nearby.

8. When the urine stops flowing, ease out the catheter 2 to 3cm. If the urine restarts, wait a few seconds before easing out the catheter another 1 to 2cm, when the bladder is completely empty, remove the catheter slowly and dispose of it in a waste bin.

CSF01326_v3
02-26

Clinisupplies

Clinisupplies Ltd. 1 Blackmoor Lane, Croxley Park, Watford, Hertfordshire WD18 8GA, UK
Tel: +44 (0)20 8863 4168 Email: info@clinisupplies.co.uk **www.clinisupplies.co.uk**

EU REP Yuhang Jinsuhi Manufacturing Limited, 1 Jinsuhi Road, Jinsuhi Development Street Yuhang District, 311121 Hangzhou, Zhejiang Province. People's Republic of China

EU REP SHANGHAI INTERNATIONAL HOLDING CORP. (GMBH) EUROPE
Eiffestraße 80, 20537 Hamburg, GERMANY

CE 0344