Liberty® Leg Bag Sleeve

EN INSTRUCTIONS FOR USE

Description
The **Liberty® Leg Bag Sleeve** is a medical device designed to provide comfortable and discreet support for urine collection bags.
It is made from a soft, skin-friendly, breathable fabric consisting of:
- **Polyamide (92%)**
- **Elastane (7%)**
- **Viscose (1%)**

The sleeve helps reduce sweating, provides a secure fit, and allows discreet wearing under clothing.

Intended Use
The **Liberty® Leg Bag Sleeve** is intended to comfortably secure a urine collection device (urine leg bag) to the calf or thigh. It prevents "drag" on the catheter tubing, improving comfort and discretion for the user.

Product Features
- Compatible with **Prosys®, Libra® and Urilite® Leg Bags**
- Made of comfortable and skin-friendly material that allows the skin to breathe
- Cut-out at the bottom to accommodate the drainage tap
- Washable and reusable
- Available in **Small, Medium, Large** sizes

Selecting the Correct Size
To select the correct size of **Liberty® Leg Bag Sleeve**:
- Measure the **circumference (Ø) of your thigh or calf** (depending on where you intend to wear the sleeve).
- Select the appropriate size according to the product size table.

* All sizes are approximate and may vary dependent on user.

Fitting Instructions
1. Pull the **Liberty® Leg Bag Sleeve** over the foot like a sock.
2. Slide it up onto the **thigh or calf**, ensuring the **large slit at the top faces outwards.**
3. Insert your **Leg Bag** through the top slit.
4. Push the **drainage tap** through the pre-cut hole at the bottom of the sleeve.
5. Adjust the sleeve until it feels secure and comfortable.

Warnings & Precautions
- **Do not use** if the packaging or product is damaged. Inspect carefully before use.
- **Check the expiration date** on the product packaging. Do not use beyond this date.
- Follow local infection control procedures if the sleeve comes into contact with contaminated material.
- If you experience discomfort, skin irritation, or poor circulation, discontinue use and consult a healthcare professional.
- **Any serious incident** related to this device must be reported to the manufacturer and the competent authority.

Cleaning & Maintenance
- The **Liberty® Leg Bag Sleeve** is **machine washable.**
- Follow washing instructions provided on the care label.
- Ensure the sleeve is completely dry before reuse.

WASHING INSTRUCTIONS

Disposal
- Dispose of used or contaminated sleeves according to **local general waste regulations.**
- Devices contaminated with potentially infectious substances must be handled in accordance with infection control guidance.

REF	SIZE	Ø	QTY
PLS3881	Small/Klein/Petite	24-40cm	4
PLS3904	Medium/ Mittelgroß/Moyenne	36-50cm	4
PLS3928	Large/Groß/Grande	40-65cm	4

Non-sterile	Keep dry	Batch code	Catalogue number
Consult instructions for use	Keep away from sunlight	Date of manufacture	Medical device
Caution		Use by date	Unique device identifier
Not made with natural rubber latex	Do not use if package is damaged	Manufacturer	UK conformity mark
Authorized representative in the European community/European Union		European conformity mark	

Clinisupplies

Clinisupplies Ltd, 1 Blackmoor Lane, Croxley Park, Watford, Hertfordshire WD18 8GA, UK
Tel: +44 (0)20 8863 4168 Email: info@clinisupplies.co.uk **www.clinisupplies.co.uk**

EU REP Advena Ltd. Tower Business Centre, 2nd floor, Tower Street, Swatar, BKR 4013. Malta.

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