EN INSTRUCTIONS FOR USE

PLEASE READ CAREFULLY

1. Wash your hands before use
2. Please ensure you are using the correct size sheath – a measuring guide is enclosed for your reference
3. If necessary, trim any pubic hair or use the hair protectors supplied. Do not shave as this can cause skin irritation. Wash penis and groin area using a mild soap without any moisturiser and ensure the penis and the area around it is completely dry before proceeding.
4. The sheath is supplied individually in a case
5. The case has an opening clasp to allow the sheath to be removed easily
6. Unhook the opening clasp to unlock the case and to access the sheath **(Diagram 1)**
7. Remove the sheath from the pack by the spout (Diagram 2). Pop On style sheaths have a paper liner to protect the sheath in storage. Take the sheath off the paper liner. (Diagram 2)
8. Position the sheath over the penis
9. When applying the sheath leave a small space between the end of the penis and the base of the spout
10. Unroll the sheath along the shaft of the penis ensuring the end of the sheath is fully unrolled **(Diagram 3)**
11. Place your hand over the shaft of the penis and sheath and gently squeeze – this will allow the adhesive to warm up and ensures a comfortable fit **(Diagram 4)**
12. Fit drainage bag to leg first before connecting to the sheath to avoid any twisting as this may prevent drainage
13. Connect your drainage bag to the spout of the sheath making sure the tap on the bag is in the closed position **(Diagram 5)**

Removal Instructions

1. To remove the sheath, release the base of the sheath all the way around
2. Use warm soapy water and wash-cloth or use an Adhesive Remover Spray to help ease the sheath off **(Diagram 6)**
3. Dispose of your sheath and wash your hands
4. The **Livfria® Sheath** should be changed every 24 hours for hygienic purposes but this may vary between users. Advice can be sought from your healthcare professional
5. **Please note:** if your skin becomes irritated whilst using this product, do not use until you have spoken to your healthcare professional

Cautions

- Use of any medical device carries associated risks. Device failures are rare, however in the event of a device failure of this product, you should contact Clinisupplies GB Ltd using the details provided within this leaflet
- If a failure has resulted in a serious incident, such as significant deterioration of health or death, you should contact Clinisupplies GB Ltd and your local Competent Authority as soon as possible. Details of your local Competent Authority can be found on your local government website

Disposal

Sheath: should be disposed of in normal household waste

Case: Compostable material made from sugarcane pulp

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Clinisupplies GB Ltd. 1 Blackmoor Lane, Croxley Park, Watford, Hertfordshire, WD18 8GA
EU REP **MD Reps Ltd**
Clinisupplies Deutschland GmbH,