EN Libra® Leg Bag

INTENDED USE: Urine drainage leg bags are designed to collect urine from either sheaths or catheters.
Leg bag straps are intended for securing the leg bags to the calf/thigh of the leg.

Libra® Leg Bags are for the daytime storage and drainage of urine. **Libra®** Leg Bags can also be connected to a Night Bag for overnight drainage. **Libra®** Leg Bags should be changed every 5-7 days.

INSTRUCTIONS FOR USE

PLEASE READ CAREFULLY

Key

1. Protective Cap
2. Ridged Inlet Connector
3. Needle Free Sample Port
4. Inlet Tube
5. Eyelets
6. Graduated Scale*
7. Outlet Tap
8. Date Fitted Box
9. Silicone Sleeve

*The graduated scale is an approximation and should not be used for accurate measurement of urine output.

Before Use

1. Wash Hands
2. Open peel pouch and remove **Libra® Leg Bag** from its packaging
3. Ensure the lever on the outlet tap **(7)** is in the closed position **(Diagram 1)**.
4. Remove the grey protective cap **(1)** from the **Libra®** inlet connector **(2)**

Diagram 1
Close Tap

Diagram 2
Open Tap

Direct Connection

1. Insert the ridged inlet connector **(2)** into the end of the Indwelling Catheter or Sheath.
2. Insert the connector fully to allow a secure connection
3. The **Libra® Leg Bag** should be positioned on the thigh or calf, depending on the tube length, and secured appropriately with the leg bag straps provided.
4. Feed the **Libra® Leg Bag Straps** through the eyelets **(5)** on the bag and ensure the wavy silicone lines face onto the leg
5. Make sure the inlet tubing **(4)** isn't kinked which can restrict the flow of urine

Emptying

1. Empty the **Libra® Leg Bag** by pushing the lever on the outlet tap **(7)** downwards **(Diagram 2)**
2. After emptying, wipe the outlet tap with a tissue and pat a number of times before gently wiping around the inside of the valve. This will dislodge any residual urine
3. Remember to close the tap after emptying **(Diagram 1)**
4. The outlet tap **(7)** can be twisted to avoid accidental emptying

Overnight Drainage

If you require additional drainage for overnight use, a night bag can be connected as part of an overnight link system

1. When you go to bed, do not disconnect your **Libra® Leg Bag** from your Catheter or Sheath
2. Insert the connector of your Night Bag to the silicone sleeve **(9)** on your **Libra® Leg Bag**. Insert the connector fully to allow a secure connection
3. Once connected, turn the **Libra® Leg Bag** tap to the open position **(Diagram 2)** to allow the overnight drainage into the Night Bag
4. Before disconnecting the Night Bag in the morning, ensure the **Libra® Leg Bag** tap is in the closed position **(Diagram 1)**

Libra® Leg Bags with adjustable tubing

Adjustable **Libra® Leg Bag** have a 45cm tube **(4)**

If connecting directly to a Sheath, these can be cut to the desired length and fitted with the spare connector included in the pouch:

1. Place the **Libra® Leg Bag** on your leg in the desired position and secure as instructed above
2. Measure how long you need the tubing and cut to the desired length
3. Remove the **Libra® Leg Bag** from your leg and insert the spare connector into the tube ensuring it is secure
4. Connect the **Libra® Leg Bag** as per previously instructed

Caution: Do not cut the tube length yourself if you are connecting to an Indwelling Catheter as the sterility of the product may become compromised.

Storage and Disposal

1. The **Libra® Leg Bag** should be stored at an ambient temperature and allowed to return to room temperature before use
2. Dispose with normal household waste – do not flush down the toilet

Others

- Ensure you fill in the Date Fitted Box **(8)** to allow you to monitor how long the **Libra® Leg Bag** has been in use. Change the Leg Bag every 5-7 days.
- The **Libra® Leg Bag** should not be reused after disconnection from a catheter due to increased risk of infection
- Keep out of the reach of children
- Use of any medical device carries associated risks. Device failures are rare, however in the event of a device failure of this product, you should contact Clinisupplies using the details provided within this leaflet
- If a failure has resulted in a serious incident, such as significant deterioration of health or death, you should contact Clinisupplies and your local Competent Authority as soon as possible. Details of your local Competent Authority can be found on your local government website

Sterilization:

Leg bags are sterilized with EO gas. They are not suitable for resterilization

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