Libra® Cath Kit

Instructions For Use

This pack contains

1 x Removal Pack

1 x Insertion Pack

1 x **Liberty®** Gel Combi (12ml)

1 x **Libra®** All Silicone Foley Catheter

1 x Prefilled Syringe of Sterile Water for Balloon Inflation (10ml)


1 x **Libra®** Leg Bag

1 x **Libra®** Leg Bag Strap

1 x **Liberty®** Fix-it Catheter Retaining Strap – Long (80cm)

2 x Saline Solution (20ml)


EN Instructions for Use page 2
NL Gebruiksaanwijzing page 6
FR Mode d'emploi page 10
DE Gebrauchsanweisung page 14
DA Brugsanvisning page 18
FI Käyttöohje page 22
NO Brukerveiledning page 26
IT Istruzioni per l'uso page 30
PT Instruções de utilização page 34
ES Instrucciones de uso page 38
SV Bruksanvisning page 42

EN Libra® Cath Kit

Instructions for Use

Removal and Insertion Pack

Catheter Removal Pack
This pack contains all the necessary equipment to remove a catheter already in-situ.

Item	Quantity
Green Drape	1
Disposable Apron	1
Disposal Bag	1
Gloves	1 pair
Empty Syringe	1
Gauze Squares	3
Rectangular Gallipot	1

1. Wash hands
2. Open the 'Removal Pack' peel pouch and remove the green drape from the top
3. Open the green drape, taking care only to touch the green side, and place on a flat surface. The white side is your sterile field
4. Wash hands and apply the gloves and apron. Place the green modesty drape over the patient
5. Use the empty syringe provided to remove the amount of water instilled in the catheter balloon, and remove the catheter as per the catheter instructions and local clinical guidelines

Catheter Insertion Pack
This pack contains all the necessary equipment to insert a catheter.

Item	Quantity
Green Drape	1
Gloves	2 pairs
Fenestrated Drape	1
Disposable Apron	1
Disposal Bag	1
Round Gallipot	1
Non-Woven Balls	5
Gauze Squares	3
Rectangular Gallipot	1

1. Wash hands and open the 'Insertion Pack'. Remove the green drape from the top
2. Open the green drape, taking care only to touch the green side and place on a flat surface. The white side of the drape is your sterile field
3. Remove the following items from the **Libra® Cath Kit** Box:
 – **Libra® All Silicone Foley Catheter**
 – Prefilled syringe of sterile water
 – **Liberty® Gel Combi**
 – **Libra® Leg Bag**
 – Saline solution

 Do not place the **Libra® Leg Bag Strap** or **Libra® Fix-it Strap** onto the sterile field
4. Open each item into your sterile field, maintaining your sterile environment as per local clinical guidelines. Dispense the saline solution into the round gallipot
5. Wash hands and apply the gloves and apron. Place the fenestrated drape over the patient, exposing the catheter insertion site
6. Use saturated gauze balls to clean the catheter insertion site. Dry the patient with the gauze squares
7. Dispense the lubricating gel into the urethra, as per local clinical guidelines
8. Once the Lidocaine is active, wash hands and apply clean gloves
9. Connect the **Libra® Leg Bag** to the catheter funnel and proceed with catheterisation according to local clinical guidelines
10. To prevent trauma and/or discomfort, secure the catheter using the **Liberty® Fix-it Catheter Retaining Strap** provided
11. Secure the **Libra® Leg Bag** to the patient's leg using the **Libra® Leg Bag Straps**
12. Dispose of waste as per local clinical guidelines

Libra® All Silicone Foley Catheter

Product description
Libra® All Silicone Foley Catheter is a sterile, single-use product. The catheter has a drainage lumen and inflation channel. It has an integrated balloon towards the tip of the catheter and a funnel with an inflation valve at the bottom of the catheter. The capacity of the balloon can be found on the funnel and on the catheter pack label. **Libra® All Silicone Foley Catheter** is made of 100% silicone. The catheter has a nelaton (cylindrical) tip with two opposite drainage eyes. The valve of the catheter is colour coded and can be operated with a luer syringe tip. Each catheter is sterile and double-wrapped to aid in non-touch insertion.

Indications
Libra® All Silicone Foley Catheter is suitable for drainage of urine from the bladder.

Contra-indications
Libra® All Silicone Foley Catheter should not be used if the patient has a known allergic reaction to silicone.

Important
Catheterisation with a **Libra® All Silicone Foley Catheter** can only be performed by a suitably qualified and trained healthcare professional. The catheter is strictly single-use and must not be reused.

Precautions
- Recent urethral or bladder surgery
- Prior transurethral resection of the prostate with large tissue defect
- Significant symptoms of urinary obstruction prior to treatment
- History of abdominoperineal resection for rectal cancer, rectal stenosis or other major rectal pathology

Additional warnings and precautions
- Do not use ointments or lubricants with a petroleum base
- Do not use sharp instruments with the **Libra® All Silicone Foley Catheter** – sharp instruments could damage the catheter and cause it to malfunction
- It is recommended to drink more water when catheterised – please follow local clinical policies on fluid intake
- Instructions for use shall be read and abided by carefully
- Medical advice shall always be followed depending on patient's clinical symptoms
- Patients with indwelling catheters should be monitored in accordance with local and national clinical policies
- Strict adherence to aseptic technique and hand hygiene are essential to minimise the risk of Catheter Associated Urinary Tract Infections (CAUTI)

CSF00830_1_v1 LCK IFU Multi.indd 2 24/04/2025 15:55

Pre-Testing

Before inserting the catheter, use the prefilled syringe to inject the fluid volume into the inflation port, and then pull back on the syringe plunger to remove the fluid. This is to ensure the balloon has not been damaged.

Instructions for use

The catheter should be inserted using an aseptic technique

1. Lubricate the catheter using a water-based gel and proceed with catheterisation in line with the normal urethral catheterisation procedure
2. Insert the catheter into the meatus until urine starts to flow
3. Guide the catheter gently 5-8cm beyond the point at which urine begins to flow. The rationale for inserting the catheter further into the bladder ensures the catheter is beyond the neck of the bladder
4. Use the prefilled water syringe and inflate the balloon with the recommended volume into the injection port
5. Retract the catheter until you feel resistance
6. Secure the catheter with a catheter fixation device
7. Attach a drainage bag or catheter valve

The catheter is licensed to be left in-situ for up to 12 weeks. Regular monitoring of the patient should be undertaken to ensure there are no side effects and that the catheter is functioning correctly. Change the catheter according to local clinical policies. This should be no longer than 12 weeks.

Please note

- Only use water-based lubricants or gels
- Use a Luer Tip syringe to inflate the balloon with sterile water
- Fill the balloon slowly, taking care only to instil the amount of water specified on the catheter funnel and packaging

Possible complications

- Should balloon rupture occur, care should be taken to ensure no fragments are left inside the bladder
- In case the balloon does not deflate, please refer to your local Clinical Policies to resolve the issue. Do not cut the balloon inflation valve
- Use of any medical devices carries associated risks. Device failures are rare, however in the event of a device failure following the use of this product, you should contact Clinisupplies GB using the details provided within this leaflet
- If a failure has resulted in a serious incident, such as significant deterioration of health or death, you should contact Clinisupplies GB and your local Competent Authority as soon as possible. Details of your local Competent Authority can be found on your local government website

Sterilisation

The Libra® All Silicone Foley Catheter has been sterilised by Ethylene Oxide (EtO). The expiry date can be found on the sterile pouch. Do not use the catheter if the packaging has been damaged or accidentally opened.

Warning

Libra® All Silicone Foley Catheters are for single-use only. Do not re-sterilise or re-use.

MRI safety – This device is MR Conditional

Non-clinical testing demonstrated that this device is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla only
- Maximum spatial gradient magnetic field of 1,850-Gauss/cm (18.5-T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg (Normal Operating Mode).

Storage

Libra® All Silicone Foley Catheters should be stored in their original packaging in a cool, dry place away from direct sunlight.

Disposal

Libra® All Silicone Foley Catheters should be disposed of in line with your local Clinical Policies.

Liberty® Gel Combi

Liberty® Gel Combi is a sterile, water-soluble, non-irritating lubricant gel with supportive antiseptic and anaesthetic properties used to facilitate the insertion of catheters and other medical instruments (catheterisation, cystoscopy) into the urethra and bladder.

Liberty® Gel Combi helps prevent trauma by lubricating the interface between the urethral mucosa and catheter during catheterisation. The antiseptic property of chlorhexidine reduces the risk of infection and the local anesthetic property of lidocaine reduces perioperative pain associated with the procedure. It minimises patient discomfort during the catheterisation procedure.

Ingredients

100 g of **Liberty® Gel Combi** contains; Deionized Water, Hydroxyethyl Cellulose (Lubricative), Mono Propylene Glycol (Moisturizer), 2g Lidocaine (Local anesthetic), 0.05g Chlorhexidine Gluconate (Antiseptic), Methyl Hydroxybenzoate (Preservative), Propyl Hydroxybenzoate (Preservative).

Warnings and precautions for use

- Under anesthesia a lubricant without lidocaine should be used
- Instillation should only be performed by a physician or by qualified medical personnel
- Patients with glucose-6-phosphate dehydrogenase deficiency or congenital or idiopathic methaemoglobinaemia are more susceptible to drug induced methaemiglobinaemia
- When applying **Liberty® Gel Combi**, a temporary burning sensation may occur. Patients should be warned in this regard. In this case, cold compress to the applied region can alleviate the burning sensation
- Allergic reactions to lidocaine and/or chlorhexidine may rarely occur. Please consult your doctor if redness, burning and stinging occurs in the applied site
- Irritates eyes and/or skin in cases of direct contact. Wash with plenty of water in cases of exposure
- If swallowed, seek medical advice
- Keep away from children
- Store in dry and cool places
- Keep away from direct sunlight
- For single use only
- Do not re-sterilize
- Do not use products with damaged packages because of the risk of sterilisation deterioration

Contra-indications

Liberty® Gel Combi should not be used in patients with hypersensitivity to any of the components or patients with marked bradycardial rhythm disorders.

Liberty® Gel Combi should not be used in patients with severe cardiac insufficiency, severe hepatic and renal dysfunction, traumatic mucosa, inflammation, sepsis at the application site, and in patients with a tendency to convulsion (epilepsy, severe shocks). Do not use in children under 2 years.

Pregnancy and lactation

Liberty® Gel Combi can be used during pregnancy and lactation. The risk/benefit assessment should be performed by the physician. Lidocaine passes to placenta and enters breast milk in small quantities, however, therapeutic dose used in **Liberty® Gel Combi** does not generally pose a risk to the fetus and infants.

CSF00830_1_v1 LCK IFU Multi.indd 3 24/04/2025 15:55

Drug interaction

Interaction due to absorption of lidocaine may occur when used in combination with the following drugs;

- Propranolol: Reduction in plasma clearance of lidocaine
- Cimetidine: Reduction in plasma clearance of lidocaine
- Antiarrhythmic products: Increase in toxicity of lidocaine
- Phenytoin or barbiturates: Reduction in plasma levels of lidocaine

The specified interactions can be seen in long-term use and repeated high doses. No clinically significant interactions have been reported when administered at the recommended doses.

Due to the risk of interaction depending on chlorhexidine, use with the following is not recommended;

- Chlorhexidine Gluconate (CHX) + Sodium Hypochlorite (NaOCl)
- Chlorhexidine Gluconate (CHX) + Hydrogen Peroxide (H2O2)
- Chlorhexidine Gluconate (CHX) + Povidone-Iodine (BTD)

According to scientific literatures, CHX + NaOCl interaction is known to produce toxic by-products for humans.

Instructions for use

1. Wash your hands and prepare an area that is aseptic as possible to perform the operation. Prepare a clean serving tray and place required materials
2. Remove the syringe from its sterile package
3. Lightly press on the plunger before removing the blue cap to provide more even and gentle application. Remove cap
4. Carefully insert the cone of the syringe into the urethral opening and slowly inject the gel by applying gentle and even pressure
5. Apply the remaining gel to the surface of the catheter or to the surface of another medical device
6. **Liberty® Gel Combi** shows lubricative effect as soon as the application starts. The anaesthetic effect starts 3-5 minutes later
7. If more gel has to be used, condition of the applied tissue, age of the patient and any adverse condition of the patient which may be related to sensitivity of chlorhexidine must be taken into care by the physician
8. When the whole product is not used, the syringe and the remaining gel should be evaluated and discarded like other medical wastes
9. **Liberty® Gel Combi** should only be used under the supervision of specialist health personnel

Side effects

Allergic reactions to lidocaine and chlorhexidine are extremely rare. In the case of hypersensitivity, skin reactions are possible. Skin reactions depending on anaphylaxis (hypersensitivity) may occur due to the use of lidocaine in the damaged mucosa. Systemic adverse reactions to lidocaine are rare and may result from hypersensitivity or diminished tolerance.

Sterilisation

Liberty® Gel Combi is sterilised with radiation after packaging. In the case that the packaging is damaged, do not use as sterility may be compromised.

Cautions

Use of any medical device carries associated risks. Device failures are rare, however in the event of a device failure, you should contact Clinisupplies GB using the details provided in this leaflet.

If a failure has resulted in a serious incident, such as significant deterioration of health or death, you should contact Clinisupplies GB and your local Competent Authority as soon as possible. Details of your Competent Authority can be found on your local government website.

EasyPod Saline

Instructions for use

Sterile saline pods for topical irrigation and wound cleansing. (Sodium Chloride 0.9% w/v)

1. Hold pod upright and open it by twisting off the top
2. For directional irrigation position the pod and squeeze firmly, taking care not to touch the nozzle. Discard any unused solution

Additional warnings and precautions

- For irrigation only
- Not for injection
- Do not use unless the solution is clear and the unit is intact


Liberty® Fix-it Retaining Strap

The **Liberty® Fix-it Catheter Retaining Strap** is designed to secure indwelling and suprapubic catheters.

Selecting the appropriate strap length

Use a measuring tape to measure the circumference of the thigh at the widest point, or around the abdomen at the area the strap is to be worn. Choose a Fix-it Strap that is longer than the measured circumference to ensure proper fit and secure placement.


Instructions for use

Remove **Liberty® Fix-it** and place in the desired position with the two silicone lines facing the skin. Use the hook and eye fastening to fix into place – ensuring a comfortable and secure grip. Cut strap to size if necessary

At the tapered end lift the **Liberty® Fix-it** tab and place the tube on to the strap. Once in position thread the tapered end through the eyelet

Gently pull the **Liberty® Fix-it** tab to tighten and fasten the tapered end to the fabric using the velcro strip. Ensure there is no restriction to blood flow

Washing instructions

For best results, place the **Liberty® Fix-it Strap** in a garment bag or pillow case and wash at 30°C on a delicate setting.

CSF00830_1_v1 LCK IFU Multi.indd 4 24/04/2025 15:55


Libra® Leg Bag

Libra® Leg Bags are for the daytime storage and drainage of urine

Designed to be connected directly to a catheter for continuous drainage


Libra® Leg Bags can also be connected to a night bag for overnight drainage (See 'Overnight Drainage')


Libra® Leg Bags should be changed every 5-7 days

Before Use

1. Wash hands
2. Open peel pouch and remove the **Libra® Leg Bag** from its packaging
3. Ensure the lever on the outlet tap **(7)** is in the closed position **(Diagram 1)**
4. Remove the grey protective cap **(1)** from the **Libra®** inlet connector **(2)**

Direct Connection

1. Insert the ridged inlet connector **(2)** into the end of the Indwelling Catheter. Insert the connector fully to allow a secure connection. To reduce the risk of introducing bacteria, do not touch connector whilst connecting
2. The **Libra® Leg Bag** should be positioned on the thigh or calf, depending on the tube length, and secured appropriately with the leg bag straps provided
3. Feed the leg bag straps through the eyelets **(5)** on the bag and ensure the wavy silicone lines face onto the leg
4. Make sure the inlet tubing **(4)** isn't kinked which can restrict the flow of urine

Emptying

1. Empty the **Libra® Leg Bag** by pushing the lever on the outlet tap **(7)** downwards **(Diagram 2)**
2. After emptying, wrap the outlet tap with a tissue and pat a number of times before gently wiping around the inside of the valve, this will dislodge any residual urine
3. Close the tap after emptying **(Diagram 1)**
4. The outlet tap **(7)** can be twisted to avoid accidental emptying

Overnight Drainage

1. If you require additional drainage for overnight use, a night bag can be connected as part of an overnight link system. When you go to bed, do not disconnect your leg bag from your catheter
2. Insert the connector of your night bag to the silicone sleeve **(9)** on your leg bag. Insert the connector fully to allow a secure connection
3. Once connected, turn the **Libra® Leg Bag** tap to the open position **(Diagram 2)** to allow the overnight drainage into the night bag
4. Before disconnecting the night bag in the morning, ensure the **Libra® Leg Bag** tap is in the closed position **(Diagram 1)**

Storage and Disposal

The **Libra® Leg Bag** should be stored at an ambient temperature and allowed to come to room temperature before use

Dispose with normal household waste – do not flush down the toilet

Cautions

Ensure you fill in the Data Fitted Box **(8)** to allow you to monitor how long the **Libra® Leg Bag** has been in use. Change the **Libra® Leg Bag** every 5-7 days.

The **Libra® Leg Bag** should not be reused after disconnection from a catheter due to increased risk of infection.

Keep out of reach of children.

Use of any medical device carries associated risks. Device failures are rare, however in the event of a device failure, you should contact Clinisupplies GB using the details provided in this leaflet.


If a failure has resulted in a serious incident, such as significant deterioration of health or death, you should contact Clinisupplies GB and your local Competent Authority as soon as possible. Details of your Competent Authority can be found on your local government website.

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CSF00830_1_v1 LCK IFU Multi.indd 45 24/04/2025 15:55

Manufacturer's details

Component	Manufacturer's details	Applicable labelling
Libra® Leg Bag	**Clinisupplies GB Ltd.** 1 Blackmoor Lane, Croxley Park, Watford, Hertfordshire, WD18 8GA, UK +44 (0) 20 8863 4168 EC REP **MD Reps Ltd** 12-16 Bridge Street, Belfast, Northern Ireland, BT1 1LU	
Libra® Leg Bag Strap		
Liberty® Fix-it Catheter Retaining Strap		
Liberty® Gel Combi	**Turkuaz Medikal** Akçaburgaz Mah. Muhsin Yazıcıoğlu Cad No:45/5 34522 Esenyurt, Istanbul/Türkiye +90 (212) 428 68 48	
Libra® All Silicone Foley Catheter	**Fortune Medical Instrument Corp.** 6FL, No.29, Sec.2, Jhongjheng E. Rd, Danshuei Dist., New Taipei City, 251, Taiwan +(00)886 2 2624 2233 EC REP **Emergo Europe** Prinsessegracht 20, 2514 AP The Hague, The Netherlands	
Removal Pack	**Rocialle** Cwm Cynon Business Park (North), Mountain Ash, Rhondda Cynon Taff, South Wales, CF45 4ER, UK. +44 (0)1433 471 300 EC REP **Advena Limited** Lower Business Centre, 2nd Flr, Tower Street, Swatar, BKR 4013, Malta	
Insertion Pack		
Prefilled Water Syringe	**Well Lead Medical Co. Ltd** C-4 Jinhu Industrial Estate Hualong, 511434 Panyu, Guangzhou, People's Republic of China. +(00)86 20 8475 8878 EC REP **Shanghai International Holding Corp.** GmmH (Europe), Eiffestrasse 80, 20537 Hamburg, GERMANY	
EasyPod Saline	**Amanta Healthcare Ltd.** Plot no: 876, N. H. no. 8, Hariyala, Taluka: Matar, Dist: Kheda-387411 Gujarat, India. +91 2694 287700 EC REP **CMC Medical devices & Drugs S.L.,** Spain	

CSF00830_1_v1 LCK IFU Multi.indd 46 24/04/2025 15:55

Symbol explanation

EN: Do not use if package is damaged
NL: Niet gebruiken indien de verpakking beschadigd is
FR: Ne pas utiliser si l'emballage est endommagé
DE: Nicht verwenden, wenn die Verpackung beschädigt ist
DA: Må ikke anvendes, hvis emballagen er beskadiget
FI: Älä käytä, jos pakkaus on vaurioitunut
IT: Non utilizzare se la confezione è danneggiata
NO: Skal ikke brukes hvis pakken er skadet
PT: Não utilizar se a embalagem estiver danificada
ES: No usar si el envase presenta algún daño
SV: Får inte användas om förpackningen är skadad

EN: Consult instructions for use
NL: Raadpleeg instructies voor gebruik
FR: Lire la notice avant utilisation
DE: Gebrauchsanweisung beachten
DA: Se brugsanvisningen før brug
FI: Lue käyttöohjeet
IT: Consultare le istruzioni per l'uso
NO: Se bruksanvisningen
PT: Consultar as instruções de utilização
ES: Consultar las instrucciones de uso
SV: Läs bruksanvisningen

EN: Box and instructions for use are recyclable
NL: Doos en gebruiksinstructies zijn recyclebaar
FR: La boîte et le mode d'emploi sont recyclables
DE: Schachtel und Gebrauchsanweisung sind recycelbar
DA: Æske og brugsanvisning kan genbruges
FI: Laatikko ja käyttöohjeet ovat kierrätettäviä
IT: La scatola e le istruzioni per l'uso sono riciclabili
NO: Eske og bruksanvisning er resirkulerbare
PT: A caixa e as instruções de utilização são recicláveis
ES: La caja y las instrucciones de uso son reciclables
SV: Förpackning och bruksanvisning är återvinningsbara

EN: Use-by-date
NL: Te gebruiken tot
FR: Date limite d'utilisation
DE: Verfallsdatum
DA: Sidste anvendelsesdato
FI: Viimeinen käyttöpäivä
IT: Da utilizzarsi entro
NO: Holdbarhetsdato
PT: Data de validade
ES: Fecha de caducidad
SV: Utgångsdatum

EN: Keep away from sunlight
NL: Beschermen tegen zonlicht
FR: À conserver à l'abri de la lumière directe du soleil
DE: Von Sonnenlicht fernhalten
DA: Holdes væk fra sollys
FI: Pidä poissa suorasta auringonvalosta
IT: Tenere lontano dalla luce del sole
NO: Må ikke utsettes for sollys
PT: Manter afastado da luz solar
ES: Mantener fuera de la luz del sol
SV: Förvara skyddad från direkt solljus

UK CA
EN: UK conformity mark
NL: UK conformiteitsmerkteken
FR: Marque de conformité brittanique
DE: GB-Konformitätszeichen
DA: UK-overensstemmelsesmærke
FI: Yhdistyneen kuningaskunnan vaatimustenmukaisuusmerkintä
IT: Contrassegno di conformità Regno Unito
NO: Britisk samsvarsmerke
PT: Marca de conformidade do Reino Unido
ES: Marca de conformidad del Reino Unido
SV: Brittiskt överensstämmelsemärke

EN: Importer
NL: Importeur
FR: Importateur
DE: Importeur
DA: Importør
FI: Maahantuoja
IT: Importatore
NO: Importør
PT: Importador
ES: Importador
SV: Importör

EN: Distributor
NL: Distributeur
FR: Distributeur
DE: Lieferant
DA: Distributør
FI: Jakelija
IT: Distributore
NO: Distributør
PT: Distribuidor
ES: Distribuidor
SV: Distributör

EN: Caution
NL: Let op
FR: Attention
DE: Achtung
DA: Forsigtig
FI: Varoitus
IT: Attenzione
NO: Forsiktig
PT: Cuidado
ES: Advertencia
SV: Försiktighet

EN: Double sterile barrier system
NL: Dubbel steriel barrièresysteem
FR: Système de barrière stérile double
DE: Doppeltes Sterilbarrieresystem
DA: Dobbelt sterilt barrieresystem
FI: Kaksinkertainen steriili suojajärjestelmä
IT: Sistema a barriera doppia sterile
NO: Dobbelt sterilt barrieresystem
PT: Sistema de barreira estéril duplo
ES: Sistema doble de barrera estéril
SV: System med dubbel steril barriär

STERILE EO
EN: Sterilised Ethylene oxide
NL: Gesteriliseerd met ethyleenoxide
FR: Stérilisé à l'oxyde d'éthylène
DE: Sterilisiert mit Ethylenoxid
DA: Steriliseret ved hjælp af bestråling
FI: Steriloitu etyleenioksidilla
IT: Sterilizzato con ossido di etilene
NO: Sterilisert ved hjelp av etylenoksid
PT: Esterilizado com óxido de etileno
ES: Esterilizado usando óxido de etileno
SV: Steriliserad med etylenoxid

REF
EN: Catalogue number
NL: Catalogusnummer
FR: Numéro de catalogue
DE: Katalognummer
DA: Katalognummer
FI: Luettelonumero
IT: Numero di catalogo
NO: Katalognummer
PT: Número do catálogo
ES: Número de catálogo
SV: Katalognummer

EN: Not made with natural rubber latex
NL: Dit apparaat bevat geen natuurlatex
FR: Ce dispositif ne contient pas latex naturel
DE: Die Einmalkatheter sind frei von Naturkautschuklatex
DA: Denne enhed er fri for naturlig latex
FI: Tämä laite ei sisällä luonnonkumilateksia
IT: Dispositivo privo di lattice naturale
NO: Denne enheten er fri for naturlig lateks
PT: Este dispositivo não contém látex natural
ES: Este producto no contiene látex natural

EN: Do not re-use
NL: Niet hergebruiken
FR: Ne pas réutiliser
DE: Nicht wiederverwenden
DA: Må ikke genbruges
FI: Älä käytä uudelleen
IT: Non riutilizzare
NO: Må ikke brukes om igjen
PT: Não reutilizar
ES: No reutilizar
SV: Får inte återanvändas

EN: Temperature limits
NL: Temperatuurgrenzwerte
FR: Limites de température
DE: Temperaturgrenzwerte
DA: Temperaturgrænser
FI: Lämpötilan raja-arvot
IT: Limiti di temperatura
NO Temperaturgrenser
PT: Limites de temperatura
ES: Límites de temperatura
SV: Temperaturgränser

EN: Keep dry
NL: Droog bewaren
FR: À conserver au sec
DE: Trocken aufbewahren
DA: Holdes tør
FI: Pidä kuivana
IT: Tenere asciutto
NO: Må holdes tørt
PT: Manter seco
ES: Mantener seco
SV: Förvaras torrt

LOT
EN: Batch code
NL: Batchcode
FR: Numéro de lot
DE: Chargencode
IT: Lotto di produzione
NO: Batch-kode
DA: Batchkode
FI: Sarjanumero
PT: Código do lote
ES: Código de lote
SV: Satskod

EN: Do not use sharp instruments
NL: Gebruik geen scherpe instrumenten
FR: Ne pas utiliser d'instruments pointus
DE: Verwenden Sie keine scharfen Gegenstände
DA: Brug ikke skarpe instrumenter
FI: Älä käytä teräviä instrumentteja
IT: Non utilizzare strumenti affilati
NO: Ikke bruk skarpe instrumenter
PT: Não utilize instrumentos pontiagudos
ES: No usar instrumentos afilados
SV: Använd inte vassa instrument

STERILE R
EN: Sterilized using irradiation
NL: Gesteriliseerd door middel van bestraling
FR: Stérilisé par irradiation
DE: Sterilisiert durch Bestrahlung
DA: Steriliseret ved hjælp af bestråling
FI: Steriloitu säteilyttämällä
IT: Sterilizzato tramite irradiazione
NO: Sterilisert ved hjelp av bestråling
PT: Esterilizado por irradiação
ES: Esterilizado usando irradiación
SV: Steriliserad med strålning

EN: Single sterile barrier system
NL: Enkel steriel barrièresysteem
FR: Système de barrière stérile unique
DE: Einfaches Sterilbarrieresystem
DA: Enkelt sterilt barrieresystem
FI: Yksinkertainen steriili suojajärjestelmä
IT: Sistema a barriera sterile singola
NO: Enkelt sterilt barrieresystem
PT: Sistema de barreira estéril único
ES: Sistema de protección estéril único
SV: System med en steril barriär

MR
EN: MR conditional
NL: MR voorwaardelijk
FR: MR conditionnel
DE: Bedingt MR-sicher
DA: MR-kompatibel
FI: MR-ehdollinen
IT: MR condizionale
NO: MR-betinget
PT: Condicionado para RM
ES: Condicional para RM
SV: MR Conditional

EN: Wash hands
NL: Hände waschen
FR: Lavez-vous les mains
DE: Hände waschen
DA: Vask hænder
FI: Pese kädet
IT: Lavarsi le mani
NO Vask hendene
PT: Lave as mãos
ES: Lávese las manos
SV: Tvätta händerna

EN: Do not resterilize
NL: Niet opnieuw sterilise
FR: Ne pas restériliser
DE: Nicht erneut Sterilisieren
DA: Må ikke resteriliseres
FI: Älä steriloi uudelleen
IT: Non sterilizzare di nuovo
NO: Må ikke resteriliseres
PT: Não reesterilizar
ES: No reesterilizar
SV: Får inte omsteriliseras ren

EN: Date of manufacture
NL: Datum van fabricage
FR: Date de fabrication
DE: Herstellungsdatum
DA: Fremstillingsdato
FI: Valmistuspäivämäärä
IT: Data di produzione
NO: Produksjonsdato
PT: Data de fabrico
ES: Fecha de facturación
SV: Tillverkningsdatum

EN: Manufacturer
NL: Fabrikant
FR: Fabricant
DE: Hersteller
DA: Producent
FI: Valmistaja
IT: Produttore
NO: Produsent
PT: Fabricante
ES: Fabricante
SV: Tillverkare

MD
EN: Medical device
NL: Medisch hulpmiddel
FR: Appareil médical
DE: Medizinisches Gerät
DA: Medicinsk enhed
FI: Lääketieteellinen laite
IT: Dispositivo medico
NO: Medisinsk enhet
PT: Dispositivo médico
ES: Dispositivo médico
SV: Medicintekniskt produkt

UDI
EN: Unique device identifier
NL: Unieke code voor hulpmiddelidentificatie
FR: Identifiant d'appareil unique
DE: Eindeutige Gerätekennung
DA: Entydig enhedsidentifikator
FI: Ainutlaatuinen laitetunniste
IT: Identificatore unico del dispositivo
NO: Unik enhetsidentifikator
ES: Identificador único del dispositivo
PT: Identificador único do dispositivo
SV: Unik produktidentifiering (UDI)

CE
EN: European conformity mark
NL: Europees conformiteitsteken
FR: Marque de conformité européenne
DE: Europäisches Konformitätszeichen
DA: Europæisk overensstemmelsesmærke
FI: Eurooppalainen vaatimustenmukaisuusmerkki
IT: Contrassegno di conformità europeo
NO: Europeisk samsvarsmerke
PT: Marca europeia de conformidade
ES: Marca de conformidad europea
SV: Europeiskt överensstämmelsemärke

EC REP
EN: Authorized representative in the European community/ European union
NL: Gemachtigd vertegenwoordiger in de Europese Gemeenschap/Europese Unie
FR: Mandataire agréé dans la communauté européenne/l'union européenne
DE: Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union
DA: Bemyndiget repræsentant i Det Europæiske Fællesskab/ Den Europæiske Union
FI: Valtuutettu edustaja Euroopan yhteisössä/Euroopan unionissa
IT: Rappresentante autorizzato nella comunità Europea/ Unione Europea
NO: Autorisert representant i EU
PT: Representante autorizado na comunidade Europeia/União Europeia
ES: Representante autorizado en la Comunidad Europea/Unión Europea
SV: Auktoriserad representant i Europeiska gemenskapen/ Europeiska unionen

STERILE A
EN: Sterilised using aseptic processing techniques
NL: Gesteriliseerd met aseptische verwerkingstechnieken
FR: Stérilisé par techniques de traitements aseptiques
DE: Sterilisiert durch aseptische Aufbereitungstechniken
DA: Steriliseret ved hjælp af aseptiske behandlingsteknikker
FI: Steriloitu aseptisia käsittelytekniikoita käyttäen
IT: Sterilizzato con tecniche di trattamento asettiche
NO: Sterilisert ved hjelp av aseptiske behandlingsteknikker
PT: Esterilizado através de técnicas de processamento asséptico
ES: Esterilizado usando técnicas de procesamiento asépticas
SV: Steriliserad med aseptisk teknik

EN: Single sterile barrier system with protective packaging inside
NL: Enkel steriel barrièresysteem met beschermende verpakking binnenin
FR: Système de barrière stérile simple avec emballage protecteur interne
DE: Einzelnes steriles Barrieresystem mit Schutzverpackung im Inneren
DA: Enkelt sterilt barrieresystem med beskyttende indvendig emballage
FI: Yksinkertainen steriili suojajärjestelmä, jonka sisällä on suojapakkaus
IT: Sistema a barriera sterile singola con involucro protettivo esterno
NO: Enkelt sterilt barrieresystem med beskyttende indre emballasje
PT: Sistema a barriera singola sterile con confezione protettiva interna
ES: Sistema sencillo de barrera estéril con embalaje protector en el interior
SV: System med en steril barriär med skyddsförpackning på insidan

CSF00830_1_v1 LCK IFU Multi.indd 47 24/04/2025 15:55

Scan the QR code
to see our complete guide to catheter removal and insertion

Clinisupplies

Assembled by: **Clinisupplies GB Ltd.**
1 Blackmoor Lane, Croxley Park, Watford, Hertfordshire, WD18 8GA, UK
Tel: +44 (0) 20 8863 4168 Email: info@clinisupplies.co.uk **www.clinisupplies.co.uk**

EC REP **MD Reps Ltd**
12-16 Bridge Street, Belfast, Northern Ireland, BT1 1LU

Clinisupplies GB Ltd is formerly Great Bear Healthcare Ltd.,
1 Lambourne Crescent, Cardiff Business Park,
Cardiff, CF14 5GF, UK.

CSF00830_1_v1
02/25

CSF00830_1_v1 LCK IFU Multi.indd 48 24/04/2025 15:55