


# INSTRUMENT REPROCESSING CYCLE

From the Point-of-Use



Remove Visible Gross Soil



Transport Instruments to Decontamination Area




Thoroughly Clean Instruments




Remove Protective Equipment and Perform Hand Hygiene

Transport Instruments to Clean Area




Dry Instruments and Follow Procedural Instructions



Function Test, Inspect, and Assemble



Package Instruments




Label and Attach External Indicator


Transport Instruments to Sterilization Area



Document Using Log Sheet



Sterilize Instruments



Cool Down and Inspect



Store Sterile Instruments

Transport Sterilized Instruments to the Point-of-Use

For more information visit [www.spectrust.org](http://www.spectrust.org)  
or get in touch at [hello@spectrust.org](mailto:hello@spectrust.org)

© SPECT 2021

