



Total Solution for Sinus Graft Surgery

Collagen


Membrane

- Highly pure type 1 collagen
- Thin membrane (300µm) with multiple layers
- 6 months resorption period provides enough time for stabilizing graft materials and supporting bone to grow.
- Multiple-layered structure enables more effective bone regeneration by sparing enough space for hard tissue formation and facilitates proliferation of osteoblasts.
- Easy manipulation
- Dual-sided usage
- Size: 10x20, 15x20, 20x30, 30x40 (mm)
- Thickness: 0.3mm

Rabbit Calvaria Model - 6-12 weeks




Biodegradable barrier membrane for guided bone / tissue regeneration



• Clinical Evidence



Implantation




Collagen Membrane



4 months healing



Final prosthesis



SuperLine 4510



OSTEON™II



Collagen Membrane



5 months healing



DentiumUSA
Developed by Clinicians for Clinicians

6731 Katella Ave. Cypress, CA 90630
Tel. 1-877-304-6752 Fax. 1-714-226-0019
cs@dentiumusa.com



PN3041-1