

This is an Appendix to the NGPod System Instructions for use (IFU). Please refer to the full IFU before operating the device.

Prior to use of the NGPOD system, compatibility with the nasogastric tube(s) (NGT) to be used should be established.

The table below lists NGPOD sensor code to be used with NGTs that NGPOD has been validated with.

If the nasogastric tube you wish to use is not included in the below list, please contact NGPOD Global www.ngpodglobal.com, info@ngpodglobal.com for advice or to arrange a validation.

Nasogastric Tube Compatibility Reference Table

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Avanos/Halyard	CORFLO	42-8226	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	56	6	NGPOD SENSOR 01-3L9	5060644021255
Avanos/Halyard	CORFLO	40-7438	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	109	8	NGPOD SENSOR 01-6LE	5060644020487
Avanos/Halyard	CORFLO	42-7368	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	91	8	NGPOD SENSOR 01-5LC	5060644020302
Avanos/Halyard	CORFLO	42-9368	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	91	8	NGPOD SENSOR 01-5LC	5060644020302

Avanos/Halyard	CORFLO	42-7361	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	91	10	NGPOD SENSOR 01-5LC	5060644020302
Avanos/Halyard	CORFLO	42-9361	Nasogastric/Nasointestinal Feeding tube with Stylet	91	10	NGPOD SENSOR 01-5LC	5060644020302
MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Bexen Medical		052.12	Nasogastric Tube for Nutrition and Aspiration	110	14	NGPOD SENSOR 01-7L1	5060644020517

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
B. Braun	Nutritub Gastral	9246550	Gastric Polyurethane tube	100	15	NGPOD SENSOR 01-5LC	5060644020302
MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
CAIR L. G. L	NUTRICAIR™ ENTERAL-ENFIT®	NCE106G	Adult tube with ENFit® + Weight and Guide, in PU	110	6	NGPOD SENSOR 01-6L7	5060644020418
CAIR L. G. L	NUTRICAIR™ ENTERAL-ENFIT®	NCE108G	Adult tube with ENFit® + Weight and Guide, in PU	110	8	NGPOD SENSOR 01-6L7	5060644020418
CAIR L. G. L	NUTRICAIR™ ENTERAL-ENFIT®	NCE110G	Adult tube with ENFit® + Weight and Guide, in PU	110	10	NGPOD SENSOR 01-6L7	5060644020418

CAIR L. G. L	NUTRICAIR™ ENTERAL-ENFIT®	NCE112G	Adult tube with ENFit® + Weight and Guide, in PU	110	12	NGPOD SENSOR 01-6L7	5060644020418
CAIR L. G. L	NUTRICAIR™ ENTERAL-ENFIT®	NCE208G	Adult tube with ENFit® + Weight and Guide, in PU	110	8	NGPOD SENSOR 01-6L8	5060644020425
CAIR L. G. L	NUTRICAIR™ ENTERAL-ENFIT®	NCE210G	Adult tube with ENFit® + Weight and Guide, in PU	110	10	NGPOD SENSOR 01-6L7	5060644020418
CAIR L. G. L	NUTRICAIR™ ENTERAL-ENFIT®	NCE212G	Adult tube with ENFit® + Weight and Guide, in PU	110	12	NGPOD SENSOR 01-6L7	5060644020418
CAIR L. G. L	Cair LGL	AL114G	Nasogastric tube, Polyurethane with guidewire	110	14	NGPOD SENSOR 01-6L9	5060644020340
CAIR L. G. L	Cair LGL	AL116G	Nasogastric tube, Polyurethane with guidewire	110	16	NGPOD SENSOR 01-6L8	5060644020425

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
COMPAT	Compat Soft Junior	12362546	Enteral feeding tube	90	6	NGPOD SENSOR 01-5L5	5060644020234
COMPAT	Compat Soft	12363210	Enteral feeding tube	120	8	NGPOD SENSOR 01-7L2	5060644020524
COMPAT	Compat Soft	12362842	Enteral feeding tube	120	10	NGPOD SENSOR 01-7L2	5060644020524
COMPAT	Compat Soft	12362605	Enteral feeding tube	120	12	NGPOD SENSOR 01-7L2	5060644020524
COMPAT	Compat Soft	12363179	Enteral feeding tube	120	14	NGPOD SENSOR 01-7L2	5060644020524

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
COVIDIEN™	Salem Sump™	266122	Stomach Tube dual lumen with anti-reflux valve	Unspecified	14	NGPOD SENSOR 01-6LD	5060644020470
COVIDIEN™	Salem Sump™	266130	Stomach Tube dual lumen with anti-reflux valve	Unspecified	16	NGPOD SENSOR 01-6LB	5060644020456
COVIDIEN™	Salem Sump™	266148	Stomach Tube dual lumen with anti-reflux valve	Unspecified	18	NGPOD SENSOR 01-6LC	5060644020463
COVIDIEN™	Salem Sump™	264945	Stomach Tube – dual lumen	Unspecified	14	NGPOD SENSOR 01-6LB	5060644020456
COVIDIEN™	Salem Sump™	264960	Stomach Tube -dual lumen	Unspecified	16	NGPOD SENSOR 01-6LB	5060644020456
COVIDIEN™	Salem Sump™	1180264531	Polyurethane Gastroduodenal tube	114	16	NGPOD SENSOR 01-7L4	5060644020548
COVIDIEN™	Salem Sump™	1180264432	PVC Gastroduodenal tube	114	16	NGPOD SENSOR 01-7L2	506044020524

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Fresenius Kabi	Freka Tube	7981839	Nasogastric tube with ENfit and guide wire	80cm	8	NGPOD SENSOR 01-4LD	5060644020555
Fresenius Kabi	Freka Tube	7981840	Nasogastric tube with ENfit and guide wire	120cm	8	NGPOD SENSOR 01-7L5	5060644020555
Fresenius Kabi	Freka Tube	7981842	Nasogastric tube with ENfit and guide wire	120cm	10	NGPOD SENSOR 01-7L5	5060644020555

Fresenius Kabi	Freka Tube	7981843	Nasogastric tube with ENfit and guide wire	120cm	12	NGPOD SENSOR 01-7L7	5060644020579

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
GB UK Enteral	Nutricare	NC06.050ISOSAF	Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	50	6	NGPOD SENSOR 01-3L3	5060644021194
GB UK Enteral	Nutricare	NC08.050ISOSAF	Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	50	8	NGPOD SENSOR 01-3L3	5060644021194
GB UK Enteral	Nutricare	NC08.092ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	92	8	NGPOD SENSOR 01-092	5060644020005
GB UK Enteral	Nutricare	NC10.092ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	92	10	NGPOD SENSOR 01-092	5060644020005
GB UK Enteral	Nutricare	NC12.092ISOSAF/ NC12.092ISOSAF. ROW	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	92	12	NGPOD SENSOR 01-092	5060644020005
GB UK Enteral	Nutricare	NC08.080ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	80	8	NGPOD SENSOR 01-5L1	5050644020197

GB UK Enteral	Nutricare	NC06.080ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	80	6	NGPOD SENSOR 01- 5L1	5050644020197
GB UK Enteral	Nutricare	NC12-110ISOSAF-ROW	Nasogastric feeding tube, polyurethane with ENFIT connector	110	12	NGPOD SENSOR 01- 6LF	5060644020494
GB UK Enteral	Nutricare	NC08.120ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guide wire	120	8	NGPOD SENSOR 01- 7L9	5060644020609
GB UK Enteral	Nutricare	NC08.092ISOSAF.ROW	Nasogastric feeding tube, polyurethane	92	8	NGPOD SENSOR 01- 092	5060644020005
GB UK Enteral	Nutricare	NC14.092ISOSAF.ROW	Nasogastric feeding tube, polyurethane	92	14	NGPOD SENSOR 01- 5L7	5060644020258
GB UK Enteral	Nutricare	NC08.120ISOSAF.ROW	Nasogastric feeding tube, with guidewire	120	8	NGPOD SENSOR 01- 7L9	5060644020609
GB UK Enteral	Nutricare Infant Clear	NCI06.075ISOSAF	Nasogastric feeding tube, polyurethane with ENFIT connector	75	6	NGPOD SENSOR-01- 4LA	5060644021125
GB UK Enteral	Nutricare Infant Clear	NCI10.075ISOSAF	Nasogastric feeding tube, polyurethane with ENFIT connector	75	10	NGPOD SENS0-01- 4LA	5060644021125
GB UK Enteral	Nutricare Infant Clear	NCI10.075ISOSAF	Nasogastric feeding tube, polyurethane, with Enfit connector	75	10	NGPOD SENSOR 01- 4LA	5060644021125
GB UK Enteral	Carefeed	CF01.10.092ISOSAF	Feeding tube with ENFit connector	92	10	NGPOD SENSOR 01- 5LA	5060644020289
GB UK Enteral	Carefeed	CF01.12.092ISOSAF	Feeding tube with ENFit connector	92	12	NGPOD SENSOR 01- 5L9	5060644020272
GB UK Enteral	Carefeed	CF01.14.092ISOSAF	Feeding tube with ENFit connector	92	14	NGPOD SENSOR 01- 5L9	5060644020272

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
HMC Premedical S.p.A	Medicina®	NGP6/55L	Enfit Nasogastric feeding tube with guidewire	55	6	NGPOD SENSO-01-3L6	5060644021224
HMC Premedical S.p.A	Medicina®	NGP8/55L	Enfit Nasogastric feeding tube with guidewire	55	8	NGPOD SENSO-01-3L6	5060644021224
HMC Premedical S.p.A	Medicina®	LG6/80	Enfit Nasogastric feeding tube	80	6	NGPOD SENSO-01-4LF	5060644021064
HMC Premedical S.p.A	Medicina®	LG8/80	Enfit nasogastric feeding tube	80	8	NGPOD SENSOR 01-4LE	5060644021163
HMC Premedical S.p.A	Medicina®	LG12/80	Enfit nasogastric feeding tube	80	12	NGPOD SENSOR 01-4LE	5060644021163
HMC Premedical S.p.A	Medicina®	NGP8/85L	PUR Feeding tube	85	8	NGPOD SENSOR 01-5L2	5060644020104
HMC Premedical S.p.A	Medicina®	NGP10/85L	PUR Feeding tube	85	10	NGPOD SENSOR 01-5L5	5060644020234
HMC Premedical S.p.A	Medicina®	RT10/100L	ENFit Short Term Feeding/Aspiration Tube	100	10	NGPOD SENSOR 01-6L0	5060644020340
HMC Premedical S.p.A	Medicina®	RT12/100L	ENFit Short Term Feeding/Aspiration Tube	100	12	NGPOD SENSOR 01-6L0	5060644020340
HMC Premedical S.p.A	Medicina®	RT14/100L	ENFit Short Term Feeding/Aspiration Tube	100	14	NGPOD SENSOR 01-6L0	5060644020340
HMC Premedical S.p.A	Medicina®	RT16/100L	ENFit Short Term Feeding/Aspiration Tube	100	16	NGPOD SENSOR 01-6L0	5060644020340
HMC Premedical S.p.A	Medicina®	RT18/100L	ENFit Short Term Feeding/Aspiration Tube	100	18	NGPOD SENSOR 01-6L0	5060644020340

HMC Premedical S.p.A	Medicina®	205101ML	Enfit Nasogastric feeding tube weighted with guidewire	110	10	NGPOD SENSOR 01- 6L4	5060644020388
HMC Premedical S.p.A	Medicina®	205121ML	Enfit Nasogastric feeding tube weighted with guidewire	110	12	NGPOD SENSOR 01- 6L3	5060644020371
HMC Premedical S.p.A	Medicina®	NGP12/120L	Enfit Nasogastric feeding tube with guidewire	120	12	NGPOD SENSOR 01- 7L6	506064402056
HMC Premedical S.p.A	Medicina®	NGP10/120L	Enfit Nasogastric feeding tube with guidewire	120	10	NGPOD SENSOR 01- 7L7	5060644020579

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Hoist Medical	STARFIT	NGGV 14-120	Nasogastric Duodenal Feeding Tube Guidewire PUR ENFit	120cm	14	NGPOD SENSOR 01- 7L6	5060644020562

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
--------------	-------	----------	-------------	-------------	----------------------	------------------------------	--------

Kangaroo™	Kangaroo	8884711253E	Nasogastric Feeding Tube - 7 g Weighted Dobbhoff Tip, Rigid Port, Stylet	109	12	NGPOD SENSOR 01- 6LE	5060644020487
Kangaroo™	Kangaroo	8884720825E	Nasogastric Feeding Tube - 7 g Weighted Dobbhoff Tip, Rigid Port, Stylet	91	8	NGPOD SENSOR 01- 092	5060644020005

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Maxter-Catheters	ENTRAL™	SFT885	Nasogastric Feeding Tube	85	8	NGPOD SENSOR 01- 5L4	5060644020227
Maxter-Catheters	ENTRAL™	DFT12-100	Nasogastric Feeding Tube	100	12	NGPOD SENSOR 01- 5LF	5060644020333

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
NeoMed®	NeoConnect™	PFTL6.5P- NC	Enfit compliant Neoconnect feeding tubes	90	6.5	NGPOD SENSOR 01- 5L8	5060644020265
NeoMed®	NeoConnect™	PFTL8.0P - NC	Enfit compliant Neoconnect feeding tubes	90	8	NGPOD SENSOR 01- 5L8	5060644020265
NeoMed®	NeoConnect™	PFTM10.0P-NC	ENFit Compliant NeoConnect Feeding Tubes	109	10	NGPOD SENSOR 01- 6LB	5060644020456

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Nutricia	Flocare®	594569	Flocare® PUR tube with stylet	110	8	NGPOD SENSOR 01-7L1	5060644020517
Nutricia	Flocare®	594570	Flocare® PUR tube with stylet	110	10	NGPOD SENSOR 01-7L1	5060644020517
Nutricia	Flocare®	594571	Flocare® PUR tube with stylet	110	12	NGPOD SENSOR 01-7L1	5060644020517
Nutricia	Flocare®	594572	Flocare® PUR tube with stylet	110	14	NGPOD SENSOR 01-7L1	5060644020517
Nutricia	Flocare®	624856	Flocare® PUR tube with Guidewire/Suction	110	14	NGPOD SENSOR 01-7L9	5060644020609

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
REDA	SMOOTH- <i>flo</i>	83-1290	Enteral feeding tube	90cm	12	NGPOD SENSOR 01-5LA	5060644020289
REDA	SMOOTH- <i>flo</i>	83-1490	Enteral feeding tube	90cm	14	NGPOD SENSOR 01-5LA	5060644020289

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
NB: Ryles tubes are compatible with the NGPOD System when being used for feeding only (not drainage)							
Ryles	Ryles	IC 4032-12	Nasogastric tube for enteral feeding/aspiration catheter with integrated colour coded dual mode mount	105cm	12	NGPOD SENSOR 01-5LE	5060644020326

Ryles	Ryles	IC 4032-14	Nasogastric tube for enteral feeding/aspiration catheter with integrated colour coded dual mode mount	105cm	14	NGPOD SENSOR 01-6LO	5060644020340
Ryles	Ryles	IC 4032-16	Nasogastric tube for enteral feeding/aspiration catheter with integrated colour coded dual mode mount	105cm	16	NGPOD SENSOR 01-6LO	5060644020340

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
TIK	TIK	54872	Stomach tube	80cm	16	NGPOD SENSOR 01-4LA	5060644021125
TIK	TIK	54873	Stomach tube	80cm	18	NGPOD SENSOR 01-4LA	5060644021125
TIK	TIK	54874	Stomach tube	80cm	20	NGPOD SENSOR 01-4LA	5060644021125

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Unomedical	Unomedical	11003182	Ryles tube	125cm	12	NGPOD SENSOR 01-7L4	5060644020548
Unomedical	Unomedical	11004182	Ryles tube	125cm	14	NGPOD SENSOR 01-7L4	5060644020548
Unomedical	Unomedical	11005182	Ryles tube	125cm	16	NGPOD SENSOR 01-7L4	5060644020548
Unomedical	Unomedical	10001182	Duodenal tube	125cm	8	NGPOD SENSOR 01-7L4	5060644020548
Unomedical	Unomedical	10019182	Duodenal tube	125cm	10	NGPOD SENSOR 01-7L4	5060644020548
Unomedical	Unomedical	10021182	Duodenal tube	125cm	14	NGPOD SENSOR 01-7L4	5060644020548
Unomedical	Unomedical	10022182	Duodenal tube	125cm	16	NGPOD SENSOR 01-7L4	5060644020548
Unomedical	Unomedical	10023182	Duodenal tube	125cm	18	NGPOD SENSOR 01-7L4	5060644020548

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
VYGON	Vygon – Value Life	391.14	Gastro-duodenal feeding tube Levin type	125	14	NGPOD SENSOR 01- 7L6	5060644020562
VYGON	Vygon – Value Life	392.12	Gastro-duodenal tube with weighted tip -PVC	125	12	NGPOD SENSOR 01- 7LA	5060644020616
VYGON	Vygon – Value Life	392.14	Gastro-duodenal tube with weighted tip -PVC	125	14	NGPOD SENSOR 01- 7LA	5060644020616