

This is an Appendix to the NGPod System Instructions for use (IFU). Please refer to the full IFU before operating the device.

Prior to use of the NGPOD system, compatibility with the nasogastric tube(s) (NGT) to be used should be established.

The table below lists NGPOD sensor code to be used with NGTs that NGPOD has been validated with.

If the nasogastric tube you wish to use is not included in the below list, please contact NGPOD Global www.ngpodglobal.com, info@ngpodglobal.com for advice or to arrange a validation.

Nasogastric Tube Compatibility Reference Table

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Avanos/Halyard	CORFLO	42-7368	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	91	8	NGPOD SENSOR 01-5LC	5060644020302
Avanos/Halyard	CORFLO	42-9368	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	91	8	NGPOD SENSOR 01-5LC	5060644020302
Avanos/Halyard	CORFLO	42-7361	CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet	91	10	NGPOD SENSOR 01-5LC	5060644020302
Avanos/Halyard	CORFLO	42-9361	Nasogastric/Nasointestinal Feeding tube with Stylet	91	10	NGPOD SENSOR 01-5LC	5060644020302

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Fresenius Kabi	Freka Tube	7981840	Nasogastric tube with ENfit and guide wire	120cm	8	NGPOD SENSOR 01-7L5	5060644020555
Fresenius Kabi	Freka Tube	7981842	Nasogastric tube with ENfit and guide wire	120cm	10	NGPOD SENSOR 01-7L5	5060644020555

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
GB UK Enteral	Nutricare	NC08.092ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	92	8	NGPOD SENSOR 01-092	5060644020005
GB UK Enteral	Nutricare	NC10.092ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	92	10	NGPOD SENSOR 01-092	5060644020005
GB UK Enteral	Nutricare	NC12.092ISOSAF/ NC12.092ISOSAF. ROW	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	92	12	NGPOD SENSOR 01-092	5060644020005
GB UK Enteral	Nutricare	NC08.092ISOSAF. ROW	Nasogastric feeding tube, polyurethane	92	8	NGPOD SENSOR 01-092	5060644020005

GB UK Enteral	Nutricare	NC08.080ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	80	8	NGPOD SENSOR 01-5L1	5050644020197
GB UK Enteral	Nutricare	NC06.080ISOSAF	Long term Nasogastric feeding tube, polyurethane, with ENFIT connector and guidewire	80	6	NGPOD SENSOR 01-5L1	5050644020197

MANUFACTURER	BRAND	MFR CODE	Description	Length (cm)	French Gauge (FR)	Compatible Sensor Variant	UDI-DI
Nutricia	Flocare®	594569	Flocare® PUR tube with stylet	110	8	NGPOD SENSOR 01-7L1	5060644020517
Nutricia	Flocare®	594570	Flocare® PUR tube with stylet	110	10	NGPOD SENSOR 01-7L1	5060644020517
Nutricia	Flocare®	594571	Flocare® PUR tube with stylet	110	12	NGPOD SENSOR 01-7L1	5060644020517
Nutricia	Flocare®	594572	Flocare® PUR tube with stylet	110	14	NGPOD SENSOR 01-7L1	5060644020517