

# POWERFUL + PORTABLE

Mobile Blood Clot Prevention



The VenaPro (TM) DVT prevention system is specifically designed with the patient in mind. The completely tubeless and cordless system allows for use in both the home and clinical environments to help prevent the onset of a DVT during and after surgery.

The unique VenaPro calf compression system stimulates blood flow in the lower extremities by mimicking natural muscle contractions of the calf.

#### VenaPro Can Be Used To:

- Aid in the prevention of DVT
- Enhance blood circulation
- Diminish post-operative pain and swelling
- Reduce wound healing time



TUBELESS &  
CORDLESS  
OPERATION

EASY TO USE

UP TO 20 HOUR  
RECHARGEABLE  
BATTERY

PORTABLE &  
LIGHTWEIGHT

USB COMPLIANCE  
MONITORING

## HOW IT WORKS




#### Place on Calf

Open wrap and place around your calf with inflatable part of wrap on calf.



#### Turn On

Press and hold power button on both units for 3 seconds and release. Device will start.



#### Pressure & Charging

The device will gently fill with air creating pressure on your calf. Recharge device when needed. Scan QR to see product tutorial.