


POWERFUL + PORTABLE

Mobile Blood Clot Prevention



CIRCUL8[®]

— PRO —

CIRCUL8[®]

— PRO —

The CIRCUL8[®] PRO is intended to be an easy to use portable system, prescribed by a physician, for use in the home or clinical setting to help prevent the onset of DVT in patients by stimulating blood flow in the extremities (stimulating muscle contractions).

CIRCUL8 PRO Can Be Used To:

- Aid in the prevention of DVT
- Enhance blood circulation
- Diminish post-operative pain and swelling
- Reduce wound healing time

EASY TO USE

UP TO 20 HOUR
RECHARGEABLE
BATTERY LIFE

PORTABLE &
LIGHTWEIGHT

COMPLIANCE
COMMUNICATION
ADAPTER



HOW IT WORKS



Place on Calf

Open wrap and place around your calf with inflatable part of wrap on calf.



Turn On

Press and hold power button on both units for 3 seconds and release. Device will start



Pressure & Charging

The device will gently fill with air creating pressure on your calf. Recharge device when needed.