

I-CAM Chairside Conversion



Pre-Surgery



Record Supporting
Anatomy With IOS.
Use Fiducial Markers

We recommend PLY files.
STL files are still acceptable



Finish Implant
Surgery and Place
MUAs

Post-Surgery



Using the I-CAM
Flags, Record Implant
Positions



Position Scan Bodies to
MUAs & Scan With IOS



Suture Patient &
Acquire IOS of the
Scan Bodies



Send All Acquired
Data to ITXPROS



Delayed Conversion
Within 48 Hours



Same Day
Conversion Within
90 Minutes

Requires Pre-Booking




SCAN QR CODE
TO START

Updated October, 2025

I-CAM ALL-ON-ITX



We recommend PLY files.
STL files are still acceptable



SCAN QR CODE
TO START

Updated October, 2025