

The **PICASSO[®] plus**

Specimen Radiography System

**See Things Differently
in the Gross Room
to Improve Efficiency**




The Image Blender™
Combine X-ray and Optical
Images for the Perfect Blend!



Large Imaging Area
A 10" x 12" detector for larger
specimens, like mastectomies.



Voice Control
Innovative speech recognition
technology for true
hands-free navigation



Selenium Detector
Amorphous Selenium
technology provides crisp,
high quality images.



Optional Pablo™ Tablet
Convenient access to images
and annotation tools anywhere
in the gross room.



Gesture Control
Operate the optional Pablo
Tablet without touching it.

KUBTEC[®]
MEDICAL




KUBTEC@KUBTEC.COM
1.203.364.8544
WWW.KUBTEC.COM



The Next Generation of Specimen Imaging for the Gross Room.


In today's fast-paced pathology lab, PAs need X-ray imaging in the gross room to reduce turnaround time and eliminate specimen transport to radiology or mammography, while enabling accurate identification of clips, markers, and microcalcifications.

The PICASSO Plus System redefines imaging in the Pathology Lab. It includes workflow features designed specifically for the gross room as well as being the only Pathology specimen imaging system featuring Amorphous Selenium Direct Capture Imaging Technology.




The Image Blender™ See it All Together Prior to Grossing the Specimen

Blend both the X-ray and optical images to quickly locate surgical clips, margins, and micro-calcifications to best gross specimens and select tissue to place into cassettes.



+



=



The Pablo™ Tablet See Images Anywhere in the Gross Room

With the optional Pablo™ Tablet you can conveniently access your images and annotation tools anywhere in the gross room, helping ensure a smooth transition to the grossing station and beyond, so you can keep working without interruption.

The Pablo tablet features both voice and gesture control for touch free operation and stows away inside the PICASSO Plus System when not in use.