

Additional support

For more details on how to set up and use your Nuance Audio Glasses, scan the QR code



Support Contacts

Call Center
1-844-844-0919
8:00 AM - 8:00 PM
Monday through Friday
(EST time zone)

Website
www.nuanceaudio.com/support

THE COLLECTION

Styles

2 timeless models available in 3 different sizes and 2 elegant colors.




SQUARE 54/56



PANTHOS 48

Accessories & packaging

Designed to be discreet and easy to use while maintaining the same sophisticated look as the eyewear.



CHARGING PAD*



CASE



REMOTE CONTROL**

Nuance Audio, Best Paired with Transitions® Lenses



Nuance Audio, best paired with Transitions® lenses to experience seamless hearing and visual performance. With the ultra-dynamic Transitions® lenses, empower your everyday life in any light.

Available in Grey and Emerald as non-prescription lenses at purchase, and with a widest color range as on-demand prescription lenses.

Customize Nuance Audio with prescription lenses



Option to add your preferred optical lenses based on your vision needs.

Transitions and the Transitions logo are registered trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. Photochromic performance and polarization is influenced by temperature, UV exposure and lens material.

* available for purchase separately

** sold separately


nuance audio
in tune with life

The new invisible hearing solution

Our Manifesto

Empowering the world to hear, see, and connect.

Tune into a world where the beauty of sound knows no limits, where you can hear every detail and immerse yourself in the richness of life.

Because the better you hear, the more you can connect with what matters.

BREAKTHROUGH INNOVATION

Hear better in noisy environments with the all-in-one hearing glasses



Hear what you see with directional microphones

Hear comfortably with discreet open-ear speakers

The new invisible open-ear OTC hearing solution for perceived mild to moderate hearing loss.

Nuance Audio blends sound engineering and world-class smart eyewear design for a stylish, high-tech solution engineered to transform lives.

What Makes Nuance Audio Hearing Glasses Unique

Invisible hearing aid

Open-ear comfort

All-in-one solution

Directional amplification

Easy to use

USE CASES

Restaurants & Crowded Places

Struggling to follow conversations in noisy restaurants? Focus on the voices you care about with directional microphones limiting background noise interference, making dining out more enjoyable.

Workplaces

Busy meetings and bustling offices? Experience enhanced hearing that reduces listening effort, to help you stay sharp and confident throughout the day.

Family Gatherings

At lively family dinners or reunions, don't miss the shared laughter and heartfelt moments. Focus on the voices that matter most, strengthening your connections and making everyday experiences more enjoyable and accessible.

Refer to Nuance Audio training materials and terms of sale for product warranty, lens installation and return policies.

©EssilorLuxottica- April 2025- 37665058_PRO_OTH

USER EXPERIENCE

Nuance Audio Glasses use directional microphones and open-ear speakers to provide a comfortable and natural listening experience.

The Nuance Audio App is essential for setting up and personalizing the glasses to the user's hearing needs and preferences.

Audio modes

Switch between two audio modes: Frontal, to focus on face-to-face conversations, and All-around, to amplify all surrounding sounds.

Presets

Choose from the different pre-configured amplification settings (A, B, C or D) the one that best suits your hearing needs.

Background noise


Customize the background noise level you want to hear with three levels of noise reduction.

Volume level

Adjust the volume of amplification on a scale from 1 to 5.

Download for iOS

Download for Android



The Nuance Audio App is compatible with iOS devices running iOS 16 or later (release 0.3.7 onward) and Android devices running Android 12 or later with SDK 24 (release 0.4.1 onward).