



OUTBACK

ADVANCED OUTDOOR REHAB POWER

BRYAN ANDERSON
Iraq War Veteran
and Quantum®
Spokesperson



FDA Class II
Medical Device*

QUANTUM® OUTBACK
shown with Q-Logic 3e
controller and optional
10" safe seat elevation



MAKE YOUR OWN PATH

Quantum Rehab® welcomes adventurers and nature enthusiasts to unleash their wild side with the release of our cutting-edge 4x4 all-terrain power chair, the Quantum® Outback! Designed for those who dare to defy conventional limits, this rugged mobility solution promises to elevate any outdoor adventure, whether navigating wooded trails or exploring uncharted terrain.



QUANTUM®
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	Standard: 15" pneumatic
MAXIMUM SPEED ¹	Up to 5 mph (6 mph available with optional upgrade)
GROUND CLEARANCE	4.92" (at motors)
TURNING RADIUS ²	49.2"
OVERALL LENGTH ²	42"
BASE WIDTH	28.8"
SUSPENSION	Independent SRS
RANGE ³	Up to 16 miles at 300 lbs. Up to 17 miles at 220 lbs.
SEATING SIZES	
Q-Captain Seat	Width: 16", 18", 20", 22" Depth: 16", 18", 20", 22"
TRU-Balance® 3 Power Tilt	Width: 16-22" Depth: 16-22"
TRU-Balance® 3 Power Lift & Tilt	Width: 16-22" Depth: 16-22"
SEAT-TO-FLOOR HEIGHTS	
Q-Captain Seat ⁵	19.5", 20.5", 21.5"
TRU-Balance® 3 Tilt ⁵	19", 20", 21"
TRU-Balance® 3 Lift & Tilt ⁵	19.5", 20.5", 21.5"
BATTERY SIZE ³	
BATTERY WEIGHT ⁶	58.6 lbs. (each)
AVAILABLE ELECTRONICS ⁷	120A Q-Logic 3
BATTERY CHARGER	8A Off-board
BRAKES	Regenerative and Electromechanical
DRIVETRAIN	4 motor, independent drive
MAXIMUM OBSTACLE CLIMBING	4"
WEIGHT CAPACITY ⁴	300 lbs. (TRU-Balance® 3) 350 lbs. (Q-Captain Seating)
BASE WEIGHT	168.6 lbs.
HCPCS CODE ⁸	K0899 Power Mobility Device not coded by DME PDAC

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and pneumatic tires. Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

*The Quantum® Outback is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

Standard Features

- 4-pole motors
- New camo color choices
- Full LED lighting package
- SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions

Options

- 10" of power adjustable safe seat elevation at speeds of up to 3.2 mph
- Quantum® backup camera

Colors

CAMO



Black



Desert



White

MATTE



Cadet



Tan



FDA Class II Medical Device *

QUANTUM®

#1 FOR REHAB POWER

(US) 866-800-2002

(CAN) 888-570-1113

quantumrehab.com

