



CM-T Flex

A New innovative design for the treatment for APHAKIA
Sutureless Scleral Fixation Hydrophilic Foldable IOL

CM-T Flex

Features

- Single piece T shaped flexible haptic foldable intraocular lens for aphakic patients
- It is made up of biocompatible hydrophilic acrylic material
- CM T flex IOL is for single vision
- This IOL can be used for all aphakia indications such as surgical aphakia, traumatic aphakia and secondary aphakia
- Two T shaped anchor haptics present on either side
- Can be placed horizontally and vertically.
- The two tunnels must be lined-up 180° from each other for good centration
- The angulations between the optic and haptic prevent the rubbing of posterior surface of the iris thus helping in better pupil dilation
- Provides sutureless sclera
- Fixed haptics provides more stability
- Its large optic diameter can be used in pupils of various size and hence avoids lens glare during night vision

Technical Specification

Material	:	Hydrophilic 26%
Optic Size	:	6.00 mm
Overall Size	:	13.75mm
Angulation	:	10°
Refractive Index	:	1.460
Diopter Range	:	-15D to 55D (In 0.5 steps)
Sterilization	:	Steam
A-Constant	:	118.0

Manufactured by



84/1, Aavin Dairy Farm Road,
SIDCO Industrial Estate,
Ambattur, Chennai-600098, INDIA

Marketed by



APPASAMY ASSOCIATES (P) LTD.,
20, SBI Officers' Colony, First Street, Arumbakkam, Chennai - 600 106, INDIA
Tel : (91-44) 30101400, 30101401
Email: info@appasamy.com Website: www.appasamy.com