


nAspro-BBY

Blue Blocking Yellow

Nano Aspheric Precision Optics

Aspheric Foldable Posterior Chamber
Intraocular lens



Enhanced 360° Square Edge
Blue Blocking Optics



Blue Blocking Yellow

Technical Specification

Model	NASY 207
Overall length	12.50 mm
Optic Diameter	6 mm
Haptic Angulations	0°
Optics Design	Biconvex
Diopter Range	-15.0D to +53.0D (0.5 Diopter increments)
Material	Biocompatible HEMA Yellow
Refractive Index	1.46
Sterilization	Moist heat sterilization
A constant	118.0
Special Features	Aspheric Optics with 360° Square edge
Water Content	26%

Manufactured by



84/1, Aavin Dairy Farm Road,
SIDCO Industrial Estate,
Ambattur, Chennai-600098, INDIA

Marketed by



APPASAMY ASSOICATES (P) LTD.,
20, SBI Officers' Colony, First Street, Arumbakkam, Chennai - 600 106, INDIA
Tel : (91-44) 30101400, 30101401
Email: info@appasamy.com Website: www.appasamy.com