

Compatible Devices

Sway Medical continuously updates its software to ensure the best performance, security, and user experience. As a result, the minimum system requirements may change over time. Sway Medical cannot guarantee ongoing compatibility for all devices, and we recommend regularly reviewing the latest device and system specifications before use.

iOS

- **iPhones** - Requires iOS 15.0 or later
 - **iPads** - Requires iOS 16.0 or later
-

Android (Phones and Tablets)

- Requires Android 10.0 or later
 - Device must include a **gyroscope** for functional testing
-

Unsupported Devices

Any device that does not meet the above minimum software requirements are unsupported for testing. The following devices are not supported due to known compatibility issues or missing required sensors:

- **iPad Mini 4**
- **iPad Air 2**
- **Samsung Galaxy Tab A Series**
- **Samsung Galaxy S8 and older models**
- **Motorola Moto G Play**
- **Motorola Moto G 5g**
- **Samsung Galaxy A Series** - The following models and earlier generations are not supported:
 - Galaxy A15 and older (A1x series)
 - Galaxy A23 and older (A2x series)
 - Galaxy A33 and older (A3x series)
 - Galaxy A41 and older (A4x series)
 - Galaxy A52 and older (A5x series)
 - Galaxy A62 and older (A6x series)
 - Galaxy A72 and older (A7x series)
 - Galaxy A82 and older (A8x series)
- **Galaxy Z-Fold and Z-Flip** - While they might work, we cannot guarantee it at this time and therefore do not recommend their use.