

PERIPHERAL INTERVENTION




Distail™ Microcatheter

Navigating Through Selective Vessels

Comprehensive specifications

- Tip curves fit for different clinical demands
- 5 types of O.D. with 3 different length for meeting clinical requirements
- Provides well-equipped accessories for higher efficiency



Clear Visibility under Fluoroscopy

- One or two markers




Three-Layer Structure

- Excellent torquability and high pressure endurance
- Gradually decreased stiffness towards the tip facilitates passing tortuous anatomies smoothly




Thin Wall with Big Lumen



Embolic Agent Compatibility

O.D. (F)	Tip I.D.	Pressure (psi)	Compatible Embolic Agent
1.8F	0.017"	800	300µm Particles / 300µm Microspheres
2.0F	0.019"	800	500µm Particles / 500µm Microspheres
2.2F	0.021"	800	0.018" Coil / 500µm Particles / 700µm Microspheres
2.5F	0.025"	800	0.018" Coil / 900µm Particles / 900µm Microspheres
2.7F	0.027"	800	1000 µm Particles / 1200µm Microspheres

Hydrophilic Coating



Special design for radial access

150cm via radial artery for Oncology Intervention

Tip Curve	Length (cm)	Marker	1.8F	2.0F	2.2F	2.5F	2.7F
Straight	150	1	3618501ST00	3620501ST00	3622501ST00	3625501ST00	3627501ST00
45°			36185014500	36205014500	36225014500	36255014500	36275014500
70°			36185017000	36205017000	36225017000	36255017000	36275017000
J			3618501JC00	3620501JC00	3622501JC00	3625501JC00	3627501JC00
Swan Neck			3618501SN00	3620501SN00	3622501SN00	3625501SN00	3627501SN00
Cobra			3618501CB00	3620501CB00	3622501CB00	3625501CB00	3627501CB00
Straight	150	2	3618502ST00	3620502ST00	3622502ST00	3625502ST00	3627502ST00
45°			36185024500	36205024500	36225024500	36255024500	36275024500
70°			36185027000	36205027000	36225027000	36255027000	36275027000
J			3618502JC00	3620502JC00	3622502JC00	3625502JC00	3627502JC00
Swan Neck			3618502SN00	3620502SN00	3622502SN00	3625502SN00	3627502SN00
Cobra			3618502CB00	3620502CB00	3622502CB00	3625502CB00	3627502CB00

*Without Guidewire

Tip Curve	Length (cm)	Marker	1.8F	2.0F	2.2F	2.5F	2.7F
Straight	150	1	3618501ST01	3620501ST01	3622501ST01	3625501ST01	3627501ST01
45°			-	36205014501	36225014501	36255014501	36275014501
70°			-	36205017001	36225017001	36255017001	36275017001
J			-	3620501JC01	3622501JC01	3625501JC01	3627501JC01
Swan Neck			-	3620501SN01	3622501SN01	3625501SN01	3627501SN01
Cobra			-	3620501CB01	3622501CB01	3625501CB01	3627501CB01
Straight	150	2	3618502ST01	3620502ST01	3622502ST01	3625502ST01	3627502ST01
45°			-	36205024501	36225024501	36255024501	36275024501
70°			-	36205027001	36225027001	36255027001	36275027001
J			-	3620502JC01	3622502JC01	3625502JC01	3627502JC01
Swan Neck			-	3620502SN01	3622502SN01	3625502SN01	3627502SN01
Cobra			-	3620502CB01	3622502CB01	3625502CB01	3627502CB01

150cm Distail equipped with 170cm micro guidewire
All sizes equipped with: shaping needle, guidewire introducer, torque device and syringe

One Marker

Tip Curve	Length (cm)	1.8F	2.0F	2.2F	2.5F	2.7F
Straight	110	3618101ST00	3620101ST00	3622101ST00	3625101ST00	3627101ST00
45°		36181014500	36201014500	36221014500	36251014500	36271014500
70°		36181017000	36201017000	36221017000	36251017000	36271017000
J		3618101JC00	3620101JC00	3622101JC00	3625101JC00	3627101JC00
Swan Neck		3618101SN00	3620101SN00	3622101SN00	3625101SN00	3627101SN00
Cobra	130	3618101CB00	3620101CB00	3622101CB00	3625101CB00	3627101CB00
Straight		3618301ST00	3620301ST00	3622301ST00	3625301ST00	3627301ST00
45°		36183014500	36203014500	36223014500	36253014500	36273014500
70°		36183017000	36203017000	36223017000	36253017000	36273017000
J		3618301JC00	3620301JC00	3622301JC00	3625301JC00	3627301JC00
Swan Neck		3618301SN00	3620301SN00	3622301SN00	3625301SN00	3627301SN00
Cobra		3618301CB00	3620301CB00	3622301CB00	3625301CB00	3627301CB00

*Without Guidewire

Tip Curve	Length (cm)	1.8F	2.0F	2.2F	2.5F	2.7F
Straight	110	-	3620101ST01	3622101ST01	3625101ST01	3627101ST01
45°		-	36201014501	36221014501	36251014501	36271014501
70°		-	36201017001	36221017001	36251017001	36271017001
J		-	3620101JC01	3622101JC01	3625101JC01	3627101JC01
Swan Neck		-	3620101SN01	3622101SN01	3625101SN01	3627101SN01
Cobra	130	-	3620101CB01	3622101CB01	3625101CB01	3627101CB01
Straight		3618301ST01	3620301ST01	3622301ST01	3625301ST01	3627301ST01
45°		-	36203014501	36223014501	36253014501	36273014501
70°		-	36203017001	36223017001	36253017001	36273017001
J		-	3620301JC01	3622301JC01	3625301JC01	3627301JC01
Swan Neck		-	3620301SN01	3622301SN01	3625301SN01	3627301SN01
Cobra		-	3620301CB01	3622301CB01	3625301CB01	3627301CB01

130cm Distail equipped with 150cm micro guidewire

110cm Distail equipped with 130cm micro guidewire

All sizes equipped with: shaping needle, guidewire introducer, torque device and syringe

Two Markers

Tip Curve	Length (cm)	1.8F	2.0F	2.2F	2.5F	2.7F
Straight	110	3618102ST00	3620102ST00	3622102ST00	3625102ST00	3627102ST00
45°		36181024500	36201024500	36221024500	36251024500	36271024500
70°		36181027000	36201027000	36221027000	36251027000	36271027000
J		3618102JC00	3620102JC00	3622102JC00	3625102JC00	3627102JC00
Swan Neck		3618102SN00	3620102SN00	3622102SN00	3625102SN00	3627102SN00
Cobra	130	3618102CB00	3620102CB00	3622102CB00	3625102CB00	3627102CB00
Straight		3618302ST00	3620302ST00	3622302ST00	3625302ST00	3627302ST00
45°		36183024500	36203024500	36223024500	36253024500	36273024500
70°		36183027000	36203027000	36223027000	36253027000	36273027000
J		3618302JC00	3620302JC00	3622302JC00	3625302JC00	3627302JC00
Swan Neck		3618302SN00	3620302SN00	3622302SN00	3625302SN00	3627302SN00
Cobra		3618302CB00	3620302CB00	3622302CB00	3625302CB00	3627302CB00

*Without Guidewire

Tip Curve	Length (cm)	1.8F	2.0F	2.2F	2.5F	2.7F
Straight	110	-	3620102ST01	3622102ST01	3625102ST01	3627102ST01
45°		-	36201024501	36221024501	36251024501	36271024501
70°		-	36201027001	36221027001	36251027001	36271027001
J		-	3620102JC01	3622102JC01	3625102JC01	3627102JC01
Swan Neck		-	3620102SN01	3622102SN01	3625102SN01	3627102SN01
Cobra	130	-	3620102CB01	3622102CB01	3625102CB01	3627102CB01
Straight		3618302ST01	3620302ST01	3622302ST01	3625302ST01	3627302ST01
45°		-	36203024501	36223024501	36253024501	36273024501
70°		-	36203027001	36223027001	36253027001	36273027001
J		-	3620302JC01	3622302JC01	3625302JC01	3627302JC01
Swan Neck		-	3620302SN01	3622302SN01	3625302SN01	3627302SN01
Cobra		-	3620302CB01	3622302CB01	3625302CB01	3627302CB01

130cm Distail equipped with 150cm micro guidewire

110cm Distail equipped with 130cm micro guidewire

All sizes equipped with: shaping needle, guidewire introducer, torque device and syringe

NXT MED

Official Distribution Partner

Address: Fritz-Vomfelde-Straße 6 40547 Düsseldorf, Germany

Email: hello@nxt-med.com

TEL: +49 211 58678600

@NXT MED



www.nxt-med.com



APT Medical

Manufacturer

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. of China
Email: international@aptmed.com
TEL: +86 755-83591013

@APT Medical



www.aptmedical.com

CE 1639