



MYCO[®]

M E D I C A L

The**SharpChoice**[®]

NMSDC Certified Diversity Supplier

DATASHEET

RELI[®] Safety Slide[™] Blood Collection/Infusion Sets



The Reli[®] Safety Slide[™] Safety Blood Collection and Infusion Set is a sterile, single-use, winged blood collection device designed for safety and efficiency. The safety shield can be activated in-vein by sliding the shield over the needle during removal from skin helping to reduce the risk of accidental needle sticks.

- Flexible color-coded wing for improved tactile feel and easy gauge identification
- Cartridge sliding design for simple and safe use

Pictured: Reli[®]
Safety Slide[™] Blood
Collection/Infusion Set
21G x 3/4" x 12"
(SBCS21G-12")

- Not made with natural rubber latex
- OSHA compliant
- Sterile blister pack, single use only
- Siliconized needle for smooth insertion
- Audible click when needle is in the final lock position
- One- or two-handed irreversible activation
- Features single-handed safety shield activation for convenience and ease of use

Reli® Safety Slide™ Blood Collection/Infusion Sets (Sterile)

MYCO CODE	NEEDLE GAUGE	TUBE LENGTH	COLOR	UNIT OF MEASURE
SBCS21G-12	21G X ¾"	30 cm (11.8")	Green	50 ea/bx; 200 ea/cs
SBCS23G-12	23G x ¾"	30 cm (11.8")	Blue	50 ea/bx; 200 ea/cs
SBCS25G-12	25G x ¾"	30 cm (11.8")	Orange	50 ea/bx; 200 ea/cs

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to "Bringing Reliability, Innovation and Value" to the Healthcare industry with an emphasis on "clinical preference" surgical devices and disposables.



Scan to learn more

or receive complimentary evaluation samples

mycomed.co/safety-slide

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com



MYCO
MEDICAL
TheSharpChoice®

2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com

www.mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2025 MYCO Medical, Inc. All rights reserved.
MYCO® is a registered trademark of MYCO Medical, Inc.

LIT-01028 Rev 0 Eif 10.02.2025