


Sterile Disposable Skin Stapler

The **RELI® Disposable Skin Stapler** is a sterile, single-patient-use mechanical instrument designed to deliver precise, rectangular-shaped stainless steel staples for routine wound closure. Its handle is textured for a user-friendly experience, and the angled head provides a clear view to ensure precise staple placement and allows the pin to staple the tissue easily.



RELI® Disposable Skin Stapler

MYCO CODE	DESCRIPTION	STAPLE MATERIAL	SIZE	STAPLES	UNIT OF MEASURE
SST-35W	RELI® Disposable Skin Stapler	316L Stainless Steel	Wide	35 Staples	6 Units/Box, 21 Boxes/Case

The **RELI® Disposable Skin Stapler** is a sterile, single-use device supplied in a single-pack blister tray.



Scan to learn more

or receive complimentary
evaluation samples

mycomedical.co/RELI-Disposable-Skin-Stapler