



MYCO
MEDICAL

TheSharpChoice®

NMSDC Certified Diversity Supplier



QUINCKE POINT Spinal Needles

RELI® Quincke Point Spinal Needles are the value alternative to leading brands and offer superior quality in a broad range of sizes.




- Needle gauges available: 18G to 27G
- Needle lengths available: 1.5" to 8"
- Computerized grinding ensures sharp, consistent, high-quality needles
- Exact fit between stylet and cannula bevel to minimize coring and trauma
- Large, clear hubs for enhanced tactile feel and visualization of CSF
- Smooth insertion resistance provides for easy identification of tissue layers
- Available in sterile and non-sterile packaging



18G Pink



19G Beige



20G Yellow




22G Black




23G Blue



24G Purple



25G Orange



26G Brown



27G Gray



Large clear hub



Matched-point grinding



Broad range of sizes

BRAND (COMPANY)	LARGE HUBS	CLEAR HUBS	COLOR-CODED STYLETS	FITTED, CLOSE TOLERANCE STYLET	SHORT BEVEL OFFERINGS
Reli® Quincke Point Spinal Needle (MYCO Medical)	✓	✓	✓	✓	✓
Spinocan (B. Braun)		✓		✓	
Quincke Spinal Needle (BD)			✓	✓	
Quincke Spinal Needle (Kimberly-Clark)		✓	✓	✓	✓

Gauge	Length	STERILE MYCO CODE	UNIT OF MEASURE	NON-STERILE MYCO CODE	UNIT OF MEASURE	Gauge	Length	STERILE MYCO CODE	UNIT OF MEASURE	NON-STERILE MYCO CODE	UNIT OF MEASURE
18G	3"	SN18G301	25 ea/bx; 4 bx/cs	SN18G300	100 ea/bx	23G	2"	N/A	N/A	SN23G200	100 ea/bx
18G	3.5"	SN18G351	25 ea/bx; 4 bx/cs	SN18G350	100 ea/bx	23G	3.5"	SN23G351	25 ea/bx; 4 bx/cs	SN23G350	100 ea/bx
18G	5"	SN18G501	25 ea/bx; 4 bx/cs	SN18G500	100 ea/bx	23G	5"	SN23G501	25 ea/bx; 4 bx/cs	SN23G500	100 ea/bx
18G	6"	N/A	N/A	SN18G600	100 ea/bx	23G	6"	SN23G601	25 ea/bx; 4 bx/cs	SN23G600	100 ea/bx
19G	3.5"	SN19G351	25 ea/bx; 4 bx/cs	SN19G350	100 ea/bx	23G	7"	N/A	N/A	SN23G700	100 ea/bx
20G	3.5"	SN20G351	25 ea/bx; 4 bx/cs	SN20G350	100 ea/bx	24G	3.5"	SN24G351	25 ea/bx; 4 bx/cs	SN24G350	100 ea/bx
20G	5"	SN20G501	25 ea/bx; 4 bx/cs	SN20G500	100 ea/bx	25G	1.5"	SN25G151	25 ea/bx; 4 bx/cs	SN25G150	100 ea/bx
20G	6"	SN20G601	25 ea/bx; 4 bx/cs	SN20G600	100 ea/bx	25G	2.5"	SN25G251	25 ea/bx; 4 bx/cs	SN25G250	100 ea/bx
22G	1.5"	SN22G151-RW	25 ea/bx; 4 bx/cs	SN22G150-RW	100 ea/bx	25G	3"	SN25G301	25 ea/bx; 4 bx/cs	SN25G300	100 ea/bx
22G	2.5"	SN22G251-RW	25 ea/bx; 4 bx/cs	SN22G250-RW	100 ea/bx	25G	3.5"	SN25G351	25 ea/bx; 4 bx/cs	SN25G350	100 ea/bx
22G	3"	SN22G301-RW	25 ea/bx; 4 bx/cs	SN22G300-RW	100 ea/bx	25G	5"	SN25G501	25 ea/bx; 4 bx/cs	SN25G500	100 ea/bx
22G	3.5"	SN22G351	25 ea/bx; 4 bx/cs	SN22G350	100 ea/bx	25G	6"	SN25G601	25 ea/bx; 4 bx/cs	SN25G600	100 ea/bx
22G	4.75"	SN22G4341-RW	25 ea/bx; 4 bx/cs	SN22G4340-RW	100 ea/bx	25G	7"	SN25G701	25 ea/bx; 4 bx/cs	SN25G700	100 ea/bx
22G	5"	SN22G501-RW	25 ea/bx; 4 bx/cs	SN22G500-RW	100 ea/bx	25G	8"	SN25G801	25 ea/bx; 4 bx/cs	SN25G800	100 ea/bx
22G	6"	SN22G601-RW	25 ea/bx; 4 bx/cs	SN22G600-RW	100 ea/bx	26G	3.5"	SN26G351	25 ea/bx; 4 bx/cs	SN26G350	100 ea/bx
22G	7"	SN22G701-RW	25 ea/bx; 4 bx/cs	SN22G700-RW	100 ea/bx	27G	3.5"	SN27G351	25 ea/bx; 4 bx/cs	SN27G350	100 ea/bx
22G	8"	SN22G801-RW	25 ea/bx; 4 bx/cs	SN22G800-RW	100 ea/bx						

SHORT BEVEL

Gauge	Length	STERILE MYCO CODE	UNIT OF MEASURE	NON-STERILE MYCO CODE	UNIT OF MEASURE
22G	3.5"	SN22G351-SB	25 ea/bx; 4 bx/cs	SN22G350-SB	100 ea/bx
22G	5"	SN22G501-SB	25 ea/bx; 4 bx/cs	SN22G500-SB	100 ea/bx
22G	6"	SN22G601-SB	25 ea/bx; 4 bx/cs	SN22G600-SB	100 ea/bx
22G	8"	SN22G801-SB	25 ea/bx; 4 bx/cs	SN22G800-SB	100 ea/bx

METRIC MARKINGS / ECHOGENIC STYLET TIP

Gauge	Length	STERILE MYCO CODE	UNIT OF MEASURE	NON-STERILE MYCO CODE	UNIT OF MEASURE
18G	8"	SNME18G801	25 ea/bx; 4 bx/cs	SNME18G800	100 ea/bx
20G	3.5"	SNME20G351	25 ea/bx; 4 bx/cs	SNME20G350	100 ea/bx
22G	3.5"	SNME22G351	25 ea/bx; 4 bx/cs	SNME22G350	100 ea/bx
22G	6"	SNME22G601	25 ea/bx; 4 bx/cs	SNME22G600	100 ea/bx
25G	3.5"	SNME25G351	25 ea/bx; 4 bx/cs	SNME25G350	100 ea/bx



MYCO
MEDICAL

TheSharpChoice®

NMSDC Certified Diversity Supplier



1.919.460.2535



2015 Production Drive
Apex, NC 27539 USA



www.mycomedical.com



info@mycomedical.com



Scan to learn more

or receive complimentary
evaluation samples

<https://mycomed.co/quincke-spinal-needles>

Some medical devices are not licensed for sale in Canada.
© 2026 MYCO Medical, Inc. All rights reserved.
MYCO® is a registered trademark of MYCO Medical, Inc.
LIT-00174 Rev 4 Eff 08.13.2026