

MedITEX Release Notes

Version 4.0.4

Features

1) MedITEX IVF

a) Language

All MedITEX Apps are now available in Hungarian.


Each user can specify their preferred language under „Personal Settings”.

The screenshot shows the MedITEX IVF v4.0.4.0 application window. The 'System' menu is open, and the 'Personal settings' option is selected. Below the menu, the 'Edit staff data' dialog box is visible, showing the 'Person' tab. The 'Program language' dropdown menu is open, displaying a list of languages including English, Croatian, Czech, English, Estonian, French, German, Hungarian, Italian, Kazakhstan, and Polish. The 'English' option is highlighted.

b) Inbox


Lab result printouts now include the patient’s Social Security Number (NHI in New Zealand).

MedITEX Release Notes



c) Labels

The partner's Social Security Number can now be added to Label templates.




d) Lab diagnostics

2 new units, "ml/min" and "ml/s" have been added for lab parameters.

MedITEX Release Notes

e) Report Editor

The date and time can now be reported separately.



f) Cryo Storage

For documentation purposes, "Ovarian Transplantation" has been added as a reason for thawing

Oocyte / embryo cryo storage of DEMO, Mrs (1)

Currently in storage		History									
Own material											
Cryo ID	Follicle aspiration	Cryopreserved on	Storage date	Cryo technique	Cryopreserved by	Ovarian tissue	# Days	# Ooc.	# PN	# Emb.	Straw ID
38	23/04/2025	28/04/2025				<input type="checkbox"/>	0			1	
19	01/07/2014	06/07/2014	06/07/2014								
17	01/07/2014	03/07/2014	03/07/2014								
18	01/07/2014	02/07/2014	02/07/2014								
6	10/06/2010	12/06/2010	12/06/2010								
21	11/04/2008	13/04/2008	13/04/2008	Slow	Doe, Jane Dr.						

Reason for removal

Reason:


Date:

- Deceased
- Elimination (patient order)
- Elimination (termination of contract)
- Export to overseas
- Legal requirement
- Mistakenly documented
- Ovarian transplantation
- PID abnormal
- Research
- Sent to another clinic

MedITEX Release Notes



2) MedITEX Scheduler

All alerts from MedITEX IVF will now also be displayed in Scheduler. They can be opened with the “Window” button.



Also, Ultrasound appointments can be automatically added to Scheduler as a proposed appointment. To do this, the Ultrasound event can be activated under System → Configuration/Administration → Settings → System → Interfaces → Between MedITEX Programs.


As soon as a future ultrasound is documented, e.g. in a cycle, it appears under the proposed appointments in the Scheduler.



MedITEX Release Notes

3) Interfaces

- a) Sending data to external partners in MedITEX has moved. You can find the new buttons under “Send data” in all sections.



MedITEX Release Notes

Bugs

1) MedITEX IVF

a) General

The automatic filling of the last AMH value from the laboratory diagnostics for the documented cycles has been adjusted.

Only predefined users are now displayed in the user fields.

Fixed an error message that occurred when creating the cycle summary.

Incorrectly displayed special characters such as umlauts are now displayed correctly again.

b) QM-Fields

Deactivated staff are now still displayed in already documented QM fields.

c) Inbox

Documented comments in the laboratory results are now saved again as usual for the individual laboratory parameters.

d) Cycles

A display error in the treatment performed in the cycle overview has been fixed.

The automatic setting of the reason for ET not performed has been revised.

The allergy button in the treatment overview now shows all documented allergies of the couple as usual.

The freezing of embryos from the table view in the culture has been adjusted.

e) Planned Cycles

The date when creating a new entry in the planned cycle now refers to the last selected cycle day and is not shifted by 1 day anymore.

f) Semen Analysis

Error message in the semen analysis has been fixed.

g) Medical history

Incorrectly linked previous pregnancies have been updated.

The medication prescription document now contains all the necessary data again.

h) Report Editor

Macros in letters now access all dependent placeholders again and are filled as usual.

Blank lines in comments are now carried over to the placeholders of the letters again.

MedITEX Release Notes

i) Consent Forms

Linked planned treatments in the consent forms are now displayed again in the main overview.

Multiple consent forms can now be generated automatically again.

The documented comment in a consent form is now displayed in full again in the main overview.

j) Cryo Contracts

The question regarding the termination of the contract, even when clicking "Cancel" in the thawing straw window, is no longer occurring.

Suspending an invoice no longer generates an invoice number.

2) MedITEX Scheduler

Read SMS can be marked as "read" again.

3) MedITEX Accounting

Adjustments have been made to the import of the price list.

The simultaneous creation of an invoice is now prevented in order to avoid the creation of duplicate invoice numbers.

The columns "VAT" and "VAT code" are now displayed again as usual.

Fixed an error message that occurred when documenting a discount.

4) DICOM

Unsaved files after the import are now saved again as usual in the respective "Archive" and "Failed" folders.

5) HL7

Error message when importing via TCP has been fixed.

6) Email

The sending of emails for clinics with multiple locations has been adapted.

7) Data Exporter

Error message after login has been fixed.

8) Hope

The problem of missing registration after linking with the app has been fixed.

The generation of the QR code document has been revised.