

MedITEX Release Notes

Version 4.0.5

Features

1) MedITEX IVF

a) Consents

There is now a validity period, which can be predefined per Consent.

MedITEX consents [Close]

Edit consent

Name:

Valid from: Partner signature Active

Valid until:

Valid for: Valid (days) Mandatory Erasurable

Group:

One consent can be assigned to several groups. These groups can be defined/edited in System settings (System tables).

Template: [Browse] [New]

Drag and drop a new consent template here ...

Bookmark	Use	Exists in template	
Patient: ID	<input type="checkbox"/>		[Add]
Patient: ID	<input type="checkbox"/>		[Add]
Patient: ID	<input type="checkbox"/>		[Add]
Patient: ID	<input type="checkbox"/>		[Add]
Patient: ID	<input type="checkbox"/>		[Add]
Patient: name	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	[Add]
Patient: name	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	[Add]
Patient: name	<input type="checkbox"/>		[Add]
Patient: DOB	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	[Add]
Patient: address	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	[Add]

[Icon]

MedITEX Release Notes

If the created Consent is no longer valid, a red Exclamation Mark under Status will be displayed. The expiration date of the Consent can be seen under the "Valid until" column.

The screenshot shows the MedITEX IVF Internal Debug v4.0.5.0 interface. The patient information panel on the left displays 'Muster, Frau' (01.01.1990, 35 years old) with a balance of 0.00 €. The main table lists consent records with columns for 'Th. N', 'Name', 'Zuletzt geändert', 'Status', and 'Gültig bis'. A red box highlights the 'Einwilligung...' button and the 'Gültig bis' column, which shows '16.12.2024' for the last record.

Th. N	Name	Zuletzt geändert	Status	Gültig bis
1	Allgemeine Einverständniserklär	30.08.2023	!	
1	Einverständniserklärung zur Er	30.08.2023	✓	
11	Einverständniserklärung zu wei	06.12.2024	✓	
	Allgemeine Einverständniserklär	01.08.2025	✓	
	Allgemeine Einverständniserklär	01.08.2025	✓	
12	Vertrag	01.08.2025	✓	16.12.2024


Additionally, the Expiry date field can be reported on in Report editor.

The screenshot shows the 'Report editor: queries' interface. The 'Consent forms' category is expanded, and the 'Expiry date' field is selected, indicated by a red box. Other fields listed include 'Consent document (signed)', 'Consent document (template)', 'Consent name', 'Consent status', 'Consents group', 'Date (created)', 'Date (last changed)', 'Linked cycle', 'Linked cycle planning', 'Note', 'Person', and 'Personnel'.


MedITEX Release Notes

b) Cycles

The timeframe of documentation from 2 triggering medications is extended to 48 hours. This means, if in 48 hours, triggering medications are documented, no plausibility message should appear. If the timeframe is exceeded, then the cycle will be marked as not plausible.




The column “Billing category” was added to the Main Window for Cycles and Planned cycles. It can be added/removed from the Overview via the show/hide/move columns button.



MedITEX Release Notes

For better reporting, the field “Country of birth” of the baby in the “Birth” section was converted to a Dropdown from a free text field.




c) Report Editor

The following new fields were added to the Report Editor:

Table	Field
Cycle	Patient ID Time of insemination
Culture	Age at Aspiration date
Medical History	F: long term medication M: long term medication
Last medical histories of a couple	F: long term medication M: long term medication

d) Lists

The column „Billing Category” has been added to the “Anticipated aspirations” List.




MedITEX Release Notes

e) System Settings


Defined Patient Groups can now be deactivated.

*Even after deactivating them, already defined Groups will be displayed, but they cannot be assigned.




f) Donation

A new checkbox has been added.



If this is active, then the reserving of donations is permitted only when the country is specifically listed.



MedITEX Release Notes

2) MedITEX KPI

A new Cryo KPI has been added:

Name: Live births per no. of embryos transferred after thawing

Formula: (Number of births / Number of transferred embryos after thawing) * 100

Bugs

1) MedITEX IVF

a) General

Symbol to change size of columns in Main window is visible again.

Scanning documents: Path is displayed and stored correctly.

Female patient in same sex relationship can no longer be deleted when an active cycle exists.

New referrer button is always enabled.

b) Cycle

Error when clicking on cycle details from Luteal phase has been fixed.

PID will only be shown under Culture when activated.

Missing input validations/checks for fields in the Culture and Birth sections were added.

The Note field in the Cycle now changes dynamically with the tabs. (Note: The change first takes effect after changing tabs)

Thawed semen samples can no longer be deleted. If a deletion of thawed material is necessary, this must be done over Lab diagnostics and Cryo storage of the man.

c) Planned Cycles

The free entry in planned cycle templates is now displayed correctly.

d) Examinations

Follicle ultrasounds created in a cycle are now removed from the examination once it has been deleted from the associated cycle.

e) Administration

Users can search for and select the insurance company without opening a new window, when creating/editing a patient.

The tab key in the demographics window now proceeds to the next (correct) field.

f) Report Editor

The automatically generated "Date of creation" of "Medical History" is now correctly generated.

g) Cryo Contract Management

Deleting a partial payment for an invoice will now display the correct "paid amount".

When sorting by columns under "Straw selection", the straws will now be correctly sorted.

When printing the cryo storage overview, when applying a location filter, only the storage for the correct location is now shown.

MedITEX Release Notes

h) Medical History

Previously documented uses of sterility factors remain visible, even after they've been deactivated.

2) MedITEX Scheduler

The import of holidays now functions as expected.

The "Linked resources" in the Waiting Room update immediately.

3) Hope

Error message when creating a QR-Code has been fixed.