

sigvaris

2026 CATALOGUE



Medical Compression Therapy Solutions

CUSTOMER ID

		PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
MEDICAL COMPRESSION THERAPY SOLUTIONS	ESSENTIAL	COTTON 230	● calf, calf w/Grip-Top, thigh w/Grip-Top ○ calf, thigh w/Grip-Top	20–30 mmHg 30–40 mmHg	10
		THERMOREGULATING COTTON ETH	● calf, calf w/Grip-Top	CCL1 (20–30 mmHg)	12
		OPAQUE 860 (petite & plus sizes available)	● calf, petite calf, calf plus, calf w/Grip-Top, thigh w/Grip-Top, petite thigh w/Grip-Top, extra-large thigh w/Grip-Top, pantyhose, petite pantyhose, plus-sized pantyhose ○ calf, calf plus, calf w/Grip-Top, thigh w/Grip-Top, extra-large thigh w/Grip-Top, thigh w/waist attachment	20–30 mmHg 30–40 mmHg	14, 16
	NATURAL FIBER	SEA ISLAND COTTON 220	● calf	20–30 mmHg	20
		MERINO WOOL 240	● calf	20–30 mmHg	22
		LINEN 250	● calf	20–30 mmHg	24
	PLANT BASED FIBER	BAMBOO 260	● calf	20–30 mmHg	26
	STYLE	MEDIUM SHEER 750 <i>(for women only)</i>	● calf, thigh w/Grip-Top, pantyhose	20–30 mmHg	28
		SHEER 780 <i>(for women only)</i>	● calf, thigh w/Grip-Top, pantyhose ○ calf, thigh w/Grip-Top, pantyhose	15–20 mmHg 20–30 mmHg 30–40 mmHg	30
		MICROFIBER 820 <i>(for men only)</i>	● calf, calf w/Grip-Top, thigh w/Grip-Top (extended foot sizes available)	15–20 mmHg 20–30 mmHg 30–40 mmHg	32
	MICROFIBER	MICROFIBER PATTERNS 830	● calf	20–30 mmHg	34
		SOFT OPAQUE 840 <i>(for women only)</i>	● calf, calf w/Grip-Top, thigh w/Grip-Top, pantyhose, maternity pantyhose ○ calf, thigh w/Grip-Top, pantyhose	15–20 mmHg 20–30 mmHg 30–40 mmHg	36
		CUSHIONED COTTON 360	● calf	20–30 mmHg	40
	MOTION	HIGH TECH 412C	● calf	20–30 mmHg	42
		THERMOREGULATING WOOL 422	● calf	20–30 mmHg	44
		FLOW TECH® 432	● calf	20–30 mmHg	46
	SPECIALTY	NATURAL RUBBER 500	○ calf, thigh, thigh w/Grip-Top, thigh w/waist attachment, pantyhose	30–40 mmHg 40–50 mmHg 50–60 mmHg	50
		ADVANCE ARMSLEEVE 910	○ armsleeve w/Grip-Top w/gauntlet and without gauntlet	15–20 mmHg 20–30 mmHg	52
		INELASTIC WRAPS AND MEDICALLY COMPLEX EDEMA WEAR			54

	PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
CONSUMER HEALTH	CASUAL COTTON 146/186	● calf	15–20 mmHg	70
	TRAVENO TRAVEL SOCKS 146/T - 186/T	● calf	15–20 mmHg	71
	SEA ISLAND COTTON 151/191	● calf	15–20 mmHg	72
	ALL-SEASON MERINO WOOL 152/192	● calf	15–20 mmHg	73
	LINEN 153/193	● calf	15–20 mmHg	74
	MICROFIBER SHADES 143/183	● calf	15–20 mmHg	75
	SHEER FASHION 120	● calf, thigh w/Grip-Top, pantyhose ○ calf	15–20 mmHg	76
	SOFT SILHOUETTE LEGGINGS 170	○ open foot	15–20 mmHg	78
	ATHLETIC RECOVERY SOCKS 401	● calf	15–20 mmHg	80
	PERFORMANCE SLEEVES 412V	○ open foot	20–30 mmHg	81
	MERINO OUTDOOR SOCKS 421	● calf	15–20 mmHg	82
	EVERSOFT DIABETIC SOCKS 160	● calf	8–15 mmHg	83
RESOURCES	ACCESSORIES			84
	MARKETING TOOLS			92
	FITTER CERTIFICATION			93
	SYMPTOMS & INDICATIONS			94
	GLOSSARY OF TERMS			95
	HOW TO MEASURE			96
	HOW TO CODE			98
MOBILIS	 ORTHOPEDIC BRACING			100
	ELBOW			104
	WRIST & THUMB			107
	BACK			111
	KNEE			114
	ANKLE			118
	HOW TO ORDER MOBILIS			121

WE'VE GOT
YOU COVERED
FOR **ALL YOUR**
COMPRESSION
THERAPY NEEDS!

Our mission

Founded in 1864 and still family-owned, SIGVARIS GROUP is headquartered in Switzerland and has been one of the leaders in the global compression market for over 160 years.

Our product range spans Medical, Consumer Health Care, Complex Edema, and Soft Orthopedic Goods — all designed to support and improve quality of life. We're driven by a passion for innovation, craftsmanship, and style that meets the needs of today's customers. With our heritage, expertise, and commitment, we're ready to go even further.

Live your life.
Be yourself.

For more information,
please visit: sigvaris.com

sigvaris

We develop compression solutions designed to address a wide range of therapeutic needs and support patients in feeling their best every day.





The **DEALER PORTAL** is a secure online platform for easily ordering products and managing your accounts.

View Order History and Invoices

24/7 access – order anytime!

Sizing and Ordering

Online Chat



— TEXT US —

Same number, now with **TEXT!**

Text us directly from your phone for quick, convenient help — on your schedule.

1.800.363.4999

TO ORDER

ONLINE DEALER PORTAL
Sigvaris-online.com

PHONE OR TEXT
1.800.363.4999

EMAIL
canadaorders@sigvaris.com

Comfort & Support for Moms-to-Be



STYLE **SOFT OPAQUE**
in Black

6 in 10 women develop their first varicose vein during pregnancy. Medical compression socks help during pregnancy by **improving circulation, reducing swelling, easing leg heaviness, and lowering the risk of varicose veins or blood clots**, making legs feel more comfortable and supported.



STYLE **SHEER**
in Warm Sand




MOTION FLOW TECH®
in Cool Grey and Orchid



STYLE MICROFIBER PATTERNS
in Navy Hearts

At SIGVARIS, we recognize the strength, resilience, and long hours that define the healthcare profession.

That's why we are proud to support healthcare professionals with high-quality compression socks designed to **enhance comfort, promote well-being, and support** them through every shift.



STYLE MICROFIBER
PATTERNS in
Pink Stripe

**We support
Healthcare
Professionals!**



Travel Smart. Feel Better.

STYLE MICROFIBER PATTERNS
in Cranberry Chevron

Whether you're catching a flight or embarking on a long road trip, sitting for extended periods can slow your circulation and leave your legs feeling tired and heavy. SIGVARIS compression socks are designed to **support healthy blood flow, reduce swelling, and help prevent discomfort and travel-related leg issues** — so you arrive feeling refreshed and ready. Make compression part of your travel essentials — your legs will thank you.



MOTION FLOW TECH®
in Sage (left) and Teal (right)



ESSENTIAL OPAQUE
in Honey



MOTION CUSHIONED
COTTON in Black



ESSENTIAL OPAQUE
in Light Beige

Sigvaris compression stockings help prevent **venous disorders** from prolonged sitting or standing. By applying graduated pressure, they **improve circulation, reduce heavy legs and swelling, and enhance workplace comfort**. Regular use supports employee well-being and long-term performance, making them a valuable tool in high-risk professions.



Compression that Works for your Occupational Health

STYLE MICROFIBER
PATTERNS in Onyx Stripe



ESSENTIAL COTTON
in White



Support Yourself

Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

ESSENTIAL PRODUCTS

Series	Product	Page
230	Cotton	10
ETH	Thermoregulating Cotton	12
860	Opaque (Women)	14
860	Opaque (Men)	16

BEST SUITED FOR

- Containment of various types of edema, including venous edema, phlebolympheidema, and early stages of lymphedema
- Ideal for dress or casual shoes
- Suited for post-vein treatment



Featured colour:
WHITE

COTTON 230

COMPRESSION LEVELS

20–30 mmHg

30–40 mmHg

FEATURES AND BENEFITS



Sensitive Skin



Cotton



Breathable

- 100% Supima Cotton against the skin
- Ideal for sensitive skin and allergies
- Breathable and durable
- Excellent containment of various types of edema
- Classic ribbed design for everyday wear

FIBER CONTENT

- 50% Nylon, 25% Supima Cotton, 25% Spandex (content can vary slightly by style)
- Latex free

“

Love the cotton-lined interior
makes them incredibly comfortable
for all-day wear.

”



20–30 mmHg
Closed Toe Calf w/ Grip-Top
Medium Long in Black

232CMLM99/S



Silicone dotted **Grip-Top**
available for calf and thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
Leg Length: Women & Open Toe	Short	Long	Short	Long	Short	Long		
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")		
Leg Length: Men	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

Open Toe = ○

WHITE 00	NAVY 10	BLACK MIST 14	LIGHT BEIGE 66	BLACK 99
-------------	------------	---------------------	----------------------	-------------

Calf (C)*	●	●	●	●○	●
Calf w/Grip-Top (/S) Silicone dotted Grip-Top band				●**	●
Thigh w/Grip-Top (N) Silicone dotted Grip-Top band				●○	●

* 30–40 mmHg extra large calf only available in black with Grip-Top

** 20–30 mmHg only

Product Code Example 1 (see inset, p. 10):

232CMLM99/S

Series	23	Cotton
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Gender/Open Toe	M	For men
Colour	99	Black
Grip-Top	/S	Yes

Product Code Example 2: **232NLSW66**

Series	23	Cotton
mmHg Level	2	20-30 mmHg
Style	N	Thigh w/Grip-Top
Size	LS	Large Short
Gender/Open Toe	W	For women
Colour	66	Light Beige



Featured colour:
NATURE

THERMO-REGULATING COTTON ETH

COMPRESSION CLASS

CCL1 (20-30 mmHg)

FEATURES AND BENEFITS



Sensitive Skin



Thermoregulating



Opaque

- Made in Europe
- Fine cotton layer lies directly on the skin and has a thermoregulating effect
- Cotton absorbs excess moisture directly from the skin, stores it inside the fibers and quickly releases it through evaporation
- Ideal for people with sensitive skin
- High quality, reinforced heel
- Opaque, non-transparent look and feel

FIBER CONTENT


- 52% Polyamide, 34% Elastane, 14% Cotton



CCL1

Essential Thermoregulating; CCL1;
Calf; Size: Medium; Length: Long;
Closed Toe Black

ETH1D23199



Quadrinova knobby

Grip-Top Square nubs feel gentle
on the skin. Features a soft,
protective top seam for comfort.


“

These feel amazing — the compression is firm yet comfortable, and I can feel the support right away. They help keep my foot and ankle swelling down throughout the day, and my legs feel lighter and less sore by evening. I won't go without them now!

”

SIZING CHART

Circumference cm	XXS	XXS Plus	XS	XS Plus	S	S Plus	M	M Plus	L	L Plus
cD	23-30	27-34	27-34	31-38	31-38	34-41	34-41	37-44	37-44	41-48
cC	25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49
cB ₁	22-26	22-26	24-29	24-29	26-32	26-32	29-35	29-35	31-37	31-37
cB	18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29
cY	27-31	27-31	28-33	28-33	29-35	29-35	31-37	31-37	32-38	32-38
cA	18-20	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29



Length cm	
normal	35 – 40
long	41 – 45



How to Measure:

cB (Ankle): Circumference at the narrowest point just above the ankle bone. This is the reference point for the highest compression and is critical for sizing.

cD (Below Knee): Circumference just below the knee, at the transition from calf to knee.

lD (Length): From the floor to just below the knee crease.

Additional measurement points are only necessary for custom made circular knit products

STYLES AND COLOURS

Closed Toe = ●

	NATURE 76	BLACK 99
Calf (D) Without silicone Grip-Top band	●	●
Calf (D) w/Grip-Top (3) With silicone Grip-Top band	●	●

Product Code Example: **ETH1D23176**

Segment & Series Name	ETH	Essential Thermoregulating
Compression Level	1	CCL1
Style	D	Calf
Plus Size		P = Yes Blank = No
Ankle Size	2	9 = XX-Small 0 = X-Small 1 = Small 2 = Medium 3 = Large
Length	3	2 = Normal 3 = Long
Grip-Top		3 = With Grip-Top Blank = Without
Foot	1	1 = Closed 2 = Open
Colour	76	Nature



Featured colour:
LIGHT BEIGE

OPAQUE 860

COMPRESSION LEVELS

20-30 mmHg

30-40 mmHg

FEATURES AND BENEFITS



Durable



Largest Size Range



Opaque

- Largest size range, from petite to plus
- Durable construction combined with a smooth opaque look
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema
- Plus size pantyhose are an alternative to maternity pantyhose

FIBER CONTENT

- 64% Nylon, 36% Spandex
(content can vary slightly by style)
- Latex free

“

Really helped during my pregnancy. My legs felt supported all day (I usually wear them 8-10 hours), and the heaviness and swelling are way less than before.

”



20-30 mmHg
Closed Toe Calf Medium
Long in Black

862CMLW99



Silicone dotted **Grip-Top** available
for calf and thigh



See product images online
sigvaris.com

SIZING CHART

	Small			Medium			Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Circumference: Calf Plus w/Grip-Top	NA			NA	39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
Calf Length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")			48–73 cm (19"–28.5")			52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
Leg Length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")			87–127 cm (34.5"–50")			98–150.5 cm (38.5"–59.5")		NA	
Plus or Maternity Size Hip Circumference	NA	≤142 cm (56")		NA	≤157.5 cm (62")		≤178 cm (70")		NA	
Size	SP	SS	SL	MP	MS	ML	LS	LL	XS	XL
Calf Plus w/Grip-Top					MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Women's Shoe Size : Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	HONEY 33	GOLDEN 36	LIGHT BEIGE 66	BLACK 99
Calf (C)	●	●	○	●●
Calf Petite Sizes (C)* Not available with dotted Grip-Top band	●	●		●
Calf w/Grip-Top (C/S) Silicone dotted Grip-Top band			●○	●
Calf Plus w/Grip-Top (C/P) Only available with silicone dotted Grip-Top band			●○	●○*
Thigh w/Grip-Top (N) Silicone dotted Grip-Top band	●	●	○	●●
Thigh w/Grip-Top (N) XS & XL Silicone dotted Grip-Top band			●○	○
Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R)			○	
Thigh Petite Sizes (N)*	●	●		●
Pantyhose (P)**	●	●		●
Pantyhose Petite (P)*	●	●		●
Pantyhose Plus and Maternity (P/P)* **			●	●

* 20-30 mmHg only

** Not available in XS or XL sizing



Product Code Example: **862NLLW66**

Series	86	Opaque
mmHg Level	2	20-30 mmHg
Style	N	Thigh
Size	LL	Large Long
Gender/Open Toe	W	For women
Colour	66	Light Beige
Plus Size or Maternity with Grip-Top		If Yes = /P If No = Blank



Featured colour:
BLACK

OPAQUE 860

COMPRESSION LEVELS

20–30 mmHg

30–40 mmHg

FEATURES AND BENEFITS



Durable



Largest Size Range



Opaque

- Durable construction with a ribbed design on closed toe calf products
- Smooth opaque design on open-toe styles
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema

FIBER CONTENT

- 64% Nylon, 36% Spandex (content can vary slightly by style)
- Latex free

“

I've tried many compression socks before, but these are by far the easiest to put on. No struggling or tugging—just smooth, comfortable support all day long. They've made wearing compression part of my routine so much easier.

”



20–30 mmHg
Open Toe Thigh Medium
Long in Smooth Light Beige

862NML066



Silicone dotted **Grip-Top** available
for calf and thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Circumference: Calf Plus w/Grip-Top	NA		39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Size	SS	SL	MS	ML	LS	LL	XS	XL
Calf Plus w/Grip-Top			MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	RIBBED LIGHT BEIGE 66	RIBBED BLACK 99	SMOOTH LIGHT BEIGE 66	SMOOTH BLACK 99
Calf (C)	●*	●*	○	○
Calf w/Grip-Top (C/S) Silicone dotted Grip-Top band			○	
Calf Plus w/Grip-Top (C/P) Only available with silicone dotted Grip-Top band		●*	○	○*
Open Toe Thigh w/Grip-Top (N) Silicone dotted Grip-Top band			○	○
Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R)			○	

* 20-30 mmHg only

** Not available in XS or XL with waist attachment

Product Code Example: **862CMLM99**

Series	86	Opaque
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Gender/Open Toe	M	For men
Colour	99	Ribbed Black
Grip-Top		If yes = /S If no = blank



STYLE **LINEN** in Denim and
STYLE **MICROFIBER PATTERNS**
in Navy Hearts



Express Yourself

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

STYLE PRODUCTS

	Series	Product	Page
NATURAL FIBER	220	Sea Island Cotton	20
	240	Merino Wool	22
	250	Linen	24
PLANT BASED FIBER	260	Bamboo	26
SHEER	750	Medium Sheer	28
	780	Sheer	30
MICROFIBER	820	Microfiber	32
	830	Microfiber Patterns	34
	840	Soft Opaque	36

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- Pregnancy
- People who prefer fashionable options
- Ideal for business or dress attire



Featured colour:
NAVY

SEA ISLAND COTTON 220

COMPRESSION LEVEL

20–30 mmHg

Also available in 15–20 mmHg, see page 72

FEATURES AND BENEFITS



Soft



Cotton



Breathable

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

FIBER CONTENT

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free



20–30 mmHg
Closed Toe Calf Small Long
in Brown for women

222CSLW11



20–30 mmHg
Closed Toe Calf Medium Long
in Black for men

222CMLM99

“

These have helped with the pain and discomfort of swelling in my lower leg following an injury and surgeries.

”



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

	NAVY 10	BROWN 11	BLACK 99
Calf (C)	●	●	●



Product Code Example: **222CMLM10**

Series	22	Sea Island Cotton
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Gender	M	For men
Colour	10	Navy



Featured colour:
BROWN

MERINO WOOL 240

COMPRESSION LEVEL

20–30 mmHg

Also available in 15–20 mmHg, see page 73

FEATURES AND BENEFITS



Wool



Thermoregulating



Stylish

- High quality extra-fine Merino wool makes this sock soft and breathable
- Moisture wicking capabilities keep feet dry
- Naturally odour resistant
- Thermoregulating for year-round wear
- Smooth, flat toe seam

FIBER CONTENT

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free

“

They are light weight so I can wear them anytime. My legs feel better wearing them and my feet stay comfortable. I love the comfort in the toe area.

”



20–30 mmHg
Closed Toe Calf Medium
Long in Brown for men

242CMLM11



20–30 mmHg
Closed Toe Calf Small
Long in Navy for women

242CSLW10



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

	NAVY 10	BROWN 11	BLACK 99
Calf (C)	●	●	●



Product Code Example: **242CMSW99**

Series	24	Merino Wool
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	MS	Medium Short
Gender	W	For women
Colour	99	Black



Featured colour:
BROWN

LINEN 250

COMPRESSION LEVEL

20–30 mmHg

Also available in 15–20 mmHg, see page 74

FEATURES AND BENEFITS



Hypoallergenic
Natural Fibers



Anti-bacterial



Moisture
Management

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

FIBER CONTENT

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free

New
Colours
coming in
2026!



20–30 mmHg
Closed Toe Calf Medium
Long in Light Grey for men

252CMLM205



20–30 mmHg
Closed Toe Calf Medium
Long in Denim for men

252CMLM201

“ I like the material as it is hypoallergenic and breathable. Considerably easier to put on. ”



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12
Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

				
	BROWN 111	DENIM 201	LIGHT GREY 205	LAVENDER 207
Closed Toe = ●	●	●	●	●
Calf (C)	●	●	●	●
Gender	For women & men			For women



Product Code Example: **252CMLW207**

Series	25	Linen
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Gender	W	For women
Colour	207	Lavender



Featured colours:
IVORY and DUSK VIOLET

NEW Soft & Breathable BAMBOO 260

COMPRESSION LEVEL

20-30 mmHg

FEATURES AND BENEFITS



Soft



Stylish



Breathable

- Bamboo viscose yarns provide excellent breathability and moisture management for added comfort
- Ideal for year-round wear - keeps feet and legs cool in summer and warm in winter
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Incredibly soft and lightweight socks crafted from bamboo viscose yarns in four colors makes legs feel and look great
- 100% bamboo viscose against the skin for an ultra-soft feel
- Aesthetically pleasing top-of-foot design

FIBER CONTENT

- 42% Viscose from Bamboo, 40% Nylon, 18% Spandex
- Latex free

Coming
summer
2026!

UNISEX SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12

Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

STYLES AND COLOURS

Closed Toe = ●

	IVORY 217	STONE BLUE 218	DUSK VIOLET 219	BLACK 199
Calf (C)	●	●	●	●



20–30 mmHg
Closed Toe Calf in Ivory

262CML217



Aesthetically-pleasing
top-of-foot design



Product Code Example: **262CML219**

Series	26	Bamboo
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Colour	219	Dusk Violet



Featured colour:
HONEY

MEDIUM SHEER 750

COMPRESSION LEVEL

20-30 mmHg

FEATURES AND BENEFITS



Medium Sheer



Durable



Breathable

- Fashionable, durable, medium sheer hosiery for everyday wear
- Excellent breathability and moisture management
- Comfortable 2" top band on calf style

FIBER CONTENT

- 68% Nylon, 32% Spandex
(content can vary slightly by style)
- Latex free

“

I spend all day in the office, and they make my tired and achy legs feel better.

”



20-30 mmHg
Closed Toe Thigh Medium
Long in Black

752NMLW99



Silicone dotted **Grip-Top**
on thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●

	HONEY 33	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	●	●
Thigh w/Grip-Top (N) Silicone dotted Grip-Top band	●	●	●
Pantyhose (P)	●	●	●



Product Code Example: **752NMLW33**

Series	75	Medium Sheer
mmHg Level	2	20-30 mmHg
Style	N	Thigh (w/Grip-Top)
Size	ML	Medium Long
Gender/	W	For women
Colour	33	Honey



SHEER 780

COMPRESSION LEVELS

15–20 mmHg

20–30 mmHg

30–40 mmHg

FEATURES AND BENEFITS



Sheer



Elegant



Breathable

- Fashionable for professional and dress wear
- Consistent sheerness throughout the leg
- Excellent breathability and moisture management

FIBER CONTENT

- 15–20 and 20–30 mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30–40 mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free

“


I really like the color, as it matches my skin tone.

”



20–30 mmHg
Closed Toe Calf Small Short
in Café

782CSSW73



Sensinnov® Grip-Top
on thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18-21.5 cm (7"-8.5")		21.5-25.5 cm (8.5"-10")		25.5-29 cm (10"-11.5")	
Calf Circumference	28-39 cm (11"-15.5")		31.5-45 cm (12.5"-17.5")		35.5-51 cm (14"-20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5-66 cm (18"-26")		48-73 cm (19"-28.5")		52-84 cm (20.5"-33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83-112.5 cm (32.5"-44.5")		87-127 cm (34.5"-50")		98-150.5 cm (38.5"-59.5")	
Size	SS	SL	MS	ML	LS	LL


Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	DARK NAVY 08	HONEY 33	WARM SAND 34*	GOLDEN 36	TOASTED ALMOND 37*	CAFÉ 73*	MOCHA 85*	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	●○	●○	●○	●○	●○	●○	●	●○
Thigh w/Grip-Top (N) Sensinnov® Grip-Top band	●	●○	●○	●○	●○	●○	●○	●	●○
Pantyhose (P)	●	●○	●○	●○	●○	●○	●○	●	●○

* 20-30 mmHg only



Product Code Example: **782PMLW99**

Series	78	Sheer
mmHg Level	2	20-30 mmHg
Style	P	Pantyhose
Size	ML	Medium Long
Gender/Open Toe	W	For women
Colour	99	Black



Featured colour:
STEEL GREY

MICROFIBER 820

COMPRESSION LEVELS

15–20 mmHg

20–30 mmHg

30–40 mmHg

FEATURES AND BENEFITS



Microfiber



Soft



Stylish

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking
- Glides on and off with ease
- Comfortable and durable sock for professional and everyday wear
- Available in extended foot sizes
- Smooth, flat toe seam

FIBER CONTENT

- 87% Nylon , 13% Spandex
(content can vary slightly by style)
- Latex free




20–30 mmHg
Closed Toe Thigh Medium
Long in Black

822NMLM99



Silicone dotted **Grip-Top** for
calf and **Sensinno® Grip-Top**
for thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
Leg Length	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
Size Shoe Size	SS 7–11	SL 9.5–14	MS 7–11	ML 9.5–14	LS 7–11	LL 9.5–14	XS 7–11	XL 9.5–14
				MX* 12.5–17.5		LX* 12.5–17.5		XX* 12.5–17.5

Men's shoe sizes: **Short Leg 7–11, Long Leg 9.5–14** (For extended foot sizes: long leg 12.5–17.5)

*MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30 mmHg only. Extended foot sizes not available with Grip-Top.

STYLES AND COLOURS

	TAN KHAKI 32	STEEL GREY 93	BLACK 99
Closed Toe = ●	●	●***	●
Calf (C)	●	●***	●
Calf w/Grip-Top (C/S)** Silicone dotted Grip-Top band	●		●
Thigh w/Grip-Top (N)** Sensinnov® Grip-Top band	●		●

** Not available in 15–20 mmHg

*** Steel Grey only available in 15–20 mmHg and 20–30 mmHg



Product Code Example: **822CMLM93**

Series	82	Microfiber
mmHg Level	2	20–30 mmHg
Style	C	Calf
Size	ML	Medium Long
Gender	M	For men
Colour	93	Steel Grey
Grip-Top		No



Featured colour:
NAVY HEARTS


MICROFIBER PATTERNS 830

COMPRESSION LEVEL

20–30 mmHg

Also available in 15–20 mmHg, see page 75

FEATURES AND BENEFITS



Microfiber



Stylish



Patterns

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- Fashionable, durable sock for everyday wear
- Glides on and off with ease
- Smooth, flat toe seam

FIBER CONTENT

- 87% Nylon, 13% Spandex
(content can vary slightly by style)
- Latex free



I really love the variety of patterns you have available! It makes wearing the socks more fun!



20–30 mmHg

Closed Toe Calf Small Long
in Royal Blue Argyle for men

832CSLM48



20–30 mmHg

Closed Toe Calf Medium Short in
Wisteria Windows for women

832CMSW61



See product images online
sigvaris.com

SIZING CHART






	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●
For women = **W**
For men = **M**

				
Calf (C)	●	●	●	●
Gender	W	W	W/M	W
				
Calf (C)	●	●	●	●
Gender	M	W/M	W	W
				
Calf (C)	●	●	●	●
Gender	W/M	W/M	M	W
				
Calf (C)	●	●	●	●
Gender	W/M	W/M	W	W/M



Product Code Example: **832CMLM42**

Series	83	Microfiber Patterns
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Gender	M	For men
Colour	42	Black Plaid



Featured colour:
CHAI

SOFT OPAQUE 840

COMPRESSION LEVELS

15–20 mmHg

20–30 mmHg

30–40 mmHg

FEATURES AND BENEFITS



Opaque



Microfiber



Soft

- Fashionable opaque hosiery for professional and everyday wear
- Long-lasting microfiber softness provides all day comfort and allows product to easily glide on and off
- Durable and breathable
- Maternity pantyhose option available

FIBER CONTENT

- 79% Nylon, 21% Spandex
(content can vary slightly by style)
- Latex free

“


I LOVE these stockings! I am currently 33 weeks pregnant and swelling like crazy! I got the maternity style ones right now so they are perfect for when you have a bump in the way! My legs thank you so much!!!
Such a pain relief!

”



20–30 mmHg
Open Toe Thigh with Grip-Top
Medium Long in Chai

842NMLO35



Silicone dotted **Grip-Top** for calf and **Sensinno® Grip-Top** for thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Maternity hip circumference	≤142 cm (56")		≤157.5 cm (62")		≤178 cm (70")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	MIDNIGHT BLUE 09**	CHAI 35	PECAN 40	MULBERRY 60*	ESPRESSO 89*	GRAPHITE 91**	BLACK 99
Calf (C)	●	●○	●	●	●	●	●○
Calf w/Grip-Top (C/S)* Silicone dotted Grip-Top							●
Thigh w/Grip-Top (N) Sensinnov® Grip-Top band	●	●○	●	●	●	●	●○
Pantyhose (P)	●	●○	●	●	●	●	●○
Maternity Pantyhose (M)*							●

* 20–30 mmHg only

** 15–20 mmHg and 20–30 mmHg only



Product Code Example: **842CMSW09**

Series	84	Soft Opaque
mmHg Level	2	20–30 mmHg
Style	C	Calf
Size	MS	Medium Short
Gender/Open Toe	W	For women
Colour	09	Midnight Blue
Grip-Top		No



MOTION FLOW TECH®
in Sky Blue and Slate



Push Yourself

Why choose **Motion**?

- Cushioned sole for added comfort
- Moisture-wicking yarns

MOTION PRODUCTS

Series	Product	Page
360	Cushioned Cotton	40
412C	High Tech	42
422	Thermoregulating Wool	44
432	Flow Tech®	46

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Casual and occupational footwear
- Walking and everyday activities





Featured colour:
WHITE

CUSHIONED COTTON 360

COMPRESSION LEVEL

20-30 mmHg

FEATURES AND BENEFITS



Cushioned Sole



Cotton



Sensitive Skin

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Smooth, flat toe seam
- Diabetic friendly
- Ideal for people who stand for prolonged periods of time

FIBER CONTENT

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free



20-30 mmHg
Closed Toe Calf Medium
Long in Black for women

362CMLW99



Flat, low friction toe seam



The socks are very comfortable, easy to put on, wash up well and keep the fit through the day.



See product images online
sigvaris.com

SIZING CHART


	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

	 WHITE 00	 BLACK 99
Calf (C)	●	●



Product Code Example: **362CMLM00**

Series	36	Cushioned Cotton
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Gender	M	For men
Colour	00	White



Featured colour:
STEEL BLUE

HIGH TECH 412C

COMPRESSION LEVEL

20-30 mmHg

FEATURES AND BENEFITS



Moisture
Wicking



Anatomically
Contoured



Colours

- High-Tech fibers for moisture wicking and improved thermoregulation
- Cushioned footbed, toe and Achilles tendon zones for optimal comfort
- Anatomically contoured to support the natural shape of the foot and calf (left and right foot specific)
- Smooth, flat toe seam
- Perfect for everyday activities and active individuals

FIBER CONTENT

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free

“

They stay in place all day without slipping, thicker material is great for steel toe work boots in 12 hour shifts. Great moisture wicking. This product is meant for construction workers and built to last!

”



20-30 mmHg
Closed Toe Calf Medium
Medium in White

412CMM00



LEFT and RIGHT
foot specific



See product images online
sigvaris.com

UNISEX SIZING CHART

	Small			Medium			Large			X Large		
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")			29–33 cm (11.5"–13")		
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")			40.5–58 cm (16"–23")		
	Short	Medium	Long	Short	Medium	Long	Short	Medium	Long	Short	Medium	Long
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

STYLES AND COLOURS

Closed Toe = ●	WHITE 00	BLUE 50	LIME 54	PINK 56
Calf (C)	●	●	●	●
	RED 59	STEEL BLUE 62	BLACK 99	
Calf (C)	●	●	●	

Product Code Example: **412CML56**

Series	41	High Tech
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Colour	56	Pink



THERMO-REGULATING WOOL 422

COMPRESSION LEVEL

20–30 mmHg

Also available in 15–20 mmHg, see page 82

FEATURES AND BENEFITS



Wool



Cushioned Sole



Thermoregulating

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned soles help prevent blisters
- Soft and comfortable
- Smooth, flat toe seam
- Perfect for outdoor work environments and activities

FIBER CONTENT

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free

“

These have been a life changer. At my age with 3 year old twins, these compression socks give me a feeling of hope and ease. Knowing my legs won't be affected while seeing my kids grow means the world to me.

”



See product images online
sigvaris.com



20–30 mmHg
Closed Toe Calf Medium
Medium in Charcoal

422CMM12



20–30 mmHg
Closed Toe Calf Medium
Short in Olive

422CMS31

UNISEX SIZING CHART

	Small			Medium			Large		
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		
	Short	Medium	Long	Short	Medium	Long	Short	Medium	Long
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

STYLES AND COLOURS

Closed Toe = ●

	CHARCOAL 12	OLIVE 31
Calf (C)	●	●

Product Code Example: **422CSL12**

Series	42	Thermoregulating Wool
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	SL	Small Long
Colour	12	Charcoal



Featured colour:
SLATE

FLOW TECH® 432

COMPRESSION LEVEL

20–30 mmHg

FEATURES AND BENEFITS



Breathable



Moisture Wicking



Cushioned Sole

- MeshLink™ channels are designed to deliver optimal airflow to the parts of the foot and leg that generate the most heat, linking the bottom of the foot to the top of the calf and allowing heat to escape
- Cushioned footbed and toe with higher Achilles tendon zones for optimal comfort
- SENSIL® Aquarius fibers keep your legs and feet dry with next-generation moisture management
- Mesh flex zone across the instep for donning ease
- Smooth, flat toe seam


FIBER CONTENT

- 62% Nylon, 24% SENSIL Aquarius Nylon, 14% Spandex
- Latex free



20–30 mmHg
Closed toe Calf Medium
Long in Deep Sea

432CML212



MeshLink™ channels,
our unique pattern design
for optimal airflow

“The moment I put them on, I could feel a cooling sensation, which I recognized as increased blood flow and circulation. Within 30 minutes my feet and legs had returned to normal size and I was able to go for a walk for the first time in almost a year without any pain.”



See product images online
sigvaris.com



UNISEX SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●	 BLACK 199	 ORCHID 210	 COOL GREY 211	 DEEP SEA 212
Calf (C)	●	●	●	●
	 SAGE 213	 TEAL 214	 SLATE 215	 SKY BLUE 216
Calf (C)	●	●	●	●



Product Code Example: **432CML210**

Series	43	Flow Tech®
mmHg Level	2	20-30 mmHg
Style	C	Calf
Size	ML	Medium Long
Colour	210	Orchid





Take care of yourself

Why choose Specialty?

- Specific medical products
- High content

SPECIALTY PRODUCTS

Serie	Produit	Page
500	Calatone Pulver	50
910	Advan Advanced	52

BEST SUITED FOR

- Severe varicose veins
- Chronic venous insufficiency
- Severe venous stasis
- Early stages of lymphedema
- Active venous ulcers and prevention of recurrence
- Patients with early stages of lymphedema
- Post-operative care



Featured colour:
BEIGE

NATURAL RUBBER 500

COMPRESSION LEVELS

30–40 mmHg

40–50 mmHg

50–60 mmHg

FEATURES AND BENEFITS



Natural Rubber



Durable



Largest Size Range

- Durable, long-term containment solution with minimal fiber fatigue throughout the day
- Extended sizing for a custom-like fit
- Can be paired with **Cotton underliner** (see page 87)
- Ideal for patients with more serious indications, including:
 - Severe varicosities
 - Chronic edema post DVT
 - Severe venous stasis
 - Prevention and treatment of Active Venous Ulcers
 - Lymphedema

FIBER CONTENT


- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)

Caution: contains natural latex rubber that may cause allergic reaction.



30–40 mmHg
Open Toe Thigh Medium
Average Long in Beige

503NM2077



Silicone dotted **Grip-Top** available
on thigh



Cotton underliner
Sold separately,
see page 87



See product
images online
sigvaris.com

SIZING CHART

	Small				Medium				Large				X Large			
Ankle Circumference	19–23 cm (7.5"–9")				23–30 cm (9"–12")				30–40 cm (12"–15.5")				40–46 cm (15.5"–18")			
Calf Circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	28–38 cm (11"–15")		≤43 cm (17")		32–43 cm (12.5"–17")		≤48 cm (19")		36–48 cm (14.5"–19")		≤57 cm (22.5")		45–57 cm (18"–22.5")		≤68 cm (27")	
Calf Circumference for 505C (Full only)	505C: full calf only		32–43 cm (12.5"–17")		505C: full calf only		36–48 cm (14.5"–19")		505C: full calf only		45–57 cm (18"–22.5")		505C: not available in XL			
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")
Thigh Circumference	Average		Full		Average		Full		Average		Full		N/A			
	41–54 cm (16"–21")		≤61 cm (24")		48–61 cm (19"–24")		≤71 cm (28")		58–71 cm (23"–28")		≤81 cm (32")					
Thigh Circumference for 505T (Full only)	505T: full thigh only		≤61 cm (24")		505T: full thigh only		≤71 cm (28")		505T: full thigh only		≤81 cm (32")					
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")				
Hip Circumference	<122 cm (48")				<137 cm (54")				<155 cm (61")							
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

STYLES AND COLOURS

Open Toe = ○



Calf (C)	○	Calf only available in S, M, and L full circumference for 50–60 mmHg
Thigh w/out Grip-Top (T)	○	
Thigh w/Grip-Top (N) Silicone dotted Grip-Top band	○	Not available in 50–60 mmHg or XL
Pantyhose (P)	○	Pantyhose only available in 30–40 mmHg and only in full thigh circumference
Thigh w/Waist Attachment (W) Indicate Leg LEFT (/L) or RIGHT (/R)	○	Thigh w/waist attachment is available in average and full thigh circumferences in 30–40 mmHg and in full thigh circumferences only in 40–50 mmHg

Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 94 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions.



Product Code Example: **503NS3077**

Series	50	Natural Rubber
mmHg Level	3	30–40 mmHg
Style	N	Thigh w/Grip-Top
Size	S3	Small Full Short
Gender/Open Toe	O	Open Toe
Colour	77	Beige



Featured colour:
BEIGE

ADVANCE ARMSLEEVE 910

COMPRESSION LEVELS

15–20 mmHg

20–30 mmHg

FEATURES AND BENEFITS



**Moisture
Wicking**



**Sensinno®
Grip-Top Band**



Breathable

- Ideal for patients with arm lymphedema or post-mastectomy care
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable and soft
- Sensinno® grip-top band ensures armsleeve stays in place all day long
- Available with or without gauntlet

FIBER CONTENT

- 50% Nylon, 28% Rayon, 22% Spandex
- Latex free



20–30 mmHg
Armsleeve With Gauntlet
Small Long in Beige

912ASLG77/S




Sensinno® Grip-Top band

“
Supportive, easy to wear, and a big help
after my mastectomy.
”

SIZING CHART

	Small		Medium				Large				X Large			
Palm (A)	17–21 cm (6.5"–8")		18–22 cm (7"–8.5")				19–23 cm (7.5"–9")				20–24 cm (8"–9.5")			
Wrist (C)	14–16 cm (5.5"–6.5")		16–18 cm (6.5"–7")				18–20 cm (7"–8")				20–23 cm (8"–9")			
Upper Mid-Arm (F)	Regular		Regular		Plus		Regular		Plus		Regular		Plus	
	26–34 cm (10"–13.5")		29–37 cm (11.5"–14.5")		37–43 cm (14.5"–17")		32–40 cm (12.5"–15.5")		40–46 cm (15.5"–18")		34–42 cm (13.5"–16.5")		42–48 cm (16.5"–19")	
Arm Length (C-G)*	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")
Size	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P

* Should have slight bend in arm when measuring length



STYLES AND COLOURS

Open Sleeve = ○



Armsleeve with Gauntlet (G)

Sensinnov® Grip-Top band

○

Armsleeve without Gauntlet (O)

Sensinnov® Grip-Top band

○

Note: all armsleeves come with Grip-Top (/S)
15–20 mmHg not available with gauntlet or plus sizing



Product Code Example: **912AMRO77/S**

Series	91	Advance Armsleeve
mmHg Level	2	20–30 mmHg
Style	A	Armsleeve
Size	MR	Medium Regular
Gauntlet	O	Without Gauntlet
Colour	77	Beige
Grip-Top	/S	With Grip-Top



INELASTIC WRAPS AND MEDICALLY COMPLEX EDEMA WEAR




INELASTIC WRAPS

Product	Page
Why wraps?	56
Which wrap do I choose?	57
Compreflex®	58
Coolflex®	60
Foot Sizing Charts	62
Foot Compression Garments	63
Liners	64
Accessories	65
Compreshort®	66
Compreshape™	67

Why wraps?

- Faster & Easier to apply
- No multilayer wrapping technique needed
- Consistent compression levels
- Easier to achieve and teach to patients/caregivers
- Time-saving
- Less need for frequent re-wrapping
- Cost-effective over time, because inelastic wraps are reusable



How Inelastic Wraps Work

EFFECTIVE WORKING PRESSURE

• Low Resting Pressure:

When the patient is at rest (sitting or lying down), the garment exerts relatively little pressure because the fabric has very little stretch. This makes it more comfortable to wear for long periods, especially for those who are less mobile or sleeping.


• High Working Pressure:

When the muscles contract during movement (such as walking), they push outward against the firm, inelastic material. Because the garment doesn't stretch much, it resists that expansion and creates a surge of pressure against the limb.

Working Pressure

Walking, stair climbing, sit to stand, toe raises

Resting Pressure
Donning Compreflex®



Resting Pressure
Laying Down, Sitting Still

Which wrap do I choose?

COMPREFLEX® Calf



Accutab® Devices



Accutabs® – Patented and adjustable pressure monitoring tabs for easy, accurate and consistent control

Complete Leg System – Option to order **Compreflex® Standard Knee** and **Compreflex® Standard Thigh**

Strap Extenders – Option to add 10cm per strap when needed

Premium materials – Soft, breathable and offers a comfortable fit.



COOLFLEX® Standard - Calf



EZ-Close™ Pressure System

One-Handed Donning – **EZ-Close™** Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues.

Stay cool – Lightweight, breathable mesh keeps skin comfortable.

Low-profile – Slim design fits easily under clothing and with footwear.

Moisture regulating – Mesh textile allows wearer to feel cool and comfortable.



COMPREFLEX®

COMPRESSION LEVEL

20–50 mmHg



PREMIUM MATERIALS

Soft, conforming material provides a comfortable fit

INTELLIGENT DESIGN


Garment features a contoured calf, providing a superior fit

SUPERIOR QUALITY

Each garment is double stitched to ensure durability

EASY APPLICATION

With a front stretch panel to secure the garment in place and straps that roll back, the garment is easy to don and adjust, improving compliance



AccuTabs® allows the patient and HCP to accurately measure the compression required

FEATURES AND BENEFITS



Easy donning
& doffing



Containment




Durable

FIBER CONTENT

- 37% Polyamide, 27% Polyurethane Foam, 26% Cotton, 6% Elastane, 4% Polyester
- Latex free

IDEAL FOR

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including ulcer due to venous stasis
- Lipedema



COMPREFLEX® CALF SIZING CHART (A, C)

	Small	Medium	Large	X Large	XX Large
C	29–39 cm	34–44 cm	40–50 cm	46–56 cm	54–64 cm
A	19–26 cm	21–30 cm	26–36 cm	31–41 cm	36–46 cm

G Length:

Regular: 30–33.9 cm

Tall: 34–42 cm

Foot Size

I: Measure circumference of ankle bend and heel

J: Measure circumference across metatarsal heads

K: Measure **medial** length from heel to 1st metatarsal head
(See page 62 for foot sizing chart)

COMPREFLEX® OPTIONS

CALF & FOOT


20–50 mmHg



For those who require a calf and foot garment

Includes:

- Compreflex® Wrap
- Compreboot® Standard
- One pair of non-compressive Basic Liners
- One set of Accutab® Devices



PRODUCT CODES

		Small	Medium	Large	X Large	XX Large
Beige	Reg	1411-BKR	1412-BKR	1413-BKR	1414-BKR	1415-BKR
	Tall	1411-BKT	1412-BKT	1413-BKT	1414-BKT	1415-BKT
Black	Reg	1401-BKR	1402-BKR	1403-BKR	1404-BKR	1405-BKR
	Tall	1401-BKT	1402-BKT	1403-BKT	1404-BKT	1405-BKT

- 1401-BKR, 1411-BKR, 1401-BKT, and 1411-BKT include a Regular Length Boot
- 1402-BKR, 1412-BKR, 1402-BKT, 1412-BKT, 1403-BKR, 1413-BKR, 1403-BKT, and 1413-BKT include a Medium/Large Regular Length Boot
- 1404-BKR, 1414-BKR, 1404-BKT, 1414-BKT, 1405-BKR, 1415-BKR, 1405-BKT, and 1415-BKT include a X-Large/XX-Large Regular Length Boot

Standard CALF

20–50 mmHg



For those who require an alternative foot garment

Includes:

- Compreflex® Wrap
- One pair of non-compressive Basic Liners
- One set of Accutab® Devices
- Can be paired with alternative foot pieces* (p 63)

* When paired with a basic liner, a foot piece is highly recommended



		Small	Medium	Large	X Large	XX Large
Beige	Reg	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
	Tall	1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT
Black	Reg	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
	Tall	1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT

Transition CALF

20–50 mmHg



For those who require a compressive liner for the foot and ankle

Includes:

- Compreflex® Wrap
- One pair of Transition Liners, which provide 10-20 mmHg in the foot and ankle
- One set of Accutab® Devices
- Ability to wear regular shoe over liner



		Small	Medium	Large	X Large	XX Large
Beige	Reg	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
	Tall	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT
Black	Reg	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
	Tall	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT

Complete CALF

20–50 mmHg



For those who require a compressive liner with compression in the foot and ankle, and graduated compression in the calf

Includes:

- Compreflex® Wrap
- One pair of Complete Liners, which provide 10-20 mmHg in the foot and ankle, and graduated compression in the calf
- One set of Accutab® Devices
- Ability to wear regular shoe over liner



		Small	Medium	Large	X Large	XX Large
Beige	Reg	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
	Tall	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT
Black	Reg	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
	Tall	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT

COOLFLEX®

COMPRESSION LEVEL

20-50 mmHg




ONE-HANDED DONNING

EZ-Close™ Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues

BREATHABLE GARMENT

The Spacer Fabric™ is designed to be highly breathable for cool-wearing comfort



INNOVATIVE DESIGN

Durable, metal hook closures for maximum strength and durability

Left and Right Specific

FEATURES AND BENEFITS



Containment



Breathable



Moisture Wicking

FIBER CONTENT

- 47% Polyester, 25% Polyamide, 20% Cotton, 8% Elastane, 1% Stainless Steel
- Latex free


IDEAL FOR

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Well-managed Lymphedema
- Non-specific or post-surgical edema
- Patients with limited dexterity, arthritis, and mobility limitations

HOW COOLFLEX® IS PROVEN TO ACCURATELY SET COMPRESSION LEVELS

EZ-Close™

The integrated **EZ-Close™** Pressure System has visual indicators for the accurate delivery of 20-30 mmHg, 30-40 mmHg, or 40-50 mmHg on each tab.



COOLFLEX® STANDARD

CALF

20–50 mmHg



INCLUDES


- Coolflex® Standard Calf
- Built-in application sleeve for the independent use and wear
- **One pair of non-compressive Basic Liners**
- You can substitute the **Basic Liner** for a **Transition Liner** (see page 64) by contacting **customer care**.



Basic Liner (included)

When pairing with a **Basic Liner**, a foot piece is highly recommended.

See **Coolflex® Foot Piece** on page 63; sold separately.




Transition Liner

When pairing with a **Transition Liner**, a foot piece may not be required.



COOLFLEX® STANDARD CALF SIZE CHART (A–C1)

		Small	Small Max	Medium / Large	Medium / Large Max	X Large / XX Large
	C1	28–39 cm	35–50 cm	34–50 cm	45–64 cm	45–64 cm
	C	28–39 cm	35–50 cm	34–50 cm	45–64 cm	45–64 cm
	B	23–33 cm	27–38 cm	29–43 cm	35–49 cm	37–54 cm
	A	16–26 cm	16–26 cm	23–33 cm	23–33 cm	30–42 cm
Left	Regular	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
	Tall	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
Right	Regular	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
	Tall	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R




G Length:

Regular: 30–35.9 cm

Tall: 36–42 cm

FOOT SIZING CHARTS



Foot Size:

I: Measure circumference of ankle bend and heel

J: Measure circumference across metatarsal heads

K: Measure **medial** length from heel to 1st metatarsal head

Gripper Patches

Available for **Compreboot® Standard FOOT** and **Compreboot® Plus FOOT**



COMPREBOOT® LITE (I, J, K)

	Small / Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long
I	28–36 cm	28–36 cm	≤44 cm	≤44 cm
J	22–26 cm	22–26 cm	≤33 cm	≤33 cm
K	14–18 cm	19–23 cm	20–22 cm	25–27 cm

PRODUCT CODE

1st set	1800-GP
---------	---------

COMPREBOOT® STANDARD AND PLUS (I, J, K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36cm	28–36 cm	≤39 cm	≤39 cm	≤44 cm	<44 cm
J	22–26 cm	22–26 cm	≤30 cm	≤30 cm	≤33 cm	≤33 cm
K	14–18 cm	19–23 cm	18–20 cm	23–25 cm	20–22 cm	25–27 cm

COOLFLEX® STANDARD FOOT (I, J, K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–39 cm	N/A	30–43 cm	N/A	34–49 cm	N/A
J	22–26 cm	N/A	25–30 cm	N/A	28–33 cm	N/A
K	14–18 cm	19–21 cm	18–20 cm	20–22 cm	20–22 cm	23–25 cm

FOOT COMPRESSION GARMENTS

Foot compression garments below are pictured with the Basic Liner (not included). Basic Liner sold separately (see page 64) or with the purchase of an Inelastic Wrap.

Compreboot® Lite FOOT

Variable compression



- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design may be trimmed to fit
- Provides mild, therapeutic compression



PRODUCT CODE (Available in Black Only)

		Small / Medium / Large	X Large / XX Large
Black	Regular	1803-BLR	1805-BLR
	Long	1803-BLL	1805-BLL

Compreboot® Standard FOOT

Variable compression



- Four-way stretch for firm compression around ankle
- Ideal for mild to moderate chronic edema
- Low stretch material conforms to foot shape



PRODUCT CODE

		Small	Medium / Large	X Large / XX Large
Beige	Regular	1811-BTR	1813-BTR	1815-BTR
	Long	1811-BTL	1813-BTL	1815-BTL
Black	Regular	1801-BTR	1803-BTR	1805-BTR
	Long	1801-BTL	1803-BTL	1805-BTL

Compreboot® Plus FOOT

Variable compression



- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate, therapeutic compression
- Used for Lymphedema stage I, II, III



PRODUCT CODE (Available in Black Only)

		Small	Medium / Large	X Large / XX Large
Black	Regular	1801-BPR	1803-BPR	1805-BPR
	Long	1801-BPL	1803-BPL	1805-BPL

Coolflex® Standard FOOT

Variable compression



- Low profile design allows boot to fit in most footwear
- Spacer Fabric™ is highly breathable for cool-wearing comfort
- Durable metal clasps, instead of Velcro closure



PRODUCT CODE (Available in Black Only)

		Small	Medium / Large	X Large / XX Large
Right (Black)	Regular	2101-BTR-R	2103-BTR-R	2105-BTR-R
	Long	2101-BTL-R	2103-BTL-R	2105-BTL-R
Left (Black)	Regular	2101-BTR-L	2103-BTR-L	2105-BTR-L
	Long	2101-BTL-L	2103-BTL-L	2105-BTL-L

LINER OPTIONS

Basic Liner

No compression

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- drirelease® moisture wicking and odour control
- Flat seams for added comfort
- Available in Heathered Silver colour only



PRODUCT CODE

(Included with Compreflex® Standard Calf & Foot, Compreflex® Standard Calf, and Coolflex® Standard Calf)

Size	Pair (gray)	10 pairs (gray)
One size	1700 - NK1	1700 - NKC

Transition Liner

10–20 mmHg

- Mild compression in foot and ankle only
- drirelease® moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Available in Black and Beige



PRODUCT CODE

(Included with Compreflex® Transition)

Size	Pair (black)	Pair (beige)
Small	1601-CK	1611-CK
Medium	1602-CK	1612-CK
Large	1603-CK	1613-CK
X Large	1604-CK	1614-CK
XX Large	1605-CK	1615-CK

* When ordering a Coolflex® standard without a foot piece, please order a transition liner.

Complete Liner

10–20 mmHg

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- drirelease® moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Available in Black and Beige




PRODUCT CODE

(Included with Compreflex® Complete)

Size	Pair (black)	Pair (beige)
Small	1601-CK-F1	1611-CK-F1
Medium	1602-CK-F1	1612-CK-F1
Large	1603-CK-F1	1613-CK-F1
X Large	1604-CK-F1	1614-CK-F1
XX Large	1605-CK-F1	1615-CK-F1

TRANSITION AND COMPLETE LINERS SIZING CHART

Size	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe
Small	19–26 cm	29–39 cm	<7½	<6½
Medium	21–30 cm	34–44 cm	7.5–9.5	6.5–8.5
Large	26–36 cm	40–50 cm	10–12	9–11
X Large	31–41 cm	46–56 cm	12.5–15	11.5–14
XX Large	36–46 cm	54–64 cm	12.5–15	11.5–14



ACCESSORIES

Compreflex® Strap Extender

- Provides up to 10 cm of additional circumferential size range to **Compreflex®** calf wraps
- Detaches easily
- Potential solution for patients outside of standard size ranges of **Compreflex®** calf wraps or needing custom size
- Max 1 strap extender per strap

PRODUCT CODE		
Qty	Black	Beige
1 unit	1401-SE1	1411-SE1
10 Unit Pack	1401-SEB	1411-SEB



Accutab® Devices



Easy to use

- Position the **Compreflex®** garment on the limb and gently position the straps without stretching or applying compression
- Attach one **Accutab®** at the end of each strap position
- Pull each strap to the prescribed pressure indicated by the lines on the **Accutab®**: 20-30, 30-40, or 40-50 mmHg
- **Accutab®** are provided with the **Compreflex®** garment

PRODUCT CODE				
Small	Medium	Large	X Large	XX Large
1491-TB1	1492-TB1	1493-TB1	1494-TB1	1495-TB1



COMPRESHORT®

SHORT & CAPRI

10–15 mmHg



IDEAL FOR

- Lymphedema
- Edema
- Lipedema


COLOUR



Black

FEATURES AND BENEFITS

- The Compreshort® offers mild, supportive compression for the treatment of abdominal or truncal and upper thigh edema. Patients with lymphedema and lipedema benefit from the therapeutic compression and support the garment provides.
- For custom sizes, please contact customer care on our full product portfolio for custom garments.
- Available in custom sizing



Product available
in Short, Capri and
custom lengths



COMPRESHORT® SIZING CHART

(Waist, Hip, Midthigh, Calf)

PRODUCT CODE

	Waist	Hip	Midthigh	Calf	Product SKU (Short)	Product SKU (Capri)
Small	69–76 cm	84–91 cm	41–46 cm	28–32 cm	1301-CS	1301-CP
Small Max	76–86 cm	91–102 cm	41–46 cm	28–32 cm	1301-CM	1301-CPM
Medium	76–86 cm	91–102 cm	46–53 cm	34–40 cm	1302-CS	1302-CP
Medium Max	86–97 cm	102–112 cm	46–53 cm	34–40 cm	1302-CM	1302-CPM
Large	86–97 cm	102–112 cm	53–61 cm	40–46 cm	1303-CS	1303-CP
Large Max	97–112 cm	112–127 cm	53–61 cm	40–46 cm	1303-CM	1303-CPM
X Large	97–112 cm	112–127 cm	61–71 cm	46–53 cm	1304-CS	1304-CP
X Large Max	112–132 cm	127–147 cm	61–71 cm	46–53 cm	1304-CM	1304-CPM
2X Large	112–132 cm	127–147 cm	71–81 cm	53–61 cm	1305-CS	1305-CP
2X Large Max	132–152 cm	147–168 cm	71–81 cm	53–61 cm	1305-CM	1305-CPM
3X Large	132–152 cm	147–168 cm	81–91 cm	61–71 cm	1306-CS	1306-CP
3X Large Max	152–172 cm	168–188 cm	81–91 cm	61–71 cm	1306-CM	1306-CPM
4X Large	152–172 cm	168–188 cm	91–101 cm	71–81 cm	1307-CS	1307-CP
4X Large Max	172–193 cm	188–208 cm	91–101 cm	71–81 cm	1307-CM	1307-CPM
5X Large	172–193 cm	188–208 cm	101–112 cm	81–91 cm	1308-CS	1308-CP

COMPRESHAPE™

SHORT & CAPRI

10–15 mmHg



IDEAL FOR

- Lymphedema
- Edema
- Lipedema

COLOUR




Black

FEATURES AND BENEFITS

- Uses premium material that is designed to be worn as outerwear
- Offers mild, supportive compression for the treatment of abdominal or truncal and upper thigh edema
- The stylish waistband provides abdominal coverage
- The smooth seams give an attractive appearance



Flat seams for a smooth appearance



Product available in Short and Capri (Model is wearing **2303-CS** in Shorts Style)




Model is wearing **2303-CP** in Capris Style

COMPRESHAPE™ SIZING CHART

(Waist, Hip, Midthigh, Calf)

PRODUCT CODE

	Waist	Hip	Midthigh	Calf	Product SKU (Short)	Product SKU (Capri)
Small	69–76 cm	84–91 cm	41–46 cm	28–32 cm	2301-CS	2301-CP
Small Max	76–86 cm	91–102 cm	41–46 cm	28–32 cm	2301-CM	2301-CPM
Medium	76–86 cm	91–102 cm	46–53 cm	34–40 cm	2302-CS	2302-CP
Medium Max	86–97 cm	102–112 cm	46–53 cm	34–40 cm	2302-CM	2302-CPM
Large	86–97 cm	102–112 cm	53–61 cm	40–46 cm	2303-CS	2303-CP
Large Max	97–112 cm	112–127 cm	53–61 cm	40–46 cm	2303-CM	2303-CPM
X Large	97–112 cm	112–127 cm	61–71 cm	46–53 cm	2304-CS	2304-CP
X Large Max	112–132 cm	127–147 cm	61–71 cm	46–53 cm	2304-CM	2304-CPM
2X Large	112–132 cm	127–147 cm	71–81 cm	53–61 cm	2305-CS	2305-CP
2X Large Max	132–152 cm	147–168 cm	71–81 cm	53–61 cm	2305-CM	2305-CPM
3X Large	132–152 cm	147–168 cm	81–91 cm	61–71 cm	2306-CS	2306-CP
3X Large Max	152–172 cm	168–188 cm	81–91 cm	61–71 cm	2306-CM	2306-CPM
4X Large	152–172 cm	168–188 cm	91–101 cm	71–81 cm	2307-CS	2307-CP
4X Large Max	172–193 cm	188–208 cm	91–101 cm	71–81 cm	2307-CM	2307-CPM
5X Large	172–193 cm	188–208 cm	101–112 cm	81–91 cm	2308-CS	2308-CP





MICROFIBER SHADES
in Pink Stripe, Royal Blue Argyle
and Onyx Stripe



Care for Yourself

Why choose **Consumer Health**?

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease, prevention of swelling and varicose veins
- Simplified sizing

CONSUMER HEALTH PRODUCTS

Series	Product	Page
146 / 186	Casual Cotton	70
146/T / 186/T	Traveno Travel Socks	71
151 / 191	Sea Island Cotton	72
152 / 192	All Season Merino Wool	73
153 / 193	Linen	74
143 / 183	Microfiber Shades	75
120	Sheer Fashion	76
170	Soft Silhouette Leggings	78
401	Athletic Recovery Socks	80
412V	Performance Sleeves	81
421	Merino Outdoor Socks	82
160	Eversoft Diabetic Socks	83



CASUAL COTTON

146/186

COMPRESSION LEVEL

15–20 mmHg

FEATURES AND BENEFITS

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation

FIBER CONTENT

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

Product Code Example: **146CA99**

Series	146	Casual Cotton
Style	C	Calf
Size	A	According to shoe size
Colour	99	Black

SIZING CHART

146	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
186	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS

Closed Toe = ●

	WHITE 00	NAVY 10	BROWN 11	KHAKI 30	BLACK 99
Calf (C)	●	●	●	●	●



TRAVENO TRAVEL SOCKS

146/T | 186/T

COMPRESSION LEVEL

15–20 mmHg

FEATURES AND BENEFITS

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation
- Box cover promotes Traveno as a Travel friendly sock

FIBER CONTENT

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

Product Code Example: **186CB99/T**

Series	186	Traveno Travel Socks
Style	C	Calf
Size	B	According to shoe size
Colour	99	Black
Travel	/T	Yes


SIZING CHART

146/T	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
186/T	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS

Closed Toe = ●

Calf (C)	●
----------	---



SEA ISLAND COTTON 151/191

COMPRESSION LEVEL

15–20 mmHg

Also available in 20–30 mmHg, see page 20

FEATURES AND BENEFITS

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

FIBER CONTENT

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free

Product Code Example: **151CA99**

Series	151	Sea Island Cotton (Women)
Style	C	Calf
Size	A	According to shoe size
Colour	99	Black

SIZING CHART

151	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
191	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
Size		A	B	C

STYLES AND COLOURS

Closed Toe = ●

	NAVY 10	BROWN 11	BLACK 99
Calf (C)	●	●	●



ALL-SEASON MERINO WOOL

152/192

COMPRESSION LEVEL

15–20 mmHg

Also available in 20–30 mmHg, see page 22

FEATURES AND BENEFITS

- High quality extra-fine Merino Wool makes this sock soft and breathable
- Moisture wicking capabilities and naturally odour resistant
- Thermoregulating for year-round wear and keeps feet dry
- Smooth, flat toe seam

FIBER CONTENT

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free

Product Code Example: **152CC99**

Series	152	All-Season Merino Wool (Women)
Style	C	Calf
Size	C	According to shoe size
Colour	99	Black

SIZING CHART

152	5–7	7.5–9.5	10–12	10–12
	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	20–27 cm (8"–10.5")
192	6.5–8.5	9–11	11.5–14	11.5–14
	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS

Closed Toe = ●

	NAVY 10	BLACK 99
Calf (C)	●	●

LINEN 153/193

COMPRESSION LEVEL

15–20 mmHg

Also available in 20–30 mmHg, see page 24

FEATURES AND BENEFITS

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

FIBER CONTENT

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free

Product Code Example: **153CB207**

Series	153	Linen (Women)
Style	C	Calf
Size	B	According to shoe size
Colour	207	Lavender

SIZING CHART

153	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
193	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
Size		A	B	C

STYLES AND COLOURS

Closed Toe = ●
For women = **W**
For men = **M**

	BROWN 111	DENIM 201	LIGHT GREY 205	LAVENDER 207
Calf (C)	●	●	●	●
Gender	W/M	W/M	W/M	W

MICROFIBER SHADES 143/183

COMPRESSION LEVEL

15–20 mmHg

Also available in 20–30 mmHg, see page 34

FEATURES AND BENEFITS

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Excellent breathability and moisture management
- Glides on and off with ease
- Smooth, flat toe seam

FIBER CONTENT

- 87% Nylon, 13% Spandex
- Latex free

Product Code Example: **183CB48**

Series	183	Microfiber Shades (Men)
Style	C	Calf
Size	B	According to shoe size
Colour	48	Royal Blue Argyle

SIZING CHART

143	Women's Shoe	5–7	7.5–9.5	10–12
	Ankle Circumference	18–23 cm (7"–9")	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")
183	Men's Shoe	6.5–8.5	9–11	11.5–14
	Ankle Circumference	19–25 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	A	B	C

STYLES AND COLOURS

Closed Toe = ●
For women = **W**
For men = **M**

	MARINER STRIPE 02	PINK STRIPE 15	GRAPHITE STRIPE 17	ONYX STRIPE 19	GRAPHITE HEATHER 27	PURPLE ARGYLE 46	GRAPHITE ARGYLE 47	ROYAL BLUE ARGYLE 48
Calf (C)	●	●	●	●	●	●	●	●
Gender	W	W	M	W/M	W/M	W	W/M	W/M



SHEER FASHION

120

COMPRESSION LEVEL

15–20 mmHg

FEATURES AND BENEFITS

- Incredible sheer design in a variety of fashionable colours
- Sheer toe and heel area
- Comfortable stretch and soft fabric, excellent for everyday wear

FIBER CONTENT

- 79% Nylon, 21% Spandex
(Content can vary slightly by style)
- Latex free

Product Code Example: **120NB33**

Series	120	Sheer Fashion
Style	N	Thigh
Size	B	See thigh sizing chart (next page)
Open toe		If yes = O; If no = blank
Colour	33	Honey



15–20 mmHg
Closed Toe Calf
Size C in Charcoal

120CC12



15–20 mmHg
Open Toe Calf
Size A in Taupe

120CAO29



15–20 mmHg
Closed Toe Pantyhose
Size A in Black

120PA99

PANTYHOSE SIZING CHART

cm / ft	Kg/Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
145 cm / 4'9"																												
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

CALF / THIGH SIZING CHART

120	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23 cm (7"-9")	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")
	Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs ; Size B: <180Lbs ; Size C: <230Lbs.

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	CHARCOAL 12	TAUPE 29	HONEY 33	GOLDEN 36	CAFÉ 73	BLACK 99
Calf (C)	●	●○	●○	●○	●	●
Thigh (N)	●	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●	●



SOFT SILHOUETTE LEGGINGS 170

COMPRESSION LEVEL

15–20 mmHg

FEATURES AND BENEFITS

- Advanced soft stretch microfiber fabric shapes legs and smooths skin
- Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts

FIBER CONTENT

- 83% Nylon, 17% Spandex
- Latex free

Product Code Example: **170LB99**

Series	170	Soft Silhouette
Style	L	Leggings
Size	B	See leggings sizing chart (next page)
Colour	99	Black



Non-binding ankle cuff
stays comfortably in place
for all day wear

STYLES AND COLOURS

Open Foot = ○

BLACK
99

Leggings (L)	○
--------------	---

LEGGING SIZING CHART

cm / ft	Kg / Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

Guidelines for Ankle Circumference: Size A: 18–23 cm (7"–9"), Size B & C: 20–26 cm (8"–10"), Size D: 22–29 cm (8.5"–11.5")

ATHLETIC RECOVERY SOCKS 401

COMPRESSION LEVEL

15–20 mmHg

FEATURES AND BENEFITS

- drirelease® yarns provide moisture wicking and odour control
- Comfortable and easy to put on and take off
- Helps reduce lactic acid buildup
- Mild cushioning in sole, heel, and toe area for added comfort
- Smooth, flat toe seam

FIBER CONTENT

- 67% drirelease® (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

Product Code Example: **401CM00**

Series	401	Athletic Recovery Socks
Style	C	Calf
Size	M	Medium
Colour	00	White

SIZING CHART

401	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X	2

STYLES AND COLOURS

Closed Toe = ●

	WHITE 00	BLACK 99
Calf (C)	●	●

15–20 mmHg

Closed Toe Calf Size M
in Black

401CM99



PERFORMANCE SLEEVES 412V

COMPRESSION LEVEL

20–30 mmHg

FEATURES AND BENEFITS

- Stabilizes and supports active muscles
- Ideal for running, triathlons, endurance sports, and gym workouts
- Fabric resists deterioration in salt, chlorine or fresh water
- Quick drying
- High-visibility reflective transfer

FIBER CONTENT

- 75% Nylon, 25% Spandex
- Latex free

Product Code Example: **412VL58**

Series	412	Performance Sleeves
Style	V	Calf Sleeve
Size	L	Large
Colour	58	Purple

SIZING CHART

412V	Ankle Circumference	18–21.5 cm (7"–8.5")	21.5–25.5 cm (8.5"–10")	25.5–29 cm (10"–11.5")	29–33 cm (11.5"–13")
	Calf Circumference	28–39 cm (11"–15.5")	31.5–45 cm (12.5"–17.5")	35.5–51 cm (14"–20")	40.5–58 cm (16"–23")
	Size	S	M	L	X

STYLES AND COLOURS

Open foot = ○

	WHITE 00	BLUE 50	PURPLE 58*	BLACK 99
Calf Sleeve (V)	○	○	○	○

* Purple is not available in X-Large

MERINO OUTDOOR SOCKS 421

COMPRESSION LEVEL

15–20 mmHg

Also available in 20–30 mmHg, see page 44

FEATURES AND BENEFITS

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned soles help prevent blisters
- Mesh flex zones for comfort and ventilation
- Smooth, flat toe seam
- Soft and comfortable
- Perfect for outdoor work environments and activities

FIBER CONTENT

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free

Product Code Example: **421CS12**

Series	421	Merino Outdoor Socks
Style	C	Calf
Size	S	Small
Colour	12	Charcoal

SIZING CHART

421	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X

STYLES AND COLOURS

Closed Toe = ●

	CHARCOAL 12	OLIVE 31
Calf (C)	●	●

EVERSOFT DIABETIC SOCKS 160

COMPRESSION LEVEL

8–15 mmHg

FEATURES AND BENEFITS

- Padded sole, heel, and toe for comfort and protection
- Smooth, flat toe seam
- Compression-free foot
- drirelease® yarns provide moisture wicking and odour control
- Suitable for individuals with limited dexterity or lower tolerance for compression
- Ideal for individuals with sensitive skin

FIBER CONTENT

- 67% drirelease® (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

Product Code Example: **160CS00**

Series	160	Eversoft Diabetic Socks
Style	C	Calf
Size	S	Small
Colour	00	White

SIZING CHART

160	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14
	Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")
	Size	S	M	L	X

STYLES AND COLOURS

Closed Toe = ●

WHITE
00

Calf (C)

●





Make Life Easier

Why choose **Accessories**?

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage
- Extra liners for added comfort

ACCESSORIES

Product	Page
Gloves	86
Cotton Underliner	87
Washing Solution	87
It Stays® Roll-On Adhesive	87
Sim-Slide®	88
Magnide On/Off®	88
Sigvaris Doff N' Donner® ROLLY®	89
Sigvaris Doff N' Donner® Cone	89
Melany®	89



GLOVES

Textile Gloves*

These high-quality gloves have a latex surface to provide a good grip when putting on and taking off compression stockings.




PACKAGE: 12 PER PACK

SIZE	PRODUCT CODE
X-Small	593R400XS
Small	593R400S
Medium	593R400M
Large	593R400L
X-Large	593R400X
ASST'D 4x XS/S/M	593R400A

Latex-Free Gloves

Gripping donning gloves recommended for use in a latex-free environment.



PACKAGE: 12 PER PACK

SIZE	PRODUCT CODE
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X

TEXTILE GLOVE SIZING CHART - Measure Length and Width of the Hand

Hand Length**	20.5–21 cm	21.5–22 cm	22.5–23 cm	23.5–24 cm	24.5–25 cm
Hand Width	7–7.5 cm	7.5–8 cm	8–8.5 cm	9–9.5 cm	10–10.5 cm
	X-Small	Small	Medium	Large	X-Large

LATEX-FREE GLOVE SIZING CHART - Measure Circumference of the Hand at the Palm

Hand Circumference	16.5–18 cm (6.5"–7")	19–20.5 cm (7.5"–8")	21.5–23 cm (8.5"–9")	24–25.5 cm (9.5"–10")
	Small	Medium	Large	X-Large

* Caution: Textile Gloves contain natural rubber latex that may cause allergic reaction.

** Hand Length: measured from the tip of the longest finger to the crease under the palm.



LINERS

Cotton Underliner

10–15 mmHg

Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10–15 mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings.

PACKAGE: 2 PER PACK

SIZE	PRODUCT CODE	NATURAL RUBBER SIZING
Small / Short	588USS	S1/S3
Small / Long	588USL	S2/S4
Medium / Short	588UMS	M1/M3
Medium / Long	588UML	M2/M4
Large-XL / Short	588ULS	L1/L3 - X1/X3
Large-XL / Long	588ULL	L2/L4 - X2/X4



WASHING SOLUTION

Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

PACKAGE OPTIONS

CONTAINS	PRODUCT CODE
Six 16 oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975



IT STAYS!®

Roll-on body adhesive

A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin.

IT STAYS!® is pliable and moves with the body.

PACKAGE OPTIONS

SIZE	PRODUCT CODE
One 2 oz Bottle (60ml)	581I350
Twelve 2 oz. Bottles (60ml)	581I300

DONNING AIDS

Sim-Slide®

(For open toe Styles)



SIZE	PRODUCT CODE
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL



See how to use
the *Sim-Slide*®
bit.ly/SIGssyt

SIZING CHART For Sim-Slide®

Women's Shoe Size	≤4-6	6.5-8	8.5-9.5	≥10
Men's Shoe Size	≤3-4.5	5-6.5	7-8	≥8.5
Size	S	M	L	XL

Magnide On/Off®

(For closed toe Styles)



SIZE	PRODUCT CODE
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL



See how to use
the *Magnide On/Off*®
bit.ly/49sxcYX

SIZING CHART For Magnide On/Off®

Women's Shoe Size	≤8.5	9-12	≥ 12.5
Men's Shoe Size	≤7.5	8-11	≥ 11.5
Size	M	L	XL

* XL suitable for extra-wide feet

DONNING AIDS

SIGVARIS **Doff N' Donner®** ROLLY®

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Latex free.



SIZE	PRODUCT CODE
One Size	358496



See how to use Sigvaris
DOFF N' DONNER® ROLLY®
bit.ly/3QSBckU

SIGVARIS **Doff N' Donner®** Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner® ROLLY®. The bottom of the base firmly attaches to most clean, flat surfaces via suction.



SIZE	PRODUCT CODE
One Size	358497

MELANY®

MELANY® from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.*



SIZE	CALF CIRCUMFERENCE	PRODUCT CODE
X-Small	25–30 cm	579SXS
Small	30–38 cm	579SSM
Medium	38–46 cm	579SMD
Large	46–56 cm	579SLG

* Individual components not sold separately



See how to use
the **MELANY®**
bit.ly/3SCQYLu





Help yourself

How can **resources** help?

- Learn about the levels of compression for various indications to better help your clients
- Understand the terms and icons that we use
- Follow specific instructions on how to measure for the perfect fit



DOWNLOAD CENTER

Download our catalogues, lookbooks and digital swatches

RESOURCES

Section	Page
Marketing Tools	92
Fitter Certification	93
Symptoms & Indications	94
Glossary of Terms	95
How to Measure	96
How to Code	98





Subscribe to our **NEWSLETTER** to have access to the **latest news and updates!**

✓ Updates


✓ New products

✓ Expert tips

✓ Exclusive offers



Our New MARKETING ASSET GUIDELINE



Get news on the latest product launches and campaigns

Download our collection of pictures

Get some caption ideas for your social media posts





Earn your Phlebology Certificate
by completing our

CERTIFIED FITTER TRAINING COURSE.

FITTER CERTIFICATION

What you'll learn :

- ✓ Venous anatomy & function
- ✓ Chronic venous conditions (CVI, DVT)
- ✓ Compression therapy principles & benefits
- ✓ Circular knit stockings & inelastic wraps
- ✓ Measuring & fitting ready-to-wear stockings
- ✓ Donning & doffing techniques
- ✓ Patient care best practices

Upon completion you will receive:

- Certificate
- Fitter pin
- Complimentary pair of compression socks



JOIN HERE



MOH - Medical Online Hub



Explore our **Medical Online Hub** for educational resources and events on lymphatic and venous disorders, and to deepen your understanding of compression therapy.



SYMPTOMS & INDICATIONS

SYMPTOMS / INDICATIONS	COMPRESSION LEVEL AT THE ANKLE
□ Tired, achy legs □ Occupational or evening edema □ Leg discomfort from long hours of standing or sitting □ Prevent varicosities & swelling during pregnancy □ Preventative for legs predisposed to risk □ Reduce swelling during travel □ Prevent DVT during travel	15–20 mmHg
□ Heavy, fatigued, aching legs □ Mild edema in lower extremities □ Mild varicosities with minimal edema □ Mild varicosities during pregnancy □ Lymphedema upper extremities □ Post-procedure of small veins □ Prevention or management of DVT <i>Note: A prescription is recommended for stockings with the indications of 20 mmHg or higher</i>	20–30 mmHg
□ Moderate varicosities with mild to moderate edema (also during pregnancy) □ Post-procedure of larger veins to maintain treatment success □ CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations) □ Management and treatment of DVT or Post-Thrombotic Syndrome □ Superficial phlebitis □ Lymphedema after decongestant therapy to maintain reduction □ Orthostatic / postural hypotension □ Post-fracture, post-traumatic edema	30–40 mmHg
□ Lymphedema after decongestant therapy to maintain reduction □ Severe Post-Thrombotic Syndrome □ Severe varicosities and/or edema	40 mmHg or higher
□ Severe arterial insufficiency, intermittent claudication, ischemia □ Uncontrolled congestive heart failure	Contraindications
□ Skin sensitivities or allergies □ Advanced neuropathy with sensory loss in leg or foot □ Diabetes with arterial compromise □ Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician □ Eczema, hypodermatitis, dermatosclerosis, dermatitis	Cautions
The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.	
<i>No liability accepted for non-observance of contraindications and cautions.</i> Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.	

TERM	DEFINITION	SIGVARIS QUALITY STANDARD
CHRONIC VENOUS DISEASE (CVD)	A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Disease (CVD), skin issues, and venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-ambulatory patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers Latex free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place.	SIGVARIS GROUP follows ISO Standard 13485:2016 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural fiber fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fiber in Soft Opaque 840, Microfiber Patterns 830, Microfiber Shades 143/183, and Microfiber 820.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in Thermoregulating Wool 422, Merino Outdoor Socks 421, Merino Wool 240, and All-Season Merino Wool 152/192.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our Sea Island 220 and Sea Island Cotton 151/191.
SENSINNOV®	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in Sheer 780, Soft Opaque 840, and Microfiber 820 to help keep thigh-high stockings in place all day long.
SUPIMA COTTON	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber.	SIGVARIS GROUP uses this premium cotton yarn in the Cotton 230, Sea Island 220, Sea Island Cotton 151/191, Casual Cotton 146/186, and Cushioned Cotton 360, making them comfortable and ideal for those with sensitive skin.

HOW TO MEASURE

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

CALF STYLE Open/Closed Toe	PANTYHOSE AND THIGH-HIGH STYLE Open/Closed Toe
1 Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement.	1 Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement.
2 Measure the calf circumference at the fullest part of the calf.	2 Measure the calf circumference at the fullest part of the calf.
3 Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg.	3 Skip Step 3 (but refer to General Guideline F on page 97).
	4 Measure the thigh at its widest circumference.
	5 Determine the leg length by measuring from the floor up to the gluteal fold.
	6 Determine the hip circumference by measuring the widest part of the hips.



Ankle circumference
directly above ankle bone



Calf circumference
at fullest part of the calf



Calf length from the floor to
the right-angle bend of the knee

GENERAL GUIDELINES

A	Measurements should always be taken as early in the morning as possible, when swelling is minimal.
B	Measure against bare skin, not over clothing (Steps 1-5).
C	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images, page 96 & below, for examples).
D	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
E	Length measurements should be taken from the ground up.
F	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.
G	Record shoe size; it may be helpful in determining the best possible fit.



Widest **circumference**
of the **thigh**



Leg length from the
floor up to the gluteal fold



Circumference
at the **hip**

HOW TO CODE

HOW TO CODE MEDICAL AND SPECIALTY PRODUCTS



Product Code Example: 232CMLM99 (page 10)

Series Number	Compression Level	Product Style	Size	Gender / Open Toe	Colour
23 = Cotton (See individual pages, use first 2 digits of series number)	1 = 15–20 mmHg 2 = 20–30 mmHg 3 = 30–40 mmHg 4 = 40–50 mmHg* 5 = 50–60 mmHg*	C = Knee-high or Calf (See individual pages for available styles)	ML = Medium Long (See individual pages for sizing)	M = Men W = Women O = Open toe	99 = Black (See individual pages for available colours)
23	2	C	ML	M	99

* 40–50 mmHg and 50–60 mmHg available in Natural Rubber only

To order Plus Calf with Grip-Top or Plus Size Pantyhose (860P for women only):
add “/P” after the colour in the ordering code

To order Calf Style with Grip-Top:
add “/S” after the colour in the ordering code (see product pages for availability)

To order Thigh with Waist Attachment:
add “/L” to indicate the Left leg or “/R” to indicate the Right leg after the colour in the ordering code (860W and 500W only)

To order Motion High Tech, Motion Thermoregulating Wool and Motion Flow Tech®:
The letter for Gender is not used in the ordering code.

Example: 412CSM59 = Motion High Tech, 20–30 mmHg, Closed toe Calf, Size Small Medium, Colour Red

HOW TO ORDER CONSUMER HEALTH PRODUCTS



Product Code Examples:

151CA99 (page 72)

191CB99 (page 72)

Series Number	Product Style	Size	Colour
151/191 = SEA ISLAND COTTON First 3 digits indicate gender 151 = woman 191 = man <i>(See individual pages for series number)</i>	C = Calf <i>(See individual pages for available styles)</i>	According to shoe size A= women's shoe sizes 5-7 B = men's shoe sizes 9-11 <i>(See individual pages for sizing)</i>	99 = Black <i>(See individual pages for available colours)</i>
151	C	A	99
191	C	B	99

ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS, and EVERSOFIT DIABETIC SOCKS are UNISEX: the first 3 digits are the same for both genders.

Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20 mmHg, Closed toe Calf, Size Large, Colour Olive

Sigvaris, SIGVARIS GROUP, PFS2, Sensinno, Compreflex, Coolflex, Flow Tech, Compreboot, Accutab, Doff N' Donner, ROLLY and MELANY are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2025 SIGVARIS, Inc.

Spacer Fabric and MeshLink are trademarks of SIGVARIS AG, CH-9014, St. Gallen/Switzerland, in many countries worldwide.

Lycra is a registered trademark of INVISTA. Supima is a registered trademark of the Supima Association of America Corporation. SENSIL is a registered trademark of Nilit. All rights reserved. drirelease is copyright ©2021 drirelease. drirelease is an Optimizer Brands company. All rights reserved. Sim-Slide and Magnide are registered trademarks of Arion Holding B.V. IT STAYS is a registered trademark of New Horizons Diversified Interests.

A full-page photograph of two men running on a paved path outdoors. The man on the left is Black, wearing a black t-shirt, white shorts, and a blue and black knee brace on his right knee. The man on the right is white, wearing a dark blue t-shirt, dark shorts, and grey knee-high socks. They are both running towards the camera. The background consists of lush green trees and a brick building. The lighting is bright, suggesting a sunny day.

mobilis
by SIGVARIS GROUP

**SOFT ORTHOPEDIC
SUPPORTS AND BRACES**

INDEX

Intro

102

Mode of action

103

How to measure

103

Elbow

104

EPICARE - Elbow Support

105

EPIBAND - Elbow Strap

105

EPIACTIVE - Elbow Brace

106

Wrist and thumb

107

MANUCARE - Wrist Support

108

MANUWRAP - Wrist and Thumb Brace

108

MANUACTIVE - Wrist Brace

109

MANUSUPPORT - Wrist Brace

109

MANUPLUS SUPPORT - Wrist and Thumb Stabilizer

110

RHIZOSUPPORT - Thumb Brace

110

Back

111

LUMBOCARE - Back Support

112

LUMBOSUPPORT - Lumbar Brace

112

LUMBOSUPPORT EVO - Lumbar Brace

113

LUMBOINTENSE - Lumbar Brace

113

Knee

114

GENUCARE - Knee Support

115

GENUWRAP - Knee Wrap

115

GENUBAND - Knee Strap

116

GENUACTIVE STAB - Stabilizing Knee Brace

116

GENUACTIVE PAD - Patellar Knee Brace

117

Ankle

118

MALLEOCARE - Ankle Support

119

MALLEOACTIVE PAD - Ankle Brace

119

MALLEOACTIVE BAND - Ankle Brace

120

MALLEOSUPPORT - Ankle Brace

120

MALLEOINTENSE - Ankle Stabilizer

121

How to order Mobilis

121

MOBILIS: SOFT ORTHOPEDIC SUPPORTS AND BRACES

A collection of high-quality supports and braces from SIGVARIS GROUP designed to provide support, stability, and pain relief while allowing individuals to stay active.




Mobilis product names are made of two parts:

1 Body part where the product is used


	LINE
Elbow	→ EPI
Hand	→ MANU
Thumb	→ RHIZO
Back	→ LUMBO
Knee	→ GENU
Ankle	→ MALLEO

2 Support level:




Example 1 + 2 = **MALLEO SUPPORT**

MODE OF ACTION




Compression effect

- Decreases swelling and pressure on joints
- Provides pain relief
- Improves mobility



Massage effect with padding

- Targeted swelling reduction
- Improves blood circulation
- Aids healing



Joint stabilization

- Rigid elements stabilize the joint and improve muscle activity



Proprioception

- Aids in perception of movements and body position



Warming effect


- Increases comfort and may provide pain relief

HOW TO MEASURE

Elbow



Measure the arm circumference at 3 cm below the elbow (two fingers)




Measure the elbow circumference with 30 degrees of inclination

Wrist and thumb




Measure wrist circumference

Back



Measure waist circumference at belly button

Knee



Measure knee circumference with 30 degrees of inclination



Measure thigh circumference at 15 cm above the knee

Ankle



Measure ankle circumference above ankle bone

Elbow



	EPI CARE	EPI BAND	EPI ACTIVE
INDICATIONS			
Arthritis	○		○
Joint fluid and swelling	○		○
Chronic irritation/ Tendomyopathy	○		○
Post trauma, Post surgery	○		○
Tennis/golfer's elbow (Epicondylopathy)	○	○	○
SUPPORT LEVEL	 Support level mild	 Support level medium	 Support level medium

EPI CARE

Elbow Support




Contains natural rubber latex that may cause allergic reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)

ACTION MODE




FEATURES

- Mild support
- Maximum comfort
- Discreet look
- Can be worn on either elbow
- Optimal adhesion due to silicone band

COMPOSITION

64% Polyamide, 27% Elastodiene, 5% Elastane, 4% Silicone



Size	Circumference	SKU
Size XS	20-23 cm	MECE-XS-66
Size S	23-26 cm	MECE-SM-66
Size M	26-29 cm	MECE-MD-66
Size L	29-32 cm	MECE-LG-66
Size XL	32-35 cm	MECE-XL-66

EPI BAND

Elbow Strap



INDICATIONS

Tennis/golfer's elbow (epicondylopathy)

ACTION MODE




FEATURES

- Adjustable strap to modify support level
- Silicone pad provides relief for tendon
- Can be worn on either elbow

COMPOSITION

89% Polyamide, 11% Elastane



Size	Circumference	SKU
One size	20 - 35 cm	MEBE-91

EPIACTIVE

Elbow Brace



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)

ACTION MODE




FEATURES

- Elbow brace with silicone inserts
- Removable elastic strap to adjust support
- Anatomical shape and 3D knitting for maximum comfort
- Can be worn on either elbow

COMPOSITION

56% Thermoplastic elastomer, 32% Polyamide, 7% Elastodiene, 3% Elastane, 2% Polyester



Watch a
product video



Size	Circumference	SKU
Size XS	20 - 23 cm	MEAE-XS-91
Size S	23 - 26 cm	MEAE-SM-91
Size M	26 - 29 cm	MEAE-MD-91
Size L	29 - 32 cm	MEAE-LG-91
Size XL	32 - 35 cm	MEAE-XL-91



Wrist & Thumb



	MANU CARE	MANU WRAP	MANU ACTIVE	MANU SUPPORT	MANU PLUS SUPPORT	RHIZO SUPPORT
INDICATIONS						
Arthritis	○	○	○	○	○	○
Chronic irritation/ Tendomyopathy	○	○	○	○	○	○
Carpal tunnel syndrome			○	○	○	
Instability			○	○	○	○
Post trauma, Post surgery		○	○	○	○	
Rhizarthrosis					○	○
Skier's thumb/Torn thumb ligament (Lateral metacarpophalangeal ligament lesion)					○	○
SUPPORT LEVEL						
	 Support level mild	 Support level moderate	 Support level medium	 Support level strong	 Support level strong	 Support level strong

MANUCARE

Wrist Support



INDICATIONS

Arthritis, chronic irritation/tendomyopathy

ACTION MODE




FEATURES

- Mild support
- Soft compression
- Ultimate comfort
- Discreet look
- Can be worn on either wrist

COMPOSITION

66% Polyamide, 29% Elastodiene,
5% Elastane



Measure wrist circumference



Size	Circumference	SKU
Size XS	13 - 15 cm	MMCW-XS-66
Size S	15 - 17 cm	MMCW-SM-66
Size M	17 - 19 cm	MMCW-MD-66
Size L	19 - 21 cm	MMCW-LG-66
Size XL	21 - 23 cm	MMCW-XL-66

MANUWRAP


Wrist and Thumb Brace



INDICATIONS

Arthritis, chronic irritation/tendomyopathy,
post trauma, post surgery

ACTION MODE




FEATURES

- Moderate adjustable support level
- Breathable fabric
- Can be worn on either hand

COMPOSITION

96% Polyamide, 4% Elastane



Measure wrist circumference



Size	Circumference	SKU
Sizes S/M	13 - 17 cm	MMWW-SM/MD-91
Sizes L/XL	17 - 23 cm	MMWW-LG/XL-91


MANUACTIVE

Wrist Brace

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

ACTION MODE




FEATURES

- Brace with silicone inserts around cubitus bone
- Metal rod to limit wrist movement
- Elastic strap allows support level adjustment
- Anatomically shaped 3D knit for comfort and fit
- Specific for left or right hand

COMPOSITION

56% Polyamide, 25% Thermoplastic Elastomer, 9% Elastodiene, 5% Polyester, 5% Elastane



Contains natural rubber latex that may cause allergic reaction



Watch a product video



Size	Circumference	SKU Left Hand	SKU Right Hand
Size XS	13 - 15 cm	MMAW-XS-L-91	MMAW-XS-R-91
Size S	15 - 17 cm	MMAW-SM-L-91	MMAW-SM-R-91
Size M	17 - 19 cm	MMAW-MD-L-91	MMAW-MD-R-91
Size L	19 - 21 cm	MMAW-LG-L-91	MMAW-LG-R-91
Size XL	21 - 23 cm	MMAW-XL-L-91	MMAW-XL-R-91


MANUSUPPORT

Wrist Brace

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

ACTION MODE



FEATURES


- Immobilization/stabilization brace: rigid palm structure, external and internal rods to restrain wrist movements
- Three hook-and-loop straps to adjust the fit
- Can be worn on either wrist

COMPOSITION

92% Polyamide, 5% Elastane, 3% Cotton



Watch a product video



Size	Circumference	SKU
Sizes S/M	13 - 17 cm	MMSWS-SM/MD-91
Sizes L/XL	17 - 23 cm	MMSWS-LG/XL-91

MANUPLUS SUPPORT


Wrist & Thumb Stabilizer



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

ACTION MODE




FEATURES

- Immobilization/stabilization brace for wrist and thumb
- The two hook-and-loop straps can be clamped by laces to adjust the fit
- Specific for right or left hand

COMPOSITION

50% Polyester, 46% Polyamide, 4% Polyurethane



Measure wrist circumference



Size	Circumference	SKU Left Hand	SKU Right Hand
Size XS	13 - 15 cm	MMPSW-XS-L-91	MMPSW-XS-R-91
Size S	15 - 17 cm	MMPSW-SM-L-91	MMPSW-SM-R-91
Size M	17 - 19 cm	MMPSW-MD-L-91	MMPSW-MD-R-91
Size L	19 - 21 cm	MMPSW-LG-L-91	MMPSW-LG-R-91
Size XL	21 - 23 cm	MMPSW-XL-L-91	MMPSW-XL-R-91

RHIZOSUPPORT


Thumb Brace



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, instability, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

ACTION MODE



FEATURES


- Immobilization/stabilization brace for thumb due to external and internal rods
- Two hook-and-loop straps to adjust the fit
- Can be worn on either hand

COMPOSITION

92% Polyamide, 8% Elastane




Watch a product video



Measure wrist circumference



Size	Circumference	SKU
Sizes S/M	13 - 17 cm	MRST-SM/MD-91
Sizes L/XL	17 - 23 cm	MRST-LG/XL-91



Back



	LUMBO CARE	LUMBO SUPPORT	LUMBO SUPPORT EVO	LUMBO INTENSE
INDICATIONS				
Muscular tension	○	○	○	○
Lower back pain (Lumbar syndrome)	○	○	○	○
Muscular weakness and loss (degeneration)		○	○	○
Post surgery		○	○	○
Sciatica		○	○	○
SUPPORT LEVEL	 Support level mild	 Support level strong	 Support level strong	 Support level extra strong

LUMBOCARE

Back Support



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome)

ACTION MODE




FEATURES

- Adjustable belt provides mild support
- Comfortable material that provides a warming effect
- Can be worn against the skin

COMPOSITION

65% Cotton, 31% Polyester, 4% Elastodiene



Size	Circumference	SKU
Size S	70 - 95 cm	MLCB-SM-66
Size M	95 - 115 cm	MLCB-MD-66
Size L	115 - 135 cm	MLCB-LG-66

LUMBOSUPPORT


Lumbar Brace



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



FEATURES


- Four stabilizing rods to support and secure
- Adjustable straps for a precise fit
- To be worn over clothing or underwear
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

54% Polyamide, 21% Elastane,
20% Polypropylene, 5% Polyester



Watch a
product video



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size S	70 - 95 cm	MLSB-SM-S-91	MLSB-SM-T-91
Size M	95 - 115 cm	MLSB-MD-S-91	MLSB-MD-T-91
Size L	115 - 135 cm	MLSB-LG-S-91	MLSB-LG-T-91

LUMBOSUPPORT EVO

Lumbar Brace



Support level
strong


Watch a
product video



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



FEATURES


- Lumbar spine stabilizing belt with four posterior rods
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

44% Polyamide, 23% Polyester, 8% Viscose, 7% Elastodiene, 6% Ethylene vinyl acetate, 5% Polyvinyl chloride, 4% Polyurethane, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size S	56 - 71 cm	MLSEB-SM-S-99	MLSEB-SM-T-99
Size M	68 - 84 cm	MLSEB-MD-S-99	MLSEB-MD-T-99
Size L	81 - 96 cm	MLSEB-LG-S-99	MLSEB-LG-T-99
Size XL	94 - 109 cm	MLSEB-XL-S-99	MLSEB-XL-T-99
Size XXL	107 - 122 cm	MLSEB-2X-S-99	MLSEB-2X-T-99

LUMBOINTENSE

Lumbar Brace



Support level
extra strong


Watch a
product video



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



FEATURES


- Stabilizing belt with four posterior rods provides stronger pressure on the lumbar spine
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

47% Polyamide, 27% Polyester, 2% Viscose, 4% Elastodiene, 6% Ethylene-acetate, 5% Polyvinyl chloride, 4% Polyurethane, 2% Polypropylene, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



Measure waist circumference at belly button



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size S	56 - 71 cm	MLIB-SM-S-99	MLIB-SM-T-99
Size M	68 - 84 cm	MLIB-MD-S-99	MLIB-MD-T-99
Size L	81 - 96 cm	MLIB-LG-S-99	MLIB-LG-T-99
Size XL	94 - 109 cm	MLIB-XL-S-99	MLIB-XL-T-99
Size XXL	107 - 122 cm	MLIB-2X-S-99	MLIB-2X-T-99

Knee



	GENU CARE	GENU WRAP	GENU BAND	GENU ACTIVE STAB	GENU ACTIVE PAD
INDICATIONS					
Arthritis	○	○		○	○
Joint fluid and swelling	○			○	○
Chronic irritation/ Tendomyopathy	○	○	○	○	○
Post trauma, Post surgery	○	○		○	○
Kneecap (patellar) pain		○	○	○	○
Deep knee pain (Femoro-patellar syndrome)		○	○	○	○
Instability, Hyperlaxity, Meniscus lesion, Gonarthrosis				○	
SUPPORT LEVEL	 Support level mild	 Support level moderate	 Support level medium	 Support level medium	 Support level medium

GENUCARE

Knee Support



Support level
mild



Contains natural rubber latex that may cause allergic reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery

ACTION MODE




FEATURES

- Mild support
- Optimal adhesion due to silicone band
- Discreet look
- Can be worn on either knee

COMPOSITION

64% Polyamide, 28% Elastodiene
5% Elastane, 3% Silicone



Size	Circumference	SKU
Size XS	30 - 34 cm	MGCK-XS-66
Size S	34 - 38 cm	MGCK-SM-66
Size M	38 - 42 cm	MGCK-MD-66
Size L	42 - 46 cm	MGCK-LG-66
Size XL	46 - 50 cm	MGCK-XL-66



GENUWRAP

Knee Wrap



Support level
moderate

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




FEATURES

- Moderate support knee wrap with three straps for optimum fit
- Adjustable size
- Breathable
- Can be worn on either knee

COMPOSITION

94% Polyamide, 6% Elastane



Size	Circumference	SKU
Size S/M	30 - 42 cm	MGWK-SM/MD-91
Size L/XL	42 - 55 cm	MGWK-LG/XL-91



GENUBAND

Knee Strap



INDICATIONS

chronic irritation/tendomyopathy, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




FEATURES

- Adjustable knee strap to modify support level
- Silicone insert is located on kneecap tendon
- Can be worn on either knee
- Anatomical shape for extra comfort

COMPOSITION

89% Polyamide, 11% Elastane



Size	Circumference	SKU
One Size	30 - 55 cm	MGBK-91

GENUACTIVE STAB


Stabilizing Knee Brace



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome), instability, hyperlaxity, meniscus lesion, gonarthrosis

ACTION MODE



FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap with three Velcro straps
- The anatomically shaped fabric is compressive
- Wearable on both knees

COMPOSITION


8% Polyamide, 23% Elastomer, 15% Elastodiene, 6% Polyester, 5% Elastane, 2% Polyvinyl chloride, 1% Silicone



Contains natural rubber latex that may cause allergic reaction



Watch a product video



Size	Knee Circumference	Thigh Circumference	SKU
Size XS	30 - 34 cm	40 - 44 cm	MGSK-XS-91
Size S	34 - 38 cm	44 - 48 cm	MGSK-SM-91
Size M	38 - 42 cm	48 - 52 cm	MGSK-MD-91
Size L	42 - 46 cm	52 - 56 cm	MGSK-LG-91
Size XL	46 - 50 cm	56 - 60 cm	MGSK-XL-91
Size 2X	50 - 54 cm	60 - 64 cm	MGSK-2X-91


GENUACTIVE PAD

Patellar Knee Brace

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/ tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE



FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap
- Anatomically shaped 3D knit for comfort and fit

COMPOSITION

52% Polyamide, 26% Thermoplastic Urethane, 14% Elastodiene, 5% Polyester, 2% Elastane, 1% Silicone




Support level
medium



Contains natural rubber latex that may cause allergic reaction



Watch a
product video



Measure knee circumference below the kneecap with 30 degrees of inclination, and thigh circumference 15 cm above the knee


Size	Knee Circumference	Thigh Circumference	SKU
Size XS	30 - 34 cm	40 - 44 cm	MGAK-XS-91
Size S	34 - 38 cm	44 - 48 cm	MGAK-SM-91
Size M	38 - 42 cm	48 - 52 cm	MGAK-MD-91
Size L	42 - 46 cm	52 - 56 cm	MGAK-LG-91
Size XL	46 - 50 cm	56 - 60 cm	MGAK-XL-91
Size 2X	50 - 54 cm	60 - 64 cm	MGAK-2X-91

3D Knit
Provides breathability and comfort for both everyday and physical activities.

Lateral Supports
Provides high stability

Silicone Insert
The anatomical shape helps the kneecap to remain stable

Bend Zone
Knitted fabric is less compressive for more comfort and flexibility





Ankle



	MALLEO CARE	MALLEO ACTIVE PAD	MALLEO ACTIVE BAND	MALLEO SUPPORT	MALLEO INTENSE
INDICATIONS					
Arthritis	○	○	○	○	○
Joint fluid and swelling	○	○	○		
Chronic irritation/ Tendomyopathy	○	○	○	○	○
Post trauma, Post surgery, Soft tissue lesion	○			○	○
Hyperlaxity				○	○
Chronic instability				○	○
Post trauma, Post surgery, Ligament injury				○	○
Moderate sprains				○	○
SUPPORT LEVEL	 Support level mild	 Support level medium	 Support level medium	 Support level strong	 Support level extra strong

MALLEOCARE

Ankle Support

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion

ACTION MODE



FEATURES

- Mild support
- Adjustable level of stability
- Can be worn on either ankle

COMPOSITION

70% Polyamide, 12% Elastane,
9% Elastodiene, 8% Polyester,
1% Polypropylene



Contains natural rubber latex that may cause allergic reaction



Taille	Circumference	SKU
Size S	18 - 22 cm	MMCA-SM-66
Size M	22 - 25 cm	MMCA-MD-66
Size L	25 - 28 cm	MMCA-LG-66


MALLEOACTIVE PAD

Ankle Brace

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy

ACTION MODE



FEATURES

- 2 lateral pads stabilize and support the ankle
- Gentle massaging effect
- Breathable 3D knit for optimal fit and comfort
- Left and right foot option
- **Includes a foot slip for easier donning**

COMPOSITION


45% Polyamide, 40% Thermoplastic
Urethane, 10% Elastodiene, 3% Elastane,
2% Polyester



Contains natural rubber latex that may cause allergic reaction



**Watch a
product video**



Size	Circumference	SKU Left Ankle	SKU Right Ankle
Size S	16 - 21 cm	MMAA-SM-L-91	MMAA-SM-R-91
Size M	21 - 26 cm	MMAA-MD-L-91	MMAA-MD-R-91
Size L	26 - 31 cm	MMAA-LG-L-91	MMAA-LG-R-91
Size XL	31 - 37 cm	MMAA-XL-L-91	MMAA-XL-R-91

MALLEOACTIVE BAND

Ankle Brace



Support level
medium



Contains natural
rubber latex that
may cause allergic
reaction




Watch a
product video

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/
tendomyopathy

ACTION MODE




FEATURES

- The elastic band allows for adjustable support
- The anatomically shaped breathable knit fabric is compressive
- Can be worn in a shoe
- Can be worn on either ankle

COMPOSITION

62% Polyamide, 18% Elastodiene, 12% Viscose,
5% Polyester, 3% Elastane



Measure ankle circumference
above ankle bone

Size	Circumference	SKU
Size XS	18 - 20 cm	MMBA-XS-91
Size S	20 - 23 cm	MMBA-SM-91
Size M	23 - 25 cm	MMBA-MD-91
Size L	25 - 27 cm	MMBA-LG-91
Size XL	27 - 30 cm	MMBA-XL-91



MALLEOSUPPORT

Ankle Brace




Support level
strong

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma,
post surgery, soft tissue lesion, hyperlaxity, chronic
instability, ligament injury, moderate sprains

ACTION MODE




FEATURES

- Stabilizing ankle brace with lateral rigid rods
- Two hook-and-loop straps for added support
- The lace closure system allows for an adjustable fit
- Can be worn on either ankle

COMPOSITION

91% Polyamide, 5% Elastane,
2% Polyester, 2% Cotton



Measure ankle circumference
above ankle bone

Size	Circumference	SKU
Size S	19 - 22 cm	MMSA-SM-91
Size M	22 - 26 cm	MMSA-MD-91
Size L	26 - 30 cm	MMSA-LG-91



MALLEOINTENSE

Ankle Stabilizer



Support level
extra strong



Contains natural
rubber latex that
may cause allergic
reaction




Watch a
product video

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion, hyperlaxity, chronic instability, ligament injury, moderate to severe sprains

ACTION MODE




Movement control | Stabilization

FEATURES

- Immobilizing ankle stabilizer based on two rigid shells
- Adaptable fit with two gel cushions to provide comfort and cooling effect
- Two Velcro straps as well as heel support strap provide ability to adjust tightness and support
- Can be worn on either ankle

COMPOSITION

50% Gel, 22% Polyvinyl chloride, 20% Polypropylene, 7% Polyamide, 1% Neoprene



Size	Height	SKU
One Size	≥ 150 cm	MMIA-91



HOW TO ORDER MOBILIS PRODUCTS

M	E	A	E	-	MD	-	-	91
Mobilis	Letter matches first part of the product name: E = Epi M = Manu MP = Manu Plus R = Rhizo L = Lumbo G = Genu M = Malleo	Letter matches line of product: C = Care W = Wrap A = Active S = Support or Stab B = Band I = Intense SE = Support Evo	Letter matches body part where the product is used: E = Elbow W = Wrist T = Thumb B = Back K = Knee A = Ankle	-	Size	-	When your product has a left or right side, L or R is indicated. When your product has a short or tall height, S or T is indicated.	Available colours: 66 = Beige 91 = Graphite 99 = Black

Add **L** to indicate the Left side or **R** to indicate the Right side after the Size in the ordering code. **L** or **R** is not specified in products that fit both left and right. **Example:** MMAW-MD-L-91.

For back braces, add **S** to indicate the Short height or **T** to indicate the Tall height after the Size in the ordering code.

You can find products without a size. **Example:** MEBE-91

You can find products with combined sizes. **Example:** MRST-SM/MD-91

Simplify
care with
WrapAbility

sigvaris

SIGVARIS
GROUP

mobilis

by SIGVARIS GROUP

Sigvaris Inelastic Wraps

The innovative alternative to bandages

Designed for daily use, Wraps are clinically proven to reduce edema and pain, enable easy self-application, and provide lasting comfort and independence.



Get more done faster:

CHAT WITH US LIVE

on the Online Dealer Portal:

sigvaris-online.com

or text us at **1.800.363.4999**



*"Wraps take less time to put on,
and you don't need anyone to help you.
The Oscar of compression goes to the Wraps!"*

Dr. Heim, General Surgeon, Switzerland



 sigvarisgroup_canada

 sigvarisgroup.canada

 SIGVARIS GROUP Canada

TO ORDER:

sigvaris-online.com
canadaorders@sigvaris.com
1.800.363.4999



DOWNLOAD CENTER
Download our **MCE Catalogue**
& **MCE Lookbook**